



Product Service

CERTIFICATE

No. Z1 08 04 42097 037

Holder of Certificate: Invacare Rea AB
Org. No. 556098-9211
Växjövägen 303
343 21 Diö
SWEDEN

Production Facility(ies): 42097

Certification Mark:



Product: Active wheelchairs

Model(s): XLT

Parameters: Technical Data: see enclosure

Tested according to: DIN EN 12182:1999
DIN EN 12183:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. The certification mark must not be altered in any way. See also notes overleaf.

Test report no.: 71334697

Date, 2008-04-18

Page 1 of 2





Product Service

Enclosure to certificate no.: Z1 08 04 42097 037

Technical data for the active wheelchair
Type: XLT

	XLT Active	XLT Dynamic	XLT Swing
total length incl. footrests:	830 – 1030 mm	730 – 880 mm	820 – 1230 mm
total length without footrests:	-	-	675 – 825 mm
max. width:	555 – 680 mm		555 – 705 mm
max. height:	610 – 1050 mm		
seat width:	355, 380, 405, 430, 455, 480 mm		355, 380, 405, 430, 455, 480, 505 mm
seat depth:	360 – 400, 410 – 450, 460 - 500 mm		
seat height:	400 – 520 mm	390 – 510 mm	390 – 520 mm
backrest height:	200 – 490 mm		
seat plane angle:	0° - 14°		
backrest angle:	-10° - 10°		
max. load:	135 kg		
brake system:	knee lever brake (user brake)		
back wheels:			
- size	22 and 24 x 1 ³ / ₈ " , 24, 25 and 26 x 1"		
- tyres	Air filled / accident proof		
- adjustable	yes		
- fixation	Quick release		
front wheels:			
- size	ø 75 x 24, ø 100 and 125 x 27, ø 140 x 37, ø 150 x 30, ø 150 x 32, ø 180 x 40 mm		
- tyres	Air filled / accident proof		
- adjustable	yes		
- fixation	rigid		
frame material:	titanium		
armrest:	height adjustable and removable		
legrest:	height adjustable	height adjustable and removable	
anti tip device	adjustable		

Department: MHS5 / me
 Date: 2008-04-18

