



Yes, you can.®



COMPATIBILITY CATALOGUE 2011



Confidence in posture solutions

Content

➤	How to read the catalogue	3
➤	Product presentation	
	Product overview	4
	All sizes	5
	All colours	7
➤	Common parts	
	Compatibility cross list	8
	Armrests	12
	Backrest	16
	Seat cushions	20
	Legrests and Hangers	22
	Castors	26
	Rearwheels	29
	Others	33
➤	Individual parts	
	<i>Rea</i> ™ Azalea ®	35
	<i>Rea</i> ™ Bellis ™	45
	<i>Rea</i> ™ Clematis ®	51
	<i>Rea</i> ™ Focus ™	55
	<i>Rea</i> ™ Spirea ² ®	61
	<i>Rea</i> ™ Spirea ³ ®	65
	<i>Rea</i> ™ Spirea ⁴ ®	69
	<i>Invacare</i> ® Spin ^x	76
	<i>Invacare</i> ® XLT ™	81

Introduction

The purpose of the catalogue is to describe the range and give an overview of the cross compatibility of our current passive Tilt & Recline wheelchairs; Rea™ **Azalea**®, Rea™ **Bellis**™ and Rea™ **Clematis**®, our medium-active wheelchairs; Rea™ **Focus**™, Rea™ **Spirea**²®, Rea™ **Spirea**³®, Rea™ **Spirea**⁴_{NG} and Invacare® **Spin**z, as well as of our new active wheelchair Invacare® **XLT**.

How to read the Compatibility Catalogue:

Product Presentation

Here you will find presentations of all wheelchairs for adults that are produced at Invacare Rea.

Common Parts

The chapter starts with a Compatibility cross list, where the whole product range is listed on a modular level, with what fits to what and which parts are compatible with each other. The order as well as the parts are according to the European Prescription Form (EPF).

All compatible parts (i.e. they fit two or more wheelchairs) are described together with article numbers

and pictures in the coming sections.

Individual Parts

This section contains a short presentation of the individual wheelchairs followed by descriptions of how the individual wheelchairs have been measured.

You will also find a complete table with all seat heights and presentations of all the parts that only fit each particular wheelchair, in other words that are not compatible with the rest of the range.

Product presentation

Product Overview



Wheelchair version		Self propelling	Transit version.	Tiltable seat	Max user weight kg	Angle adj. frame	Reclining backrest	Angle adj backrest
Rea™ Azalea ® / Azalea ® Max	Passive wheelchair with adaptable seating	×	×	×	135 / 180	-	×	-
Rea™ Bellis ™ / Bellis ™ 150	Passive wheelchair with adaptable seating	×	×	×	125 / 150	-	×	-
Rea™ Clematis ®	Passive wheelchair basic seating	×	-	×	125	-	×	-
Rea™ Azalea ® Base	Passive wheelchair	×	×	×	135	-	-	×
Rea™ Focus ™ / Focus ™ 150	Medium active wheelchair with extensive adjustability	×	-	-	135 / 150	×	× / -	×
Rea™ Spirea ²® / Rea™ Spirea ²® T	Medium active wheelchair with basic adjustments	×	×	-	125	×	-	-
Rea™ Spirea ³®	Medium active wheelchair with extensive adjustability	×	-	-	135	×	-	×
Invacare® Spin ˣ™	Medium active wheelchair with extensive adjustability	×	-	-	125 / 135 *	×	-	×
Rea™ Spirea ⁴®	Medium active wheelchair with extensive adjustability	×	×	-	135	×	-	×
Invacare® XLT ™ 75/90 Invacare® XLT ™ Swing	Medium active wheelchair with extensive adjustability	×	-	-	135	×	-	×
Invacare® XLT ™ Max	Medium active wheelchair with extensive adjustability	×	-	-	180	×	-	×

* **Spin**ˣ: seat width 48, and 50,5 = 135 kg



Size and adjustments

All sizes. For specific information, see under each product chapter, "Individual Parts"

Seat widths in mm

Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [*]	XLT _™
+0 / +50 / +100 mm adj.	± 20 mm adj.	+ 20 mm adj.	+ 25 mm adj.	-	-	-
						330 *
				355 ****	355 (spec)	355
340	410	390	380	380	380	380
390	460	440	405	405	405	405
440	510	490	430	430	430	430
490			455	455	455	455
550 **			480	480	480	480
610 **			505	505	505	505 *
			530 ^			555 ***
						605 ***

** Azalea Max

^ Focus 150

**** = Spirea⁴

* XLT Swing + XLT Max

*** XLT Max

Seat depths in mm

Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [*]	XLT _™
+70 mm adj.	+120 mm adj.	+60 mm adj.	+75 mm adj.	-50 mm adj.	-50 mm adj.	- 40 mm adj.
380	420	420	Low 375cm	400 mm	400	400
430			High 425cm	450 mm	450	450
500				Adj 400 - 500 ***	500	500 ^ ^
530 **						550 *
600 **						600 *

** Azalea Max

*** Spirea⁴

^ ^ XLT + XLT Max

* XLT Max

Possible seat heights in mm

Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [*]	XLT _™
400 / 450 / 500	360-485	450	Low 365-490	380-505 *	410-530	Swing 380-520
			High 415-540			75/90 400-520
						Max 440-520

* Spirea⁴ = 380 - 510

Backrest heights in mm

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [*]	XLT _™	XLT _™ Max
Backrest plate	610 + 75 + 50	600 + 75 + 50	610 + 70					
Tension adjustable	560 **	550						
Angle adjustable				380 - 460	350-390/ 400-480 *	360-400/ 410-490	300 - 490	400 - 490
					350-480/ 480-550 ^			
Fixed					380 - 460 ***	350 - 410	200 - 350	
Reclining				610				
Biangular				405 - 505				

** Azalea Base backrest poles, height 550 mm

* Spirea³

^ Spirea⁴

*** Spirea² + Spirea⁴

Size and adjustments

Backrest angles							
	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin[±]	XLT_±
Backrest reclining	0-30°	0-30°	-1° - +32°	0-30°	-	-	-
Backrest Angle adj.	0-30° (Base)	-	-	-6° - +13°	0 ** / +12° ***	±12°	±10°
	-	-	-	-	-6° - +13° *	-	-
Biangular Lower	-	-	-	-6° - +13°	-	-	-
Biangular Upper	-	-	-	+/- 25°	-	-	-

* Spirea⁴ ** Spirea² *** Spirea³

Seat angle							
Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin[±]	XLT_±
Tiltable	- 1° - + 20° / 25°	- 1° - + 15°	- 1° - + 19°	-	-	-	-
Adjustable	-	-	-	0 - + 6°	0 - + 6°	0 - + 6°	0 - + 14°

Legrests length in mm Seat footplate (without seat cushion). Total seat depth							
Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin[±]	XLT_±
Fixed hangers and footplates 80° with short tube	300 - 420	320 - 440	-	320 - 440	320 - 440	330 - 450	290 - 410 *
Fixed hangers and footplates 80° with long tube	390 - 510	410 - 530	-	410 - 530	410 - 530	420 - 540	390 - 510 *
Fixed hangers and footplates 90° with short tube	290 - 410	290 - 410	-	340 - 460	-	320 - 440	280 - 360 *
Fixed hangers and footplates 90° with long tube	400 - 520	400 - 520	-	450 - 570	-	420 - 540	400 - 490 *
Angle adjustable hanger and footplates							
With short tube	300-420	320-440	300-420	320-440	320-440	340-460	290-410 *
With long tube	390-510	410-530	390-510	410-530	410-530	430-550	400-510 *
Central footplate	26-53	26-53	-	-	-	-	-
Legrest amputee	250	250	-	250-320	250-330	250-330	250-310 *
Plaster legrest	-	-	-	570-730	620-780	-	-
Footbow, short	-	-	-	-	-	-	380-400 ** 360-380 ****
Footbow, long	-	-	-	-	-	-	410-480 ** 390-450 ****
Footbow, rigid XLT Max	-	-	-	-	-	-	380-480
Footrest One piece, Flip-up	-	-	-	-	-	-	380-480 *** 360-450 ****
Footrest One piece, short tube	-	-	-	-	-	-	280-400 *
Footrest One piece, long tube	-	-	-	-	-	-	400-490 *
Footrest High mounted	-	-	-	-	-	-	270-370 *** 270-350 ****

* XLT Swing ** XLT Active + XLT Max *** XLT Active **** XLT Dynamic

All Colours Upholstery

Colour	Colour code	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ^x	XLT ₊
 Grey Dartex	TR23	x	x	x	x	-	-	-
 Grey Plush	TR32	x	x	-	x	x *	-	-
 Black Jemina	TR18	-	-	-	x	x	x	x
 Black Nylon	TR 33	-	-	-	-	-	x	x

* Not for Spirea⁴

All Colours Frame

Colour	Colour code	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3} ®	Spin ^x	XLT ₊	Spirea ⁴ ®
 Pearl Grey	0515	x	x	x	x	x	x	x	-
 Shining Blue	0507	-	-	-	x	-	x	x	-
 Anthracite Bl.	0514	-	-	-	x	-	x	x	-
 Red	0501	-	-	-	-	-	x	x	-
 Lime Green	0517	-	-	-	-	-	x **	-	-
 Orange	0518	-	-	-	-	-	x *	x	-
 Bl. Titanium	0518	-	-	-	-	-	-	x	-
 Slate Grey	0615	-	-	-	-	-	-	-	x

* Special option on Spin-X ** Can only be ordered as "Special"

Cross list for common parts

Backrest type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Backrest plate	x *	x	i	i	-	-	-
Laguna backrest plate	-	-	-	i	-	-	-
Foldable backrest	-	-	-	-	-	x	x
Velcro adjustable biangular	-	-	-	i	-	-	-
Velcro adjustable	x	x	-	x	i	x	x
Sling upholstery	-	-	-	i	i	i	-
Rigid back (Special)							x **

* Not for Azalea® Max

** Not for XLT Max

Backrest Cover and Cushion Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Cover Thin	x *	x	-	x	-	-	-
Cover Lateral	x	x	-	-	-	-	-
Mistral 2	x *	x	-	x	-	x	-
Laguna 2	x	x	x	x	-	-	-
Passad 2	x *	x	-	-	-	-	-
Shoulder High	x *	x	-	-	-	-	-

* Not for Azalea® Max

Backrest options/ pushhandles	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Backrest bar	-	-	-	x	x	x	-
Push bar	x	x	-	-	-	-	-
Pushhandles fixed	-	-	-	i	i	i	i
Pushhandles external long grip	-	-	-	x	x	x	i
Pushhandles external short grip	-	-	-	x	x	x	i
Pushhandles external +10 cm	-	-	-	x	x *	-	-
Push handles braced	x	x	-	-	-	-	-
Push handles height adjustable	-	-	i	-	-	-	i

* Spirea⁴

Seat type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Firm seat, depth adjustable	i	i	i	i	-	-	-
Sliding seat	i	-	-	-	-	-	-
Sling seat	-	-	-	-	x	x	x *
Velcro adjustable seat	-	-	-	-	x	x	x

* Not for XLT Max

Seat Cushions	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Flo-Shape	x	x	x	-	-	-	-
Standard flat	-	-	-	x	x	x	x
Mistral 2 layer	-	-	-	i	-	-	-
Inlay cushion	-	i	-	i	-	-	-
Flo-tech Contour	-	-	-	-	-	-	x *

* XLT Max

Legrests / Hangers	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Central legrest, incl calf pad, footrest and attachment	x *	x	-	-	-	-	-
Hanger angle adjustable incl calf support	x	x	x	x	x	x	x **
Hanger fixed 80 deg	x *	x	-	x	x	x	x **
Hanger fixed 90 deg	x *	x	-	x	-	x	x **
Legrest amputee	x *	x	-	x	x	i	i
Plaster legrest	-	-	-	x	x	-	-
Legrest attachment fixed	i	-	i	i	i	i	i
Legrest attachment width adjustable	i	i	-	-	-	-	-

* Not for Azalea® Max

** XLT Swing

Cross list for common parts

Footrests	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Footrest one piece	x *	-	-	x	x	x	i
Fixed footplate	x *	x	x	x	x	x	x **
Adjustable footplate angle/depth	x	x	x	x	x	x	x **
Adjustable footplate angle	x *	-	x	-	-	x	x **
Footbow rigid (XLT 75/90 & XLT Max)	-	-	-	-	-	-	i
Footrest high mounted (XLT 75/90)	-	-	-	-	-	-	i
* Not for Azalea® Max ** XLT Swing							
Armrest Attachment	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Detachable	x	x	x	i	x	x	i
Flip up	-	i	-	i	i	i	i
Side support/Armrest lower part	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Height adjustable, detachable	x	x	x	x	x	x	-
Height adjustable, detachable Rio	-	-	-	x	x	x	x
Height adjustable, flip up	-	i	-	-	-	-	x *
Sideguard	-	-	-	-	-	i	i
Mudguard	-	-	-	-	-	i	i
Armrest padded	x	-	-	-	-	-	-
* Not for XLT Max							
Armrest Pads	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Short Narrow	-	x	-	x	x	x	x
Long Narrow	-	x	-	x	x	x	x **
Long Narrow padded	-	-	i	-	-	-	-
Long wide padded	x *	-	-	-	-	-	-
Long narrow Nice	-	i	-	-	-	-	-
Long wide Nice	i	-	-	-	-	-	-
Long wide	-	-	x	x	-	-	-
Height adjustable Standard (Duck tail)	i	-	-	-	-	-	-
Hemiplegics wide	-	-	-	x	x	-	-
Hemiplegics narrow	x	x	x	-	-	-	x **
* Only for Azalea Assist ** Not for XLT Max							
Seating options	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Abduction cushion in Dartex	x	x	x	-	-	-	-
Incontinence cover for seatcushion	x	x	-	-	-	-	-
Backrest Options Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3} ®	Spin ⁺	XLT ₊
Neckrest depth adjustable	x	x	x	x	x	-	i
Headrest with cheek support	x	x	x	-	-	-	i
Headrest depth / width adjustable	x	x	x	x	x	-	i
Headrest / Neckrest cover	x	x	x	x	x	-	-
Trunk support multi adjustable with standard arm	x *	x	-	-	-	-	-
Trunk support multi adjustable with extended arm	x *	x	-	-	-	-	-
Trunk support with standard arm	x *	x	-	i	-	-	-
Trunk support with extended arm	x *	x	-	-	-	-	-
Trunk support, swing away	x *	x	x	-	-	-	-
Trunk support multi adjustable swing away	x *	x	x	-	-	-	-
Lumbar support cushion	x	x	-	-	-	-	-
Wedges	x	x	-	x	-	x	x **
* Not for Azalea® Max ** Not for XLT Max							

Cross list for common parts

Legrest / Footrest Options Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±	XLT _± Max
Heel straps	x **	x	x	x	x	x	x	-
Calf pads	x	x	x	x*	x*	x*	x*	-
Calf straps	-	-	-	x	-	x	x	-
Footplate extender	x **	x	-	x	x	-	-	-

** Not for Azalea® Max

* Thinner pad

Armrest Options Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±
Hemiplegics	x *	x	-	x	x	x	x **
Side support cover	x	x	-	-	-	-	-
Cushion to sidesupport	x	x	-	-	-	-	-

* Not for Azalea® Max

** Not for XLT Max

Castor Attachment Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±
Chassis extender	i	-	-	-	-	-	-
Castor attachment fixed	x	x	x	i	i	i	i
Castor attachment quick release	-	i	-	-	-	-	-

Front Castors	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±
75 x 24 medium hard	-	-	-	-	-	-	x ***
100 x 27 medium hard	-	-	-	-	-	x	x ***
125 x 27 medium hard	-	x	-	x	x	-	x ***
140 x 37 soft	x	x	-	x	x	x	x
150 x 27 medium hard	x *	x	-	x	x	x	x
150 x 30 soft	x *	x	-	x	x	x	x
150 x 32 pneumatic	-	x	-	x	x	x	x
180 x 40 pneumatic	-	x	-	x	x **	x	x
180 x 40 soft	-	x	-	x	x **	x	x
200 x 27 medium hard	x	x	-	x	x	x	-
200 x 30 soft	x *	x	-	x	x	-	-
200 x 32 pneumatic	x *	x	-	x	x	x	-
200 x 45 medium hard	x	x	x	x	x	-	-
200 x 50 pneumatic	-	x	-	x	x	-	-

* Not for Azalea® Max

** Not for Spirea⁴

*** Not for XLT Max

Rear Wheel Axel	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±
QR	x	x	x	x	x	x	x
Camber	-	-	-	i	i *	i	i

* Spirea³ + Spirea⁴

Rear Wheels Size/Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±
16"	x	x	-	-	x *	-	-
20"	-	x	-	x	x ^	-	-
22"	x	x	-	x	x	x	x ***
24"	x	x	x	x	x	x	x
24" basic	x ^^	-	x	-	x **	-	-
25"	-	-	-	-	-	x	x ***
26"	-	-	-	-	-	-	x ***

^^ Not for Azalea® Max

* Spirea² + Spirea⁴

** Spirea³ + Spirea⁴

^ Not for Spirea⁴

*** Not for XLT Max

Tyre Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin [±]	XLT _±
Tyres pneumatic non marking	x *	x	-	x	x	x	x **
Tyres pneumatic profiled	x *	x	x	x	x	x	x **
Tyres solid	x	x	x	x	x	x	x
Tyres puncture proof	x *	x	-	x	x	x	x **
Tyres low profile	x *	x	-	x	x	x/i	x/i **
Tyres low profile Black Marathon Evol.	-	-	-	-	-	-	x

* Not for Azalea® Max

** Not for XLT Max

Cross list for common parts

Handrim Type	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Handrim aluminium	x	x	x	x	x	x	x
Handrim stainless	x */**	x *	-	x *	x *	x	x
Handrim plastic coated	x **	x	-	x	x	x	x
Handrim foam coated	x */**	x *	-	x *	x *	x	x
Handrim black grip	x */**	x *	-	x *	x *	x	x

** Not for Azalea® Max

* Not available as std from 1:st of April 2010

Castor Spokeguard	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Spokeguard Frost	x	x	x	x	x	x	x

Safety	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Standard brakes (Malta)	i	i	i	x	x	x	i
Brake lever extension	-	-	-	x	x	x	x
One arm brake	-	x	-	x	x	x	x **
Assistant operated drum brake 22" / 24"	x	x	x	x	x	x	x ***
Assistant operated drum brake 16"	x	x	-	-	x *	-	-
Pelvic belt	x	x	x	x	x	x	x
Reflectors kit (wheels and chassis)	x	x	x	x	x	x	x
Auto lock armrest	x	x	x	x	x	x	x
Anti tip device	i	x	i	x	x	i	i

* Spirea² + Spirea⁴ ** Not for XLT Max

*** Only 24" wheels for XLT Max

Accessories	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Table tray fixed	x	x	x	x	x	x	x **
Table tray swing away	x *	-	-	-	-	-	-
Pump	x	x	x	x	x	i	x
Cane holder	i	i	-	i	i	i	i

* Not for Azalea® Max

** Not for XLT Max

x	Same part
-	Not existing
i	Exists on individual product level

For the parts that are the same, see under next chapter for more info, for those parts that are individual for one wheelchair, see under specific product chapter

> Armrests

Detachable armrests



Height adjustable, light, Rio



Height adjustable, long side plate

Flip-up armrests



Height adjustable for:
Rea™ **Focus**™ and Rea™ **Spirea**^{2/3}®

Armrests to the different wheelchair models*

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Height adjustable, Rio	-	-	-	x**	x	x	x
Height adjustable, long side plate	x	x	x	x	x	x	x ***
Height adjustable flip up	-	i	-	x**	x ^	x	x ***
Height adjustable "comfort"	-	-	-	-	x ^^	-	-

** Not for Focus™ I 50.

^ Not for Spirea4

*** Not for XLT Max

^^ Only for Spirea4

* Note that all front attachments for all armrest attachments have the same insert. That means that all armrests could be inserted in all attachments. However we recommend using of the original version for the respective wheelchair product as stated below.

Armrest attachments

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Detachable	5106289			1431721	1434680		-
Flip-up, front	-	i	-				-
Flip-up, rear	-	i	-	i	1435492	i	i


* Rear attachment to Flip up is always included in the chassie



Armrest attachment
Detachable and Flip-up
Rea™ **Spirea^{2/3/4}®**
Invacare® **Spin⁺**



Armrest attachment
Detachable and Flip-up
Rea™ **Focus™**



Armrest attachment
Detachable Rea™ **Comfort™**
Rea™ **Bellis™**

> Armrests

Armrest pads including poles

I435616	Short narrow left
I435615	Short narrow right
I435618	Long narrow left
I435617	Long narrow right
I518702	"Ducktail" pad left
I518701	"Ducktail" pad right
I445787	Cushion for side support left
I445786	Cushion for side support right
I415453	Side support Cover Dartex
I493220	Hemiplegics armrest left
I493219	Hemiplegics armrest right



Short narrow



Long narrow



"Ducktail" pad



Cushion for side support
2 cm thick



Side support Cover:



Hemiplegics

Armrest, lower part

Art.no.		Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin^x	XLT^o
I415521	Height adjustable, longside plate left	x	x	x	x	x	x	x
I415520	Height adjustable, longside plate right	x	x	x	x	x	x	x
I435614	Height adjustable Flip up left	—	—	—	x	x *	x	x
I435613	Height adjustable Flip up right	—	—	—	x	x *	x	x

* Not for Spirea⁴



Height adjustable, cpl.
light weight, Rio



Height adjustable
lower part



Height adjustable, lower part
Flip up

Armrests pads including poles


Description	Art.no.	
	Left	Right
Height adjustable Standard (Duck tail)	1451563	1451562
Long Wide	1427169	1427170
Padded	1492547	



Ducktail



Long wide



Padded

What armrest pad fits to what armrest

Description	Short narrow	Long narrow	Long shaped	Duck tail
Height adjustable, short side plate Rio	x	x	-	-
Height adjustable, long side plate	-	-	x	x
Height adjustable, flip up Bellis™	-	-	x	-
Height adjustable, flip up Spirea^{2/3}® / Focus™	-	-	-	-

Armrest Rio cpl

Description	Art.no.	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin⁺	XLT₊
Height adjustable short pad left	1435461	-	-	-	x	x	x	x
Height adjustable short pad right	1435460	-	-	-	x	x	x	x
Height adjustable long pad left	1435487	-	-	-	x	x	x	x *
Height adjustable long pad right	1435486	-	-	-	x	x	x	x *

* Not for XLT Max

Armrest Flip-up, Height adjustable cpl

Description	Art.no.	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3}®	Spin⁺	XLT₊
Height adjustable short pad left	1435459	-	-	-	x	x	x	x *
Height adjustable short pad right	1435458	-	-	-	x	x	x	x *
Height adjustable long pad left	1435489	-	-	-	x	x	x	x *
Height adjustable long pad right	1435488	-	-	-	x	x	x	x *

* Not for XLT Max

Armrest, Height adjustable cpl

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3}®	Spin⁺	XLT₊
Height adjustable "ducktail" pad left			1517368 *			-	-
Height adjustable "ducktail" pad right			1517367 *			-	-

* Not for Azalea® Max

Backrests

Backrests



Backrest plate with extension



Tension adjustable backrest



Velcro adjustable for biangular backrest angle

Backrest types and versions

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Backrest plate for reclining backrest ***	x	x	i	i	-	-	-
Backrest tension adjustable for Backrest cushion Passad 2 and Laguna	x	x	-	x	-	-	-
Tension adjustable backrest for reclining backrest ***	x	x	-	-	-	-	-
Sling or Velcro adjustable for fixed backrest angle *	-	-	-	x	i	x	-
Sling or Velcro adjustable for adjustable backrest angle **	-	-	-	x	x	x	-
Velcro adjustable for biangular backrest angle	-	-	-	i	-	-	-
Velcro tension adjustable, fixed height, infinitely variable, angle adjustable, foldable and click in backrest	-	-	-	-	-	-	x
Velcro tension adjustable, fixed height, infinitely variable, angle adjustable, foldable and click in backrest	-	-	-	-	-	-	x

* = Not for Spirea³ ** = Only for Spirea³ and Spirea⁴ *** = Not for Azalea® Max

Backrest extension

Description	Art.no.	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺
Extender backrest	I440147	390	410	—	—	—	—
Extender backrest	I440148	440	460	—	—	—	—
Extender backrest	I440149	490	510	—	—	—	—

Backrests

Backrests cushions to fit below seat width (measurements in mm)

Description	Art.no			Models/Sizes					
Mistral 2	Dartex	Grey Plush	Jemima	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3} ®	Spin ⁺
small	14301668	14301669		340 ***					
				390 ***			380		380 **
medium	142460510	14305404		440 ***	410		430		405 / 430
large	14246052	14305405		490 ***	460		480		455 / 480
x-large	142460530	14305406			510		530		505
Laguna fixed shape									
Minor	1508972	15100660		340					
xx-small	1439937					390 **			
x-small	1427118	1435796		390 ***			380 ****		
small	1427119	1435797		440 ***	410	440 **	430 ****		
medium	1427120	1435798		490 ***	460	490 **	480 ****		
large					510				
Azalea Max	1525073	1525074		550					
Azalea Max	1525075	1525076		610					
Passad 2									
x-small	1429271			390 ***					
small	1424647	1444329		440 ***	410		380-405****		
medium	1424649	1444330		490 ***	460		430-455****		
large	1429278	1444331			510		480-505****		
Shoulder high									
small	1444275			390 ***					
medium	1444276			440 ***	410				
large	1444277			490 ***	460				
x-large	1444278				510				
Backrest cover Thin									
xx-small		5106739-TR32	5106739-TR18				380	380	380 **
x-small		5106740-TR32	5106740-TR18				405	405	405
		1446512 *		390					
small		5106741-TR32	5106741-TR18				430	430	430
		1446513 *		440	410				
medium		5106742-TR32	5106742-TR18				455	455	455
large		1435742	1435523				480	480	480
x-large		5106743-TR32	5106743-TR18						505
		1446522 *		490	460				
xx-large		1435743	1435524		510		530		
Backrest cover Lateral									
x-small	1508865	1444344		390					
small	1508866	1444345		440	410				
medium	1508867	1444346		490	460				
large	1508868	1444347			510				
Azalea Max	1526459	1526573		550					
Azalea Max	1526424	1526572		610					

* Includes velcro adapter 14305619 ** Not plush *** For +50 mm and +100 mm width, see individual parts **** For Focus see individual parts For XLT, see individual parts



Mistral 2



Laguna,
Fixed shape



Passad 2



Shoulder high



Backrest cover Thin



Backrest cover Lateral



Backrest cover

Backrests

What backrest cover and cushion fit what backrest base?

Description	Mistral 2	Laguna	Passad 2	Backrest cover
Backrestplate without extension	×	—	—	—
Backrestplate with extension	—	×	×	—
Plate Laguna	—	×	—	—
Velcro adjustable for fixed backrest angle	×	—	—	×
Velcro adjustable for adjustable backrest angle	×	—	—	×
Velcro adjustable for biangular backrest angle	—	—	—	×
Tension adjustable for reclining backrest	×	×	×	×

Backrest bar/Push handles (measurements in mm)

Description	Art.no.	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺
Backrest bar	510061088-C9005	-	-	-	-	355 *	-
	5106399-C9005	-	-	-	380	380	380
	5221696-C9005	-	-	-	405	405	405
	5106331-C9005	-	-	-	430	430	430
	5106190-C9005	-	-	-	455	455	455
	1435484-9005	-	-	-	480	480	480
	5106191-C9005	-	-	-	50,5	50,5	50,5
	1435490-9005	-	-	-	530 Focus 150	-	-
Push handles rear mounted + attachment long grip	1434743	-	-	-	×	×	×
Push handles rear mounted + attachment short grip	1445825	-	-	-	×	×	×
Push handles rear mounted + attachment +10 cm	1445821	-	-	-	×	×	-
Push bar	1444350	390		390	-	-	-
	5255136	440	410	440	-	-	-
	5255137	490	460	490	-	-	-
	5255138		510		-	-	-
High Push bar	1524533	550 **			-	-	-
	1524529	610 **			-	-	-
Push handles braced	1508953-0509	340			-	-	-
	1420007-0509	390		390	-	-	-
	1420008-0509	440	410	440	-	-	-
	1420009-0509	490	460	490	-	-	-
	1420010-0509		510	-	-	-	-
	1523947	550 **	-	-	-	-	-
	1523946	610 **	-	-	-	-	-

For XLT, see individual parts

** Azalea Max

* Spirea⁴



Backrest bar



Push handles



Push bar



Pushhandles braced

Backrests

Neck- and Headrests

Description	Art.no
Neckrest depth adjustable	1508848
Headrest depth / width adjustable	1420161
Headrest with cheek support	1446560

Cover Neck- and Headrests

Description	Art.no.
Headrest cover	1518516
Neckrest cover	1518517



Neckrest, depth adjustable



Headrest depth / width adjustable



Headrest with cheek support




Cover for headrest

What neck-and headrest fit to what backrest base?

Description	Backrest Plate	Backrest Tension Adjustable	Backrest Plate Laguna	Velcro	Shaped Backrest
Neckrest depth adjustable	x	x	x	x	x
Headrest with cheek support	x	x	x	x	x
Headrest depth/width adjustable	x	x	x	x	x

Lumbar support cushion

Backrest options	Art.no.	Azalea®	Bellis™	Focus™	Spin [®]	XLT™
Lumbar support cushion	5252527	x	x	-	-	-
Wedges	1435491	x	x	x	x	x *



Lumbar support cushion

* Not for XLT Max

Trunk support

Description	Art.no		Azalea®	Bellis™	Clematis®	Focus™	XLT™
	Left	Right					
Trunk support Multi adjustable with standard arm	1448692	1448693	x *	x	-	-	-
Trunk support Multi adjustable with extended arm	1484447	1484448	x *	x	-	-	-
Trunk support swingaway	1448697	1448698	x *	x	x	-	-
Trunk support swingaway	1497859	1497854	-	-	x	-	-
Trunk support with standard arm	1445798		x *	x	x	x	-
Trunk support Multi adjustable swing away	1530457	1530456	x *	x	-	-	-
	1530500	1530482	-	-	x	-	-
	1530502	1530501	-	-	-	-	x **
Trunk support with extended arm	1492608	1492609	x *	x	-	-	-

* Not for Azalea® Max

** Not for XLT Max



Trunk support Multi adjustable



Trunk support swingway




Trunk support




Trunk support with extended arm

Seat cushions

Seat cushions



Seat cushion Flo Shape



Seat cushion Standard flat

Seat cushions to fit below seat width (measurements in mm)

Description	Art.no.			Azalea [®] ****	Bellis [™]	Clematis [®]	Focus [™]	Spirea ^{2/3/4} ®	Spin ^x	XLT [™]
	Jemima	Dartex	Grey Plush							
Flo Shape		1509920	1508976	340						
		1508053	1508054	390 **		390 *				
		1508055	1508056							
		1508055	1508056	440 **	410	440 *				
		1508057	1508058							
		1508057	1508058	490 **	460	490 *				
		1508057	1508058							
		1508059	1508060		510					
Standard seat cushion seat depth 40 cm	5104255-TR18		5104317-TR32					355 ***		355*
	5104317-TR18		5103179-TR32				380	380	380*	380*
	5103179-TR18		5104262-TR32					405	405*	405*
	5104262-TR18		5103180-TR32				430	430	430*	430*
	5103180-TR18		5226509-TR32					455	455*	455*
	5226509-TR18		5104269-TR32				480	480	480*	480*
	5104269-TR18		5104319-TR32					505	505*	505*
Standard seat cushion seat depth 450 mm	5104257-TR18		5104259-TR32					355 ***		355*
	5104319-TR18		5104264-TR32				380	380	380*	380*
	5104259-TR18		5104266-TR32					405	405*	405*
	5104264-TR18		5226510-TR32				430	430	430*	430*
	5104266-TR18		5104271-TR32					455	455*	455*
	5226510-TR18		1435445				480	480	480*	480*
	5104271-TR18							505	505*	505*
Standard seat cushion seat depth 50 cm	1435429						530			
	1512809									355*
	1495830								380*	380*
	1495831								405*	405*
	1495832								430*	430*
	1495833							455 ***	455*	455*
	1495834							480 ***	480*	480*
	1495835							505 ***	505*	505*

* Not in plush

** For +50 mm and +100 mm width, see individual parts

**** For Azalea Max, see individual sheet

*** Spirea⁴

For XLT Max, see individual parts

Seat cushions

Seat cushions to the different wheelchair models

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin^x	XLT^o
Flo-Shape	x	x	x	-	-	-	-
Standard, flat	-	-	x	x	x	x	x
Flo-Tech Countour	-	-	-	-	-	-	x *

* For XLT Max

Seat cushion Flo-Shape (meassurements in mm)

Description	Grey Dartex	Grey Plush	Azalea® **	Bellis™	Clematis®
Seat cushion Flo-Shape	1509920	1508976	340*		
	1508053	1508054	390 *		390
	1508055	1508056	440 *	410	440
	1508057	1508058	490 *	460	490
	1508059	1508060		510	

* For +50 mm and +100 mm width, see individual parts

Seating options (meassurements in mm)

Description	Art.no.	Azalea® **	Bellis™	Clematis®
Abduction cushion in Dartex		1444172		1449897
Incontinence seat cover				
Minor	1510636	340	-	-
Small	5253078-TR26	390*	-	-
Medium	5253079-TR26	440*	410	-
Large	5253080-TR26	490*	460	-
X-Large	5253617-TR26	-	510	-

* For +5 cm and +10 cm width, see individual parts

** For Azalea Max, see individual sheet



Abduction cushion **Bellis™**



Abduction cushion
Clematis®



Incontinence seat cover

Legrests & Footrests

Legrests and Footrests



Fixed legrests



Angle adjustable legrests

Hangers to fit to wheelchair

Description	Art.no.	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin[*]	XLT₊
Fixed 80°	1493055-0509	x	x	-	x	x	x	x**
Fixed 90°	1443056-0509	x	x	-	x	-	x	x*
Angle adjustable left	1508945	-	-	-	x	x	x	x*
Angle adjustable right	1508948	-	-	-	x	x	x	x*
Angle adjustable left	1513007	x	x	x	-	-	-	-
Angle adjustable right	1513006	x	x	x	-	-	-	-

* Only for XLT Swing

** Only for XLT Swing, also available in chassie colour, see individual parts



Fixed hanger 80°



Fixed hanger 90°



Angle adjustable hanger

Legrests & Footrests

Footrests to fit seat width in mm

Description	Art.no.				Width						
	Left		Right								
	Short	Long	Short	Long	<i>Azalea</i> ®	<i>Bellis</i> ™	<i>Clematis</i> ®	<i>Focus</i> ™	<i>Spirea</i> ^{2/3/4} ®	<i>Spin</i> ✕	<i>XLT</i> .
Fixed footrest											
Footrest 19	I428357	I428359	I428358	I428360	***			380 / 405	380 / 405 ^	380 / 405	380 / 405
Footrest 21	I428361	I428363	I428362	I428364	390	410	390	430 / 455	430 / 455 ^	430 / 455	430 / 455
Footrest 24	I428365	I428367	I428366	I428368	440	460	440	480 / 505	480 / 505 ^	480 / 505	480 / 505
Footrest 27	I435756	I435659	I435757	I435660	490	510	490	530			
Footrest 27	I435756	I435659	I435757	I435660		510 *		530 / 55,5*			
Adjustable footrest in depth and angle											
Footrest 12	I485430	I485428	I485524	I485429					355 ^^		355
Footrest 14	I429159	I429171	I429160	I429172				380 / 405	380 / 405	380 / 405	380 / 405
Footrest 16	I429161	I429173	I429162	I429174							
Footrest 17	I429163	I429175	I429164	I429176	340 / 390		390	430 / 455	430	430 / 455	430 / 455
Footrest 18	I429165	I429177	I429166	I426178	440	410	440				
Footrest 19	I429167	I429179	I429168	I429180	490			480	480	480	480
Footrest 20	I429169	I429181	I429170	I429182		460 / 510	490	505	505	505	505
Footrest 20	I440036	I440038	I440037	I440039		510*		530 / 555*			
Adjustable footrests in depth and angle to fit 90 hangers **											
Footrest 116	I491369	I491377	I491370	I491378						405	405
Footrest 141	I491371	I492170	I491372	I492263						455	455
Footrest 166	I491375	I492264	I491376	I492265						505	505
Footrest 177	I446130	I446136	I446131	I446137				380 / 405		380	380
Footrest 202	I446132	I446138	I446133	I446139	340 / 390	410		430 / 455		430	430
Footrest 227	I446134	I446140	I446135	I446141	440	460		480* / 505		480	480
Footrest 247	I446159	I446161	I446160	I446162	490	510*					

* Heavy version 150kg

** Not for heavy version 150kg

*** For +50 mm and +100 mm width, see individual parts

^ For Spirea⁴, see individual sheet

^^ Spirea⁴



Fixed footrest



Adjustable footrest in
depth and angle



Adjustable footrest in
depth and angle
(To fit hanger 90°)

Legrests & Footrests

One piece footplate to fit seat width in mm

Description	Art.no.		Width				
	Short	Long	Azalea ®	Focus ™	Spirea ²³⁴®	Spin ˣ	XLT ™
Adjustable footrests in depth and angle to fit 80 degree hangers **	1496945	1496954		380 / 405	380	380	380*
	1496946	1496955			405	405	405*
	1496947	1496982		430 / 455	430	430	430*
	1496948	1497004			455	455	455*
	1496949	1497005		480 / 505	480	480	480*
	1496953	1497006			505	505	505*
	1495821	1495825		530** / 555**		530	
Adjustable footrests in depth and angle to fit 90 degree hangers **	1491247	1493349		380 / 405	380	380	380*
	1491248	1493350			405	405	405*
	1493267	1493351		430 / 455	430	430	430*
	1493268	1493358	340 / 390		455	455	455*
	1493269	1493359		480 / 505	480	480	480*
	1493270	1493377	440		505	505	505*
	1495821	1495825	490	530** / 555**		530	530*
	1495822	1495826	490				

** Heavy version.

*** For +50 mm and +100 mm width, see individual parts

* For XLT Swing



One piece footplate

Calf support used for angle adjustable legrests and fits in seat width mm

Description	Art.no.		Azalea ®	Bellis ™	Clematis ®	Focus ™	Spirea ²³⁴®	Spin ˣ
	Left	Right						
Calf support 36	151086106						355 *	
Calf support 38	1518617					380 / 405	380	380
Calf support 16	1518531					430 / 455	405 / 430	405 / 430
Calf support 19	1518532					480 / 505	455 / 480	455 / 480
Calf support 22	1518533					530 / 555	505	505
Calf support 16	1516313		340 / 390		390			
Calf support 19	1516314		440	410	440			
Calf support 22	1516315		490	460 / 510	490			
Calf support Azalea ® Max	1525057	1525058	550 / 610					

* Spirea⁴



Calf support **Azalea**®, **Assist**™, **Bellis**™, **Clematis**®



Calf support **Focus**™, **Spirea**²³⁴®, **XLT**™ and **Spin**ˣ



Calf support **Azalea**® **Max**

Legrests & Footrests


Footplate extender for angle adjustable footplate (Excluding Central footplate) (measurements in mm)

Description	Art.no.		Width					
	Left	Right	Azalea ®	Assist ™	Bellis ™	Clematis ®	Focus ™	Spirea ²³⁴®
Extender 135	5255297	5255296	-	-	-	380	380	380 / 405
Extender 160	5255299	5255298	340 - 490	Fit to all	Fit to all	430 / 480 / 530	405 / 430	430 / 455
							455 / 505	480 / 505


Heel straps Fits all footrests in all sizes (Excluding Central footplate)

Description	Art.no.	Azalea ®	Bellis ™	Clematis ®	Focus ™	Spirea ²³⁴®	Spin ˣ	XLT ˣ
Heel straps 40-50	5322013	x	x	x	x	x	x	x *

* Not for XLT Max



Heel strap 40-50



Footplate extender

Legrests, Footrests and others



Central legrest



Legrest amputee



Plaster legrest

Complete legrests

Description	Art.no.		Azalea ®	Bellis ™	Focus ™	Spirea ²³⁴®	Spin ˣ
	Left	Right					
Central legrest	1484821		x				
	1424640			x			
Legrest amputee	1442991	1442990	x	x	x	x *	x
Plaster legrest	1435771	1435770			x	x	

* For Spirea⁴, see individual sheet

For legrest attachment an individual adapter to respective Tilt and Recline wheelchair is needed. One version for Central Hanger and another attachment for the other Hangers. Please read more under the individual Product chapters. Note that in the article number above of Central footplate, the adapter is included

Front castors



100x27 medium hard



125x27 soft



140x37 soft




150x27 medium hard



150x30 soft



150x32 pneumatic



180x40 pneumatic



180x40 soft



200x27 medium hard



200x30 soft



200x32 pneumatic



200x45 medium hard



200x50 pneumatic

Castor attachments


Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Fixed	x	x	x	i	i	i	-
Quick release	-	x	-	-	-	-	-
Fixed with angleadjustment 0-6°	-	-	-	x	x	x	i

Castor forks


Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
75 mm - narrow	-	-	-	-	x	x	i
110 mm - narrow	x	x	-	x	x	x	i
110 mm - wide	x	x	x	x	x	x	i
150 mm - narrow	x	x	-	x	x	x	i
150 mm - wide	x	x	-	x	x	x	i

Castor attachment type


Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea ^{2/3/4} ®	Spin ⁺	XLT ₊
Castor attachment fixed	5254308			i	i	i	i




75 - narrow




110 - narrow



110 - wide



150 - narrow



150 - wide


Article number of castors in combination with forks

Description	Forks	Art.no.	
Castor (Approximate in inch")		Fixed Castor + fork	QR Castor + fork
125 x 27 medium hard (5")	75	5106784	—
	110	5320710	5322102
	150	5320727	5322104
140 x 37 soft (5,5")	110	5322135	5322137
	150	5322136	5322138
150 x 27 medium hard (6")	110	5321788	5322107
	150	5321792	5322109
150 x 30 soft (6")	110	5750808	5750809
	150	5750810	5750811
150 x 32 pneumatic (6")	110	5321788	5322106
	150	5321782	5322108
180 x 40 pneumatic (7")	110	5750789	5750790
	150	5750791	5750792
180 x 40 soft (7")	110	5750784	5750785
	150	5750786	5750787
200 x 27 medium hard (8")	110	5320717	5322112
	150	5320729	5322114
200 x 30 soft (8")	110	5750802	5750803
	150	5750804	5750805
200 x 32 pneumatic (8")	110	5104188	5322111
	150	5750889	5322113
200 x 45 medium hard (8")	110	5322143	5322144
	150	5322147	5322148
200 x 50 pneumatic (8")	110	5322145	5322146
	150	5322149	5322150

For XLT, see individual parts

For Spirea⁴, see individual parts

Rear wheels




Tyre type and size of rearwheel					
Description	16"	20"	22"	24"	24" basic
Tyres low profile (24 x 1 1/8") *	-	-	-	x	-
Tyres low profile (24 x 1")	-	-	-	x	-
Tyres pneumatic non marking	-	x	x	x	-
Tyres pneumatic profiled	x	x	x	x	x
Tyres solid	x	-	x	x	x
Tyres puncture proof	-	-	x	x	-

- Not compatible with this model


* Not standard on all models from 1:st of April 2010




low profile



pneumatic non marking



pneumatic profiled



solid



puncture proof


Rear wheels

Handrims and size of rearwheel				
Description	16"	20"	22"	24"
Handrim aluminium	-	1530271 **	1531579 **	1531577 **
Handrim stainless *	-	-	-	1531305 **
Handrim plastic coated	-	1530272 **	1531580 **	1531578 **
Handrim foam coated *	-	-	-	1531459 **
Handrim black grip *	-	-	-	1531460 **

- Not compatible with this model

* Not standard on all models from 1:st of April 2010

** Available from 1:st of April 2010



Aluminium



Stainless



Plastic coated



Foam coated



Black grip

Rear wheels

Article numbers of rearwheels in combination with tyres and handrims (for 16" see next chapter under Drum brakes)

Description	Tyres						
Handrim	Rear wheel size	Low profile	Pneumatic non marking	Pneumatic profiled	Solid	Puncture proof	Low profile XLT™ *** and Spin^x
Aluminium	20"	-	1523885 **	1523944 **	-	-	-
	22"	1524020 **	1523948 **	1523949 **	1524021 **	1524022 **	-
	24"	1524049 **	1524029 **	1524030 **	1524050 **	1524150 **	1530008 **
Aluminium (Basic)	24"	-	-	1492642	1492643	-	-
Stainless *	24"	1531063 **	1530763 **	1530764 **	1530762 **	1531129 **	-
Plastic coated	22"	1524026 **	1524023 **	1524024 **	1524027 **	1524028 **	-
	24"	1524187 **	1524163 **	1524186 **	1524188 **	1524189 **	-
Foam coated *	22"	-	-	-	-	-	-
	24"	1531533 **	1531462 **	1531508 **	1531531 **	1531532 **	-
Black grip *	24"	1531566 **	1531538 **	1531563 **	1531564 **	1531565 **	-


* Not standard on all models from 1:st of April 2010

** Available from 1:st of April 2010

*** Not for XLT Max

Spoke guards with article numbers

Description	16"	20"	22"	24"	24" basic
Spokeguard Frost	—	5102795	5102800	5102807	5102807



Frost

Article numbers of rearwheels for drum brakes in combination with tyres and handrims ***

Description	Rear wheel size	Tyres			
		Pneumatic non marking	Pneumatic profiled	Solid	Puncture proof
Handrim Aluminium	22"	1524289 **	1524290 **	1524291 **	1524292 **
	24"	1524436 **	1524437 **	1524447 **	1524512 **
Stainless *	22"	-	-	-	-
	24"	1531064 **	1531080 **	1531127 **	1531128 **
Plastic coated	22"	1524337 **	1524361 **	1524388 **	1524435 **
	24"	1524513 **	1524514 **	1524515 **	1524516 **
Black grip *	24"	1531567 **	1531568 **	1531569 **	1531570 **

* Not standard on all models from 1:st of April 2010 ** Available from 1:st of April 2010

*** For complete drum brake kit see Spare part catalogue



Rear wheel for drum brake

Rear wheels

Rear wheels for drumbrakes							
Description	Azaled [®]	Bellis [™]	Clematis [®]	Focus [™]	Spirea ^{2/3/4} [®]	Spin [×]	XLT [™]
Ass.operated drumbrake 22 / 24"	x	x	x	x	x	x	x **
Ass.operated drumbrake 16"	x	x	-	-	x *	-	-

* Spirea²® ** Only 24" wheels for XLT Max



Assistant operated drum
brake 22" / 24"



Assistant operated drum
brake 16"

Brakes incl. Attachment

Description	Art.no.		Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin^x
	Left	Right						
Standard brakes	1485423	1485424	x					
Standard brakes	1446041	1444342			x			
Standard brakes (Malta)	1415485	1415486					x	x
Standard brakes, Focus	1435800	1435799				x		
Standard brakes	5255117	5255124		x				
Brake lever extension	5203340	5203340				x	x	x
One arm brake	1492107	1492000					x	x
One arm brake	1449900	1449899				x		
One arm brake	5255418	5255426		x				
User manouvered brake 16" Bellis	1446121	1446122		x				
User manouvered brake 16" Spirea ²	1446111	1446112					x	



Standard brakes (Malta)



Brake lever extension



One arm brake

Safety

Description	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin^x	XLT™
Pelvic belt	5252276 ***				1435769 ^	1483883	5252276
Auto lock armrest	5721060						
Anti tip device	1439933	5255306	**	1435470	5255306 ^	i	i


i = see individual part
** Part of frame

* Transport version see individual part
*** For Azalea Max see individual part

^ For Spirea⁴ see individual parts



Pelvic belt



Reflectors kit



Auto lock armrest



Anti tip device

Table Tray (measurements in mm)

Description	Art.no.	Azalea® **	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin^z	XLT.
Table tray fixed small	5206727			390	380 / 405	355 ^ / 380 / 405	380 / 405	355 / 380
Table tray fixed medium	5206726		410	440	430 / 455	430 / 455	430 / 455	405 / 430
Table tray fixed large	5206725		460	490	480 / 505	480 / 505	505	455 - 505
Table tray fixed X-large	5207643		510		530 *			

^ Spirea⁴

* = Focus 150

** For Azalea see individual part



Table tray fixed

Cane holder


Description	Art.no.						
	Azalea®	Bellis™	Clematis®	Focus™	Spirea^{2/3/4}®	Spin^z	XLT.
Cane holder	5295015		-	1435764	5104964	5106390	

Step tube

Description	Art.no.						
	Azalea®	Bellis™	Clematis®	Focus™	Spirea³®	Spin^z	XLT.
Step tube	-	-	-	5103276	-	1493346	-



Cane holder



Step tube Focus, Spirea 3



Step tube XLT, Spin-X



Rea™

Azalea® – Individual Parts

Content

➤ Rea™ Azalea®	35
➤ This is how Rea™ Azalea® was measured	37
➤ Individual accessories to Rea™ Azalea®	38



Rea™

Azalea® – Individual Parts

Offers excellent seating comfort and many adjustment possibilities!

The **Azalea®** is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the user to alter the seat position, ensuring maximum comfort at all times. This combined with good driving characteristics makes **Azalea®** the ideal choice. **Azalea®** is a self propelling wheelchair with large rear wheels, and is complimented with an alternative transit version with 16" rear wheels.



Rea™ **Azalea®**



Rea™ **Azalea® Minor**



Rea™ **Azalea® Assist**



Rea™ **Azalea® Tall**



Rea™ **Azalea® Base**



Rea™ **Azalea® Max**




Rea™


Azalea® – Individual Parts

This is how Rea™ Azalea® has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



* As the dimensions are in whole cm steps it might in reality mean slightly more or less adding a couple of mm



Rea™

Azalea® – Individual Parts

Sizes	Azalea®	Azalea® 16”	Azalea® Tall	Azalea® Base	Azalea® Minor	Azalea® Max
Seat dimensions (measurements in mm)						
Seat width	390 / 440 / 490 +50 / 100 mm			340 / 390 / 440 / 490 +50 / 100 mm	340 +50 / 100 mm	550 / 610 +50 / 100 mm
Seat depth	430 - 500		480 - 550	380 - 500	380 - 450	530 - 600
Seat height	400 / 450		500	400 / 450	400 / 450	400 / 450
Cross list for common parts						
Backrest type						
Backrest plate	x	x	x	-	x	-
Velcro adjustable	x	x	x	-	-	i
Backrest tubes reclining	x	x	x	x	x	x
Backrest Cover and Cushion Type						
Cover Thin	x	x	x	-	-	-
Cover Lateral	x	x	x	-	-	x
Mistral 2	x	x	x	-	x	-
Laguna	x	x	x	-	x	x
Passad 2	x	x	x	-	-	-
Shoulder High	x	x	x	-	-	-
Backrest options/ pushhandles						
Push bar	x	x	x	-	-	i
Push handles braced	x	x	x	x	x	i
Seat type						
Firm seat plate, depth adj.	x	x	x	x	x	x
Sliding seat	x	x	x	-	-	-
Wheel base seat	-	-	-	x	-	-
H-bracket	-	-	-	x	-	-
Seat Cushions						
Flo-shape	x	x	x	-	x	-
Max seat cushion	-	-	-	-	-	x
Legrests / Hangers						
Central legrest, incl calf pad, footrest and attachment	x	x	x	x	-	-
Hanger angle adjustable	x	x	x	x	x	x
Hanger fixed 80 deg	x	x	x	x	x	-
Hanger fixed 90 deg	x	x	x	x	-	-
Legrest amputee	x	x	x	x	x	-
Legrest attachment width adj.	x	x	x	x	x	x
Footrests						
Footrest one piece	x	x	x	x	x	i
Fixed footplate	x	x	x	x	x	-
Adjustable footplate angle/ depth	x	x	x	x	x	x
Adjustable footplate angle	x	x	x	x	x	-
Armrest Pads						
Long, wide	-	x	-	-	-	x
Padded	x	x	x	x	x	-
Ducktail	x	x	x	x	x	-



Rea™

Azalea® – Individual Parts

Sizes	Azalea®	Azalea® 16"	Azalea® Tall	Azalea® Base	Azalea® Minor	Azalea® Max
Seating options						
Abduction cushion in Dartex	x	x	x	-	x	-
Incontinence cover for seat-cushion	x	x	x	-	x	x
Backrest options						
Neckrest depth adjustable	x	x	x	x	x	x
Headrest with cheek support	x	x	x	x	x	x
Headrest depth / width adjustable	x	x	x	x	x	x
Trunk support multi adjustable with standard arm	x	x	x	-	x	-
Trunk support multi adjustable with extended arm	x	x	x	-	x	-
Trunk support, detachable with standard arm	x	x	x	-	x	-
Trunk support, detachable with extended arm	x	x	x	-	x	-
Trunk support, swing away	x	x	x	-	x	-
Lumbar support cushion	x	x	x	-	x	x
Wedges	x	x	x	-	x	x
Legrest / Footrest options						
Heel straps	x	x	x	x	x	-
Calf pads	x	x	x	x	x	i
Footplate extender	x	x	x	x	x	-
Armrest options						
Hemiplegics	x	x	x	x	x	-
Side support cover	x	x	x	x	x	x
Cushion to side support	x	x	x	x	x	x
Rear Wheel size / type						
16"	-	x	x	x	x	x
22"	x	-	x	x	x	x
24"	x	-	x	x	x	x
Safety						
X-Harness	-	-	-	-	x	-
Accessories						
Table tray Fixed	x	x	x	x	-	x
Table tray Swing-away	x	x	x	x	x	-




ReaTM Azalea[®] – Individual Parts

16" Transit wheels Azalea [®] Max				
Description	Azalea [®] *	Azalea [®] Base	Azalea [®] Max	Azalea [®] Tall
16" solid DB wheels, pair	1536343		1536735	
16" pneumatic DB wheels, pair	1536342		-	1536736

* Not for Azalea Max and Azalea Tall


Seat type (measurements in mm)						
Description	Art.no.					
Seat width	340	390	440	490	550	610
Sliding seat	-	1484576	1484577	1484578	-	-
Seat plate Azalea[®] Base	1508789	1508790	1508791	1508791	1525083 *	1525082 *
H-bracket Azalea[®] Base	1508908				-	-

* Included in the chassie module



H-bracket

Seat cushion Azalea [®] Max (measurements in mm)			
Description	Seat width in mm	Art.no.	
		Grey Dartex	Grey Plush
Max seat cushion	550	1526311	1526310
	610	1526313	1526312



Max seat cushion

Seat cushion Flo-shape seat width + 50 mm and + 100 mm			
Description	Seat width in mm	Art.no.	
		Grey Dartex	Grey Plush
Seat cushion Flo-shape	340 + 50	1509920	1508976
	340 + 100	1509920	1508976
	390 + 50	1508055	1508056
	390 + 100	1508055	1508056
	440 + 50	1508057	1508058
	440 + 100	1508057	1508058
	490 + 50	1508059	1508060
	490 + 100	1508059	1508060

Backrest Velcro adjustable Azalea [®] Max		
Description	Azalea [®] Max	
	Seat width in mm	
Backrest Velcro adjustable Azalea[®] Max	550	610
	1526824	1526823



Backrest Velcro adjustable Azalea Max



Azalea® – Individual Parts

Backrest cushions seat width + 50 mm and + 10 mm

Description	Azalea®	Art.no.	
	Seat width in mm	Grey Dartex	Grey Plush
Mistral 2	340 + 50	1431668	1431669
	340 + 100	1424651	1435404
	390 + 50	1424651	1435404
	390 + 100	1424651	1435404
	4400 + 50	1424652	1435405
	440 + 100	1424652	1435405
	490 + 50	1424653	1435406
	490 + 100	1424653	1435406
Laguna 2 fixed shape	340 + 50	1508972	1510660
	340 + 100	1508972	1510660
	390 + 50	1427119	1435798
	390 + 100	1427119	1435798
	440 + 50	1427120	1444340
	440 + 100	1427120	1444340
	490 + 50	1427121	1444341
	490 + 100	1427121	1444341
Passad 2	390 + 50	1424647	1444329
	390 + 100	1424647	1444329
	440 + 50	1424649	1444330
	440 + 100	1424649	1444330
	490 + 50	1429278	1444331
	490 + 100	1429278	1444331
Shoulder high	390 + 50	1444276	-
	390 + 100	1444276	-
	440 + 50	1444277	-
	440 + 100	1444277	-
	490 + 50	1444278	-
	490 + 100	1444278	-



Pushhandle, right



Pushhandle, left



Push handles braced Azalea Max

Backrest bar / Push handles Azalea® Max (measurements in mm)

Description	Art.no.	
	Azalea® Max	
	550	610
Push handles braced	1524533	1524529
High Push bar	1523947	1523946



High Push bar Azalea Max

Seat and Backrest tilt

Description	Azalea®		Azalea® Max		Azalea® Base	
	22" / 24"	16" / 22" / 24" DB	22" / 24"	16" / 22" / 24" DB	22" / 24"	16" / 22" / 24" DB
Assistant operated seat- and backrest tilt	1536331	1536332	1536335	1536336	1536843	1536839
Electrical operated seat and backrest tilt	1484820	1483996 *	1524790		-	

* Only for 16"



ReaTM **Azalea[®]** – Individual Parts

Legrest attachment

Description	Art.no.					
Seat width in mm	340	390	440	490	550	610
Legrest attachment width adjustable	1508926-0509	1449888-0509	1449889-0509	1449890-0509	-	-
Legrest attachment width adjustable Azalea[®] Tall	-	1507966-0509	1507967-0509	1507968-0509	-	-
Legrest attachment width adjustable Azalea[®] Max	-	-	-	-	1523725 *	1523745 *

* Included in the chassie module

Legrests/hangers


Description	Art.no.	
	Left	Right
Hanger angle adjustable - alu 2008 Azalea[®] Max	1528870	1528871



Hanger angle adjustable - alu 2008
Azalea Max

Footboard converter, cpl

Description	Azalea[®] Max	
	Seat width in mm	Art.no.
Footboard converter, cpl Azalea[®] Max	550	1524731
	550 + 50	1524732
	550 + 100	1524733
	610	1524732
	610 + 50	1524733
	610 + 100	1524734



Footboard converter, cpl

Footrests

Description	Azalea[®] Max		
	Seat width in mm	Short tube	Long tube
Footrest incl. footboard converter, cpl Azalea[®] Max	550	1524977	1524721
	550 + 50	1524978	1524722
	550 + 100	1524979	1524723
	610	1524978	1524722
	610 + 50	1524979	1524723
	610 + 100	1527980	1524724



Footrest incl footboard converter, cpl



Azalea[®] – Individual Parts

Armrest attachment

Description	Art.no.
	Azalea[®] Max
Armrest attachment	1525035

Armrests pads including poles

Description	Art.no.	
	Left	Right
Wide pad, reinforced (Max)	1525067	1525066

Armrest, height adjustable and reinforced, cpl

Description	Art.no.	
	Azalea[®] Max	
	Left	Right
Height adj and reinforced, wide pad, cpl	1525069	1525068



Armrest Azalea Max

Incontinence seat cover seat width + 50 mm and + 100 mm

Description	Seat width in mm	Art.no.
Minor	340	1510636
Small	390	5253078-TR26
Medium	390 + 50	5253079-TR26
	440	
Large	390 + 100	5253080-TR26
	440 + 50	
	490	
X-Large	440 + 100	5253617-TR26
	490 + 50	
	490 + 100	
Azalea[®] Max	550	1526617
	610	1526616

Backrest options

Description	Art.no.	
	Left	Right
Trunksupport with arm extended in height	1484447	1484448
Cane holder	5295015	

Chassis extender

Description	Art.no.
Chassis extender	1484535 *

* Not for **Azalea[®] Minor**



Chassis extender



ReaTM **Azalea[®] – Individual Parts**

Safety

Description	Art.no.	
	Left	Right
Standard brake	1486411	1486410
User manouvered brake Azalea[®] Max	1524719	1524720
Tyre brake for Azalea Minor (22"+24" wheels)	1513040	
Brake kit DB cpl, pair (Azalea Max + Azalea Tall 22"+24")	1536758	1536759
Brake kit DB cpl, pair (Azalea + Azalea Base 22"+24")	1536544	1536545
Brake attachment	1491317	
Anti tip device	1439933	
Reflectors kit (wheel and chassi)	5104965	
Reflectors kit Azalea Transport	5254249	
X-Harness cpl (Minor)	1512914	
Pelvic belt Azalea[®] Max	1528872	



X-Harness

Table tray fixed

Description	Art.no.			
	Seat width in mm	Armrest wide pad	Armrest std pad	Armrest padded
Table tray fixed small	390	-	5206737	-
Table tray fixed medium	390	5206726	-	-
	390 + 50	-	5206726	-
	440	-	5206726	-
Table tray fixed large	390 + 100	-	5206725	-
	440	5206725	-	-
	440 + 50	-	5206725	-
	440	5206725	-	-
	440 + 50	-	5206725	-
	490	-	5206725	-
Table tray fixed X-large	390 + 100	5207643	-	-
	440 + 50	5207643	-	-
	440 + 100	5207643	5207643	-
	490	5207643	-	-
	490 + 50	5207643	5207643	-
	490 + 100	-	5207643	-
Table tray fixed Azalea[®] Max	550 - 710	1528869	-	-




Table tray fixed Azalea Max

Table tray swing-away *

Description	Article number	
	Seat width in mm	Article number
Table tray swing-away (small)	340	1509919
	340 + 50	
Table tray swing-away (large)	340 + 100	1449217
	390 - 590	

* Only in combination with armrest wide pad and padded armrest

Tilt lock (locking of pre-set angle)

Description	Article number
Tilt lock	1527435

Red[™] **Bellis**[™] – Individual Parts

Content

 <i>Red</i> [™] Bellis [™]	45
 This is how <i>Red</i> [™] Bellis [™] was measured	47
 Individual accessories to <i>Red</i> [™] Bellis [™]	48

Rea™ **Bellis™** – Individual Parts

Smart solutions create unique possibilities!


Rea™ **Bellis™** offers many smart solutions to provide a high level of comfort, easy multiple adjustable features and superb driving characteristics. *Rea™* **Bellis™** unique design is focused on your well being.




Rea™ **Bellis™** – Individual Parts

This is how Rea™ Bellis™ has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle




Individual Accessories to Rea™ Bellis™

Seat cushion

Description	Art.no			
	Bellis®		Bellis I50®	
	22" / 24"	16" / 22" / 24" DB	22" / 24"	22" / 24" DB
Assistant operated seat- and backrest tilt	1536429	1536430	1536431	1536432
Electrical operated seat and backrest tilt	1448691			

Seat cushion

Description	Seat width in mm	Art. no
Inlay cushion	410	5255323-TR26
	460	5255324-TR26
	510	5255325-TR26



Inlay cushion

Legrests/ hangers

Description	Seat width	Art. no
Legrest attachment	All widths	1429269-0509




Legrest attachment

ReaTM BellisTM – Individual Parts

Armrest attachment


	Article number	
	Left	Right
Flip up - rear part	5255466-0515	5255465-0515



Flip up - rear part

Armrests lower part

	Article number	
	Left	Right
Height adjustable flip up	1439751	1439750



Height adjustable flip up

Armrests pads including pole

	Article number	
	Left	Right
Long Narrow Nice	1439747	1439746
Long Wide Nice	1439749	1439748



Long Wide Nice



Long Narrow Nice

Rea™ Bellis™ – Individual Parts

16" Transit wheels Bellis

Description	Art.no.
16" solid DB wheels, pair	1536424
16" pneumatic DB wheels, pair	1536423

Castor attachment type

Description	Article number	
	Left	Right
Castor attachment quick release	5254254	5254253



Castor attachment fixed



Castor attachment quick release

Safety

Description	Article number	
	Left	Right
Standard brakes	5255117	5255124
Brake kit 22-24" drumbrake	1536403	1536422
User manouvered tyre brake	1446121	1446122
One arm brake	5255418	5255426
Reflectors kit (wheel and chassi)	5206640	



Standard brakes

Accessories

Description	Article number
Cane holder	5295015



Cane holder

Rea™ **Clematis®** – Individual Parts

Content

 Rea™ Clematis®	51
 This is how Rea™ Clematis® was measured	53
 Individual accessories to Rea™ Clematis®	54

An Improved Seat for Improved Comfort


The foundation of good posture lies in good seating. Maintaining the correct position securely and in comfort is a priority for anyone who spends long periods of time in a wheelchair or cannot move. Targeted support is now available in the new Rea™ **Clematis®** in the form of this new, ergonomically-designed support seat, so you can maintain good posture and enjoy peace of mind.




Rea™ **Clematis®** – Individual Parts

This is how Rea™ Clematis® has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



Rea™ **Clematis® – Individual Parts**

Individual Accessories to Rea™ Clematis®

Backrest Laguna 2

Description	Article number
Backrest Laguna 2 390 cpl	1518681
Backrest Laguna 2 440 cpl	1518682
Backrest Laguna 2 490 cpl	1518683

Pushhandles

Description	Article number	
	Left	Right
Pushhandles for backrest Laguna 2 390-490 mm cpl	1518504	1518503



Pushhandles, right



Pushhandles, left

Trunk support for backrest Laguna 2

Description	Article number	
	Left	Right
Trunk support "Swing-away"	1518708	1518709
Trunk support "Fixed"	1518520	1518519



Trunk support, swing-away



Trunk support, fixed

Legrests/ hangers

Legrest attachment fixed (part of frame)

Armrests pads including pole

Description	Article number	
	Left	Right
Long Narrow padded	1420019	1420018



Armrest

Safety

Description	Article number	
	Left	Right
Standard Brakes	1439932	1439923
Brake attachment (Rail brake exkl.screws)	1431607	
Anti tip device	1439933	
Reflectors kit (wheel and chassi)	5206655	
Reflectors kit including attachm.	1415574	

Rea™ **Focus™** – Individual Parts

Content

 Rea™ Focus™	55
 This is how Rea™ Focus™ was measured	57
 Individual accessories to Rea™ Focus™	58

Rea™ **Focus™** – Individual Parts

Versatility & Comfort


The Rea™ **Focus™** has rigid seat base providing excellent 'upright' postural support, combined with a broad range of body supports. The semi-rigid chassis has a folded brace construction facilitating smooth, effortless motion with negligible wheeling resistance. With state-of-the-art technology and a modular design highly compatible with other Rea™ models, the **Focus™** is very easy to maintain.




Rea™ **Focus™** – Individual Parts

This is how Rea™ Focus™ has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



* As the dimensions are in whole cm steps it might in reality mean slightly more or less adding a couple of mm

Rea™ Focus™ – Individual Parts

Individual Accessories to Rea™ Focus™

Reclined backrest and cushion

Description	Art.no		
	Seat width in mm	Grey Dartex	Grey Plush
Reclined backrest with backrest plate, cushion Laguna 2 and fixed push handles	380 - 405	I508812	I508815
	430 - 455	I508813	I508816
	480 - 505	I508814	I508817
Reclined backrest with backrest plate, cushion Passad 2 and fixed push handles	380 - 405	I518694	I518697
	430 - 455	I518695	I518698
	480 - 505	I518696	I518699
Reclined backrest, gas spring, with backrest plate, cushion Laguna 2 and fixed push handles	380 - 405	I530010	I530045
	430 - 455	I530011	I530046
	480 - 505	I530012	I530241
Reclined backrest, gas spring, with backrest plate, cushion Passad 2 and fixed push handles	380 - 405	I530242	I530266
	430 - 455	I530243	I530267
	480 - 505	I530265	I530268



Backrest plate Laguna 2



Backrest reclined with gas spring

Backrest upholstery cpl

Description	Art.no.		
	Seat width in mm	Jemima	Grey Plush
Anlge Adjustable and Biangular	380	I435649	I435653
	430	I435650	I435654
	480	I435651	I435655



Angle adjustable



Velcro Biangular


Seatcushions

Description	Art.no.			
	Seat depth in mm	Seat width in mm	Grey Dartex	Grey Plush
Mistral 2 seat cushion	375 - 425	380	I434853	I435860
Mistral 2 seat cushion	375 - 425	405 - 430	I435854	I435861
Mistral 2 seat cushion	375 - 425	455 - 480	I435855	I435862
Mistral 2 seat cushion	375 - 425	505	I435859	I435866
Mistral 2 seat cushion	450 - 500	380	I435856	I435863
Mistral 2 seat cushion	450 - 500	405 - 430	I435857	I435864
Mistral 2 seat cushion	450 - 500	455 - 480	I435858	I435865
Mistral 2 seat cushion	450 - 500	505	I4305859	I4305866

Rea™ **Focus™** – Individual Parts

Seating options

Description	Seatwidth in mm	Art.no.
Inlay cushion for seat plate	380	I434768
	430	I434769
	480	I434770
	530	I434771



Inlay cushion

Armrest attachment

Description	Article number	
	Left	Right
Flip up	I431723	



Flip up

Trunk support

Description	Article number	
	Left	Right
Trunk support, small *	5253067-TR23	
Trunk support Swing away **	I527443	I527442

* For reclined backrest without gas spring

** For reclined backrest with gas spring



Trunk support, small

Castor attachment

Description	Article number	
	Left	Right
Castor attachment low, fixed	I434733	
Castor attachment high, height adjustable	I434754	



Attachment low, fixed



Attachment high, fixed



Trunk support, swing away

Rea™ **Focus™ – Individual Parts**

Rearwheel attachment

	Article number
Camber kit	I435807



Camber kit

Safety

Description	Art.no.	
	Left	Right
One arm brake	I449900	I449899
Brake kit 22" + 24" DB	I536437	I536438
Brake kit 22" + 24" DB (narrow)	I539473	I539474
Parking leg (step tube)	5103276	
Anti tip device	I435470	
Anti tip device swing-away *	I508975	I508974
Brake attachment	524024I	5240240
Reflectors kit (wheel and chassi)	5104965	

* not in combination with drum brakes
not in combination with Focus 150 kg



Anti tip device

Accessories

Description	Art.no.
Cane holder	I435754



Cane holder

Rea[™] **Spirea²[®]** – Individual Parts

Content

 <i>Rea</i> [™] Spirea²[®]	61
 This is how <i>Rea</i> [™] Spirea²[®] was measured	63
 Individual accessories to <i>Rea</i> [™] Spirea²[®]	64

Rea™ **Spirea²**® – Individual Parts

Light-handling & high durability


Rea™ **Spirea²**® is an exciting and durable, lightweight folding wheelchair perfect for everyday use. Designed from the ground up using client input & feedback, the Rea™ **Spirea²**® is light handling, minimal adjustment requirements, high technical quality and low maintenance make it the choice for users, carers and engineers alike.




Rea™ **Spirea²®** – Individual Parts

This is how Rea™ Spirea²® has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



* As the dimensions are in whole cm steps it might in reality mean slightly more or less adding a couple of mm

Rea™ Spirea²® – Individual Parts

Backrest type

Seat width in mm	Article number		
	Sling Jemima	Velcro adjustable Jemima	Grey Plush
380	5104583-TR18	5104433-TR18	5104433-TR32
405	5104584-TR18	5104434-TR18	5104434-TR32
430	5104385-TR18	5104435-TR18	5104435-TR32
455	5104386-TR18	5104436-TR18	5104436-TR32
480	1495854	1448653	1508784
505	5104387-TR18	5104437-TR18	5104437-TR32



Sling Jemima



Velcro adjustable

Seat type

Seat width in mm	Article number			
	Depth adjustable		Velcro adjustable	
	Seat depth 400	Seat depth 450	Seat depth 400	Seat depth 450
380	5106337	5106338	5106448	5106454
405	5104005	5106217	5106449	5106455
430	5106339	5106340	5106450	5106456
455	5104006	5106218	5106451	5106457
480	1491249	1491250	1485551	1485552
505	1431755	1434717	5106452	5106458




Depth adjustable seat

Armrest attachment

Description	Article number
Flip up	1435492

Castor attachment

Description	Article number
Castor attachment fixed	1521752



Castor att. fixed

Safety

Description	Art.no	
	Left	Right
One arm brake	1492107	1492000
Brake kit 22" + 24" DB	1536462	1536463
Reflectors kit (wheel and chassi)	5104965	
Reflectors kit transport version	1435760	
Brake attachment	1492299	



Velcro adjustable

Accessories

Description	Article number
Cane holder	5104964

16" Transit wheels Spirea 2T

Description	Art.no.
16" solid DB wheels, pair	1537618
16" pneumatic DB wheels, pair	1537615



Cane holder

Rea[™] **Spirea**^{3®} – *Individual Parts*

Content

 <i>Rea</i> [™] Spirea ^{3®}	65
 This is how <i>Rea</i> [™] Spirea ^{3®} was measured	67
 Individual accessories to <i>Rea</i> [™] Spirea ^{3®}	68

Rea™ **Spirea³®** – Individual Parts

Rea™ **Spirea³®** for transport and frequent use


Rea™ **Spirea³®** is a lightweight-folding wheelchair for every day's use. **Spirea³®** features the basic adjustments for the user with upper balance. **Spirea³®** is designed for light handling in all concerns from the demand of the user and assistant to the prescriber and engineer. From requests from our customers we have made **Spirea³®** easy for refurbishment by high technical quality, low service requirements and a few but necessary and easily operated adjustment options.




Rea™ **Spirea³** – Individual Parts

This is how Rea™ Spirea³® has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



Rea™ Spirea³® – Individual Parts

Backrest type

Seat width in mm	Art.no.	
	Velcro adjustable	
	Jemima	Grey Plush
380	1435649	1435653
405	1478559	1478563
430	1435650	1435654
455	1478560	1478564
480	1435651	1435655
505	1478561	1478565



Velcro adjustable

Seat type


Seat width in mm	Art.no.			
	Depth adjustable		Velcro adjustable	
	Seat depth 40	Seat depth 45	Seat depth 40	Seat depth 45
380	5106337	5106338	5106448	5106454
405	5104005	5106217	5106449	5106455
430	5106339	5106340	5106450	5106456
455	5104006	5106218	5106451	5106457
480	1491249	1491250	1485551	1485552
505	1431755	1434717	5106452	5106458



Depth adjustable seat

Armrest attachment

Description	Art.no.
Flip up	1435492



Attachment fixed

Castor attachment

Description	Art.no.
Castor attachment fixed	1521752



Cane holder

Brake kit 22" + 24" DB, seat angle 0°

Description	Art.no.			
	Active rear wheel plate		Passive rear wheel plate	
	left	right	left	right
Brake kit 22" + 24" DB	1536375	1536376	1536462	1536463

Brake kit 22" + 24" DB, seat angle 3°

Description	Art.no.	
	pneumatic	solid
Brake kit 22" + 24" DB	1536766	1536767

One arm brake

Description	Art.no.	
Seatwidth in mm	Pneumatic	Solid
380	1513044	1513045
405	1513044	1513045
430	1513046	1513095
455	1513046	1513095
480	1513046	1513095
505	1513046	1513095

Accessories

Description	Art.no.
Cane holder	5104964

Safety

Description	Art.no.
Reflectors kit (wheel and chassis)	5104965
Brake attachment	1492299
Anti-tip device swing-away	1513004

Rea™ **Spirea⁴_{NG}** – Individual Parts

Content

 Spirea_{NG}	69
 This is how Spirea_{NG} was measured	71
 Individual accessories to Spirea_{NG}	72

Red™ **Spirea⁴_{NG}** – Individual Parts

Rediscover freedom of mobility


Red™ **Spirea⁴_{NG}** is a dynamic and durable lightweight folding wheelchair ideal for everyday use. By incorporating client input and user feedback into the design, the **Spirea⁴_{NG}** offers light handling, minimal adjustment requirements, high technical quality and low maintenance, making it a perfect choice for users, carers and engineers alike.




Rea™ Spirea⁴_{NG}™ – Individual Parts

This is how Rea™ Spirea⁴® has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



Rea™ **Spirea⁴_{NG}** – Individual Parts

Propulsion type

Description	Art.no.
	Velcro adjustable
	Jemima
Standard wheel hub	1533910
Assistant manouvered drumbrakes	1535066
16" Wheels pneumatic profiled	1537615
16" Wheels solid	1537618

Backrest base, fixed & angle adjustable backrest

Description	Art.no	
	Fixed	Angle adjustable
Backrest base	1534662	1534644



Backrest base

Backrest tube (2pcs / backrest)

Description	Art.no
Backrest tube, fixed backrest	1434742
Backrest tube, adj backrest 350 - 480	1523773
Backrest tube, adj backrest 480 - 550	1533991

Backrest upholstery

Description	Seatwidth in mm						
	355	380	405	430	455	480	505
Adjustable velcro	-	1435649	1478559	1435650	1478560	1435651	1478561
Sling upholstery	1495839	1495840	1495841	1495842	1495843	1495844	1495845

Seat type

Description	Seatdepth in mm	Seatwidth in mm						
		355	380	405	430	455	480	505
Sling seat, depth adjustable	380-400	1534848	1534850	1534851	1534852	1534863	1534870	1534940
	400-450	1534849	1437580	1437581	1437582	1437583	1437584	1437585
	450-500	-	-	-	-	1534869	1534939	1534941
Velcro tension adjustable in depth and shape	380-400	-	1533431	1533432	1533433	1533434	1533435	1533436
	400-450	-	1533438	1533439	1533440	1533441	1533442	1533443
	450-500	-	-	-	-	1533448	1533449	1533450

Rea™ **Spirea⁴_{NG}** – Individual Parts

Complete Amputee legrests

Description	Seatwidth in mm						
	355	380	405	430	455	480	505
Legrest amputee, LH	-	1429192-0509			1429211-0509		
Legrest amputee, RH	-	1429191-0509			1429212-0509		



Fixed footplate



Legrest amputee

Fixed footrest

Description	Seatwidth in mm						
	355	380	405	430	455	480	505
Footrest lower part, fixed footplate, short tube, LH	-	1530620		1530622		1530623	
Footrest lower part, fixed footplate, long tube, LH	-	1530633		1530634		1530635	
Footrest lower part, fixed footplate, short tube, RH	-	1530614		1530616		1530617	
Footrest lower part, fixed footplate, long tube, RH	-	1530630		1530631		1530632	

Armrest type

Armrest type	Art.no	
	left	right
Armrest Comfort, height adjustable, short pad	1536814	1536815
Armrest Comfort, height adjustable, long pad	1536816	1536817



Armrest Comfort, short pad

Castors

Description	Art.no		
	Solid, medium hard	Solid	Pneumatic profiled
Castors 120 mm	-	1450173	-
Castors 140 mm	1515254	-	-
Castors 150 mm	1450622	1450174	1450175
Castors 200 mm (8" x 1.1/4")	-	1431352	1431351
Castors 200 mm (8" x 2")	-	1431759	1431353



Castor with fork

Rea™ **Spirea⁴_{NG}** – Individual Parts

Castor forks

Armrest type	Art.no
Castor attachment, fork L36 (5 inch), pair	I534977
Castor attachment, fork L60 (5,5 inch), pair	I534978
Castor attachment, fork L36 (6 inch & 8 inch), pair	I535043
Castor attachment, fork L60 (6 inch & 8 inch), pair	I535044

Rearwheel attachment


Description	Art.no	
	Std wheel hub	DB
Active rear wheel position	I534828	
Standard rear wheel position	I535065	
Passive rear wheel position		
Dual amputee rear wheel position	I534828	



Rearwheel attachment
Active



Standard rearwheel position



Passive rear wheel position

Rear wheel attachment, one arm drive

Description	Art.no	
	Solid	Pneumatic profiled
Rear wheel attachment, one arm drive	I535065	I534828

One arm drive

Description	Art.no.			
Tyre type	Pneumatic profiled		Solid	
Seatwidth in mm	380 - 430	455 - 505	380 - 430	455 - 505
Dual handrim + Active rear wheel position or Dual amputee rear wheel position	I513269	I513270	I513271	I513272
Dual handrim + Standard rear wheel position or Passive rear wheel position				


Rear wheel attachment, 16" wheels

Description	Art.no	
	left	right
Rear wheel attachment, 16" wheels	I440040	I435464

Rea™ **Spirea⁴_{NG}** – Individual Parts

Transportation kit

Description	Art.no
Transportation kit Spirea4 NG	I534090



Transportation kit

Safety

Description	Art.no	
	left	right
One arm brake	I430203	I430202

Brake kit Spirea 4NG

Description	Art.no	
	left	right
Brake kit 22" + 24" DB	I536733	I536734
Brake user manouvered 16"	I533834	I533908

Anti-tip device

Description	Art.no
Anti-tip device, traditional	I484890 *
Anti-tip device, swing away	I535046 *

* In addition 2 pcs of art.no: 5716276

Pelvic belt

Description	Art.no
Belt buckle	0572-05389

Fitting kit table, Spirea 4

Armrest type	Art.no
Height adjustable Rio, long	I446177
Height adjustable Comfort, long	I446201



Anti-tip device, traditional



Anti-tip device, swing away

Content

 Spin^x	69
 This is how Spin^x was measured	71
 Individual accessories to Spin^x	72

Are you sitting comfortably?


The Invacare® **Spin^x** is designed to fit the user living an active social life and in need of body support and mobility. With an inbuilt postural frame and fully adjustable seat, every user gets individual fitting that supports them in their daily activities. Made with lightweight materials it has a dual folding mechanism making it suitable for storage and transportation.




Invacare® **Spin^x** – Individual Parts

This is how Invacare® Spin^x has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



* As the dimensions are in whole cm steps it might in reality mean slightly more or less adding a couple of mm

Individual Accessories to Invacare® Spin™
Backrest type

Description	Art.no.	
	Black nylon	Jemima
Backrest upholstery sling 380 fixed (380)	I486414	I495840
Backrest upholstery sling 410 fixed (405)	I486415	I495841
Backrest upholstery sling 430 fixed (430)	I486416	I495842
Backrest upholstery sling 460 fixed (455)	I486417	I495843
Backrest upholstery sling 480 fixed (480)	I486418	I495844
Backrest upholstery sling 500 fixed (505)	I486419	I495845
Backrest cover 380 tension adj. (380)	I486477	5106739-TR18
Backrest cover 410 tension adj. (405)	I486478	5106740-TR18
Backrest cover 430 tension adj. (430)	I486479	5106741-TR18
Backrest cover 460 tension adj. (465)	I486480	5106742-TR18
Backrest cover 480 tension adj. (480)	I486481	I435523
Backrest cover 500 tension adj. (505)	I486482	5106743-TR18



Backrest upholstery sling fixed

Backrest options/pushhandles

Description	Art.no.	
	Backrest height 360 - 400	Backrest height 410 - 490
Pushhandles fixed (foldable backrest)	I488860	I489044
Pushhandles fixed	I491207	



Backrest cover

Fixed legrests coloured as the chassie

Description	Art.no.
Hangers 80° Pearl Grey	I493055-0515
Hangers 80° Shining Blue	I493055-0507
Hangers 80° Anthracite Black	I493055-0514
Hangers 80° Red	I493055-0501
Hangers 80° Lime Green (Special)	I493055-0517
Hangers 80° Orange (Special)	I493055-0518
Hangers 90° Pearl Grey	I493056-0515
Hangers 90° Shining Blue	I493056-0507
Hangers 90° Anthracite Black	I493056-0514
Hangers 90° Red	I493056-0501
Hangers 90° Lime Green (Special)	I493056-0517
Hangers 90° Orange (Special)	I493056-0518



Fixed legrest



Sideguard

Armrest to the different wheelchair models

Description	Art.no.	
	left	right
Mudguard high	I486543	I486542
Mudguard low	I484737	I484564
Sideguard	I491238	I491237



Mudguard

Armrest attachment

Description	Art.no.	
	left	right
Flip up - rear part	1485272	1485246

Rear wheel attachment

Description	Art.no.
Camber kit	1483923

Rearwheels 24"

Description	Art.no.
24" Sporty Low profile (24x1")	1483924

Rearwheels 24" with handrims

Description	Handrims					
	Stainless		Foam coated		Black Grip	
Tyres	24" Std	24" DB	24" Std	24" DB	24" Std	24" DB
Low Profile	1531063 *	-	1531533 *	-	1531566 *	-
Pneumatic non marking	1530763 *	1531064 *	1531462 *	1531534 *	1531538 *	1531567 *
Pneumatic profiled	1530764 *	1531080 *	1531508 *	1531535 *	1531563 *	1531568 *
Solid	1530762 *	1531127 *	1531531 *	1531536 *	1531564 *	1531569 *
Puncture proof	1531129 *	1531128 *	1531532 *	1531537 *	1531565 *	1531570 *

Rearwheels 25"

Description	Art.no.
Pneumatic non marking aluminium handrim (special)	1508583

Castor attachment type

Description	Art.no.	
	left	right
Castor attachment fixed	1488858	1488857

Safety

Description	Art.no.	
	left	right
Standard brake	1415485	1415486
Brake kit 22" - 24" DB	1536465	1536466
Brake attachment	1492299	
One arm brake	1492107	1492000
Anti tip device (low 41 - 45,5)	1483885	1483884
Anti tip device (high 48 - 53)	1483999	1483998
Reflectors kit (wheel and chassi)	5104965	
Extended brake arm / pin (pair nr)	1486429	

Accessories

Description	Art.no.
Cane holder	5106390



Anti-tip device, Spin-X

Content

 XLT™	74
 This is how XLT™ was measured	76
 Individual accessories to XLT™	77

Top configurability made Simple!

The new generation *Invacare® XLT* is an active, rigid wheelchair enabling support in day to day life. Designed to be perfectly adjusted to fit the body the *Invacare® XLT* offers many features in terms of support, rigidity and quality of construction. With a frame made entirely from Titanium the *Invacare® XLT* is a light chair helping users to retain mobility as well as gain independence.



Invacare® XLT.



Invacare® XLT. Swing




Invacare® XLT. Max


Invacare® **XLT™ – Individual Parts**

This is how Invacare® XLT. has been measured:


Seat height




Seat depth




Legrest length




Seat width





Seat angle



Backrest height



Backrest angle



* As the dimensions are in whole cm steps it might in reality mean slightly more or less adding a couple of mm

Backrest type Velcro, tension adjustable, fixed height

Description	Art.no.							
	Seatwidth in mm							
Fixed backrest height in mm	330 *	355	380	405	430	455	480	505*
Fixed backrest height 200	I518650	I508716	I508720	I508724	I508728	I508732	I508736	I508740
Fixed backrest height 250	I518651	I508717	I508721	I508725	I508729	I508733	I508737	I508741
Fixed backrest height 300	I518652	I508718	I508722	I508726	I508730	I508734	I508738	I508742
Fixed backrest height 350	I518653	I508719	I508723	I508727	I508731	I508735	I508739	I508743

* XLT Swing

Backrest upholstery, fixed backrest

Description	Art.no.							
	Seatwidth in mm							
Adjustable backrest height in mm	330 *	355	380	405	430	455	480	505*
Fixed backrest height 200	I518668	I512549	I512593	I512606	I512630	I512633	I512654	I512657
Fixed backrest height 250								
Fixed backrest height 300	I518669	I512562	I512604	I512628	I512631	I512634	I512655	I512658
Fixed backrest height 350								

* XLT Swing

Backrest type Velcro, tension adjustable, adjustable height

Description	Art.no.							
	Seatwidth in mm							
Adjustable backrest height in mm	330 *	355	380	405	430	455	480	505*
Adj backrest height 300 - 320	I518654	I508716	I508720	I508724	I508728	I508732	I508736	I508740
Adj backrest height 350 - 390	I518655	I508717	I508721	I508725	I508729	I508733	I508737	I508741
Adj backrest height 400 - 440	I518656	I508718	I508722	I508726	I508730	I508734	I508738	I508742
Adj backrest height 450 - 490	I518657	I508719	I508723	I508727	I508731	I508735	I508739	I508743

* XLT Swing

Backrest upholstery, adjustable backrest

Description	Art.no.							
	Seatwidth in mm							
Adjustable backrest height in mm	330 *	355	380	405	430	455	480	505*
Adj backrest height 300 - 390	I518669	I512562	I512604	I512628	I512631	I512634	I512655	I512658
Adj backrest height 400 - 490	I518670	I512563	I512605	I512629	I512632	I512653	I512656	I512659

* XLT Swing

Narrow back *

Description	Art.no.						
	Seatwidth in mm						
Narrow back, velcro tension adj, fixed height, infinitely variable AA, foldable & click in backrest	355	380	405	430	455	480	505**
Fixed backrest height 350 mm	I516333	I516334	I516335	I516336	I516337	I516338	I516339

* Not for XLT Max

** XLT Swing


Backrest upholstery, narrow back

Description	Art.no.						
	Seatwidth in mm						
Upholstery narrow back	355	380	405	430	455	480	505*
Fixed backrest height 350 mm	I517360	I517361	I517362	I517363	I517364	I517365	I517366

* XLT Swing

Backrest type Velcro tension adjustable, adjustable height, XLT Max

Description	Art.no.		
	Seatwidth in mm		
	505	555	605
Adjustable backrest height 400 - 440 mm	I528585	I528645	I528687
Adjustable backrest height 450 - 490 mm	I528642	I528684	I528688



Narrow back

Backrest upholstery adjustable backrest, XLT Max

Description	Art.no.		
	Seatwidth in mm		
	505	555	605
Adjustable backrest height 400 - 440 cm	I535490	I535491	I535492
Adjustable backrest height 450 - 490 cm	I535490	I535491	I535492



Adjustable backrest XLT Max

Pushhandles including neckrest or headrest

Description	Art.no.							
Seat type seat width in mm	330 *	355	380	405	430	455	480	505 *
Height adjustable pushhandles external long grip incl neckrest	I530506	I530510	I512788	I512789	I512790	I512791	I512792	I512793
Height adjustable pushhandles external long grip incl headrest	I530511	I530512	I512782	I512783	I512784	I512785	I512786	I512787
Height adjustable pushhandles external long grip incl headrest with cheek support	I530535	I530536	I512743	I512744	I512745	I512746	I512747	I512748

* XLT Swing

Pushhandles

Description	Art.no.
Fixed	1511455
Foldable	1511472
Height adjustable external short grip	1512795
Height adjustable external long grip	1512797
End part backrest when no pushhandles	1511595



Fixed pushhandles



Foldable pushhandles



Height adj. pushhandles

Seat type tension adjustable, depth adjustable 40 mm

Description	Art.no.							
Seat type seat width in mm	330	355	380	405	430	455	480	505
Seatdepth 400 mm	1518674 *	1512243	1512244	1512301	1512389	1512393	1512496	1512497 *
Seatdepth 450 mm	1518675 *	1512498	1512499	1512500	1512501	1512502	1512503	1512504 *
Seatdepth 500 mm	-	1512505	1512506	1512507	1512529	1512530	1512547	1512548 *

* XLT Swing

Seat type tension adjustable, depth adjustable 40 mm

Description	Art.no		
	Tension adjustable, depth adjustable 4 cm XLT Max		
Seat type seat width in mm	SW 505	SW 555	SW 605
Seat depth 500 mm	1535453	1529661	1529754
Seat depth 550 mm	1529474	1529672	1529673
Seat depth 600 mm	1529660	1529673	1529756

Seat type tension adjustable, depth adjustable 40 mm

Description	Art.no		
	Flo-tech Contour seat cushion		
Seat type seat width in mm	SW 505	SW 555	SW 605
Seat depth 500 mm	1534690	1535692	1534694
Seat depth 550 mm	1534691	1534693	1534695
Seat depth 600 mm	-	-	-



Flo-tech Contour seat cushion

* XLT Swing

Fixed legrests coloured as the chassie (XLT Swing)

Description	Art.no.
Hangers 80° Pearl Grey	I493055-0515
Hangers 80° Shining Blue	I493055-0507
Hangers 80° Anthracite Black	I493055-0514
Hangers 80° Red	I493055-0501
Hangers 80° Orange (Special)	I493055-0518
Hangers 90° Pearl Grey	I493056-0515
Hangers 90° Shining Blue	I493056-0507
Hangers 90° Anthracite Black	I493056-0514
Hangers 90° Red	I493056-0501
Hangers 90° Orange (Special)	I493056-0518



Fixed legrest



Amputee legrest

Amputee legrest (XLT Swing)

Description	Art.no.	
	left	right
Amputee legrest	I512727	I512728

One piece footrest to fit seat width in mm

Description	Art.no.					
	Seat width in mm					
	355	380	405	430	455	480
Adjustable one piece footrest in depth and angle to fit Active and Dynamic frame	-	I512689	I512690		I512691	

Footrest High to fit seat width in mm (XLT 75/90)

Description	Art.no.					
	Seat width in mm					
	355	380	405	430	455	480
Footrests High to fit Active and Dynamic frame	I512692	I512693	I512694	I512758	I512759	I512760

Footrest High to fit seat width in mm (XLT Max)

Description	Art.no.		
	Seat width in mm		
	505	550	605
Footrests High to fit XLT Max	I536787	I536788	I536789

Footbow rigid to fit seat width in mm (XLT 75/90)

Seat width in mm	Art.no.	
	low	high
355	I512671	I512677
380	I512672	I512678
405	I512673	I512679
430	I512674	I512680
455	I512675	I512681
480	I512676	I512682



Footbow rigid

Cover for footbow rigid to fit seat width in mm

Seat width in mm	Art.no.
355	1512683
380	1512684
405	1512685
430	1512686
455	1512687
480	1512688

Footbow rigid, XLT Max, to fit seat width in mm

Seat width in mm	Art.no.
	Low
505	1528689
555	1528690
605	1528689



Footbow rigid, XLT Max

Cover for footbow rigid, XLT Max, to fit seat width in mm

Seat width in mm	Art.no.
505	1532825
555	1532826
605	1532827

Armrest attachment


Description	Art.no.	
	left	right
Flip-up Rear (pair)	1512697	
Rio XLT (pair)	1512696	
XLT Max	1533950	

Armrest to the different wheelchair models

Description	Art.no.	
	left	right
Mudguard high	1512845	1512844
Mudguard low	1512847	1512846
Sideguard	1491238	1491237



Mudguard



Sideguard

Castor attachment type

Description	Art.no.	
	Active/Dynamic frame	Swing frame / XLT Max
Castor attachment	1512844	1512845

Article number of castors in combination with forks

Castor (Approximate in inch")	Forks	Art.no. fixed castor + fork
75x24 soft (3")	75"	1512762
	110"	1512763
	150"	1512764
100x27 medium hard (4")	75"	1512765
	110"	1512766
	150"	1512767
125x27 soft (5")	110"	1512768
	150"	1512769
140x37 soft (5,5")	110"	1512770
	150"	1512771
150x27 medium hard (6")	110"	1512772
	150"	1512773
150x30 soft (6")	110"	1512774
	150"	1512775
150x32 pneumatic (6")	110"	1512776
	150"	1512777
180x40 pneumatic (7")	110"	1512780
	150"	1512781
180x40 soft (7")	110"	1512778
	150"	1512779




Castor cpl

Rearwheel attachment

Description	Art.no.
Rearwheel fixation std 0 degrees	1511869
Camber kit 3 degrees	1511930
Rear wheel extension (behind backrest frame)	1512236
Kit, DB, cpl, pair, incl, attachm.	1512718



Rearwheel fixation std XLT Max



Camber 3 deg



Rearwheel extension

Rearwheels 24"

Description	Art.no.
Sporty low profile, aluminium handrim (24 x 1")	1483924
24" Spinergy, low profile (24x1"), black marathon evolution, titanium handrims	1497787
Puma low profile, aluminium handrim (24 x 1")	5750796
Sporty low profile Black Marathon plus evolution, aluminium handrim(24 x 1")	1512822
Black Marathon plus evolution, aluminium handrim (24 x 1")	1512709
Black Marathon plus evolution, stainless handrim(24 x 1")	1512710
Black Marathon plus evolution, plastic coated handrim (24 x 1")	1512711
Black Marathon plus evolution, foam coated handrim(24 x 1")	1512712
Black Marathon plus evolution, black grip handrim (24 x 1")	1512713
Black Marathon plus evolution, no handrim, no clips (24 x 1")	1512714

Rearwheels 24" with std handrims

Description	24"		
	Handrims		
Tyres	Stainless	Foam coated	Black Grip
Low Profile	1531063	1531533	1531566
Pneumatic non marking	1530763	1531462	1531538
Pneumatic profiled	1530764	1531508	1531563
Solid	1530762	1531531	1531564
Puncture proof	1531129	1531532	1531565



Rearwheels

Rearwheels 24"DB with std handrims

Description	24"		
	Handrims		
Tyres	Stainless	Foam coated	Black Grip
Low Profile	-	-	-
Pneumatic non marking	1531064	1531534	1531567
Pneumatic profiled	1531080	1531535	1531568
Solid	1531127	1531536	1531569
Puncture proof	1531128	1531537	1531570



Handrims

Rearwheels 25" / 26"

Description	Art.no.
Pneumatic non marking, aluminium handrim 25"	I508583 *
Pneumatic non marking, aluminium handrim 26"	I508584 *

* Not for XLT Max

Rearwheels 24" XLT Max

Description	Art.no.				
	Aluminium	Stainless	Plastic coated	Foam coated	Black grip
24" Tyres black marathon evolution	I532143	I532144	I532183	I532184	I532185

Trunk support *

Description	Art.no.	
	left	right
Trunk support	I512867	I512866

* Not for XLT Max

One arm drive **

Description	Art.no.						
	Seat width in mm						
	355	380	405	430	455	480	505*
Profiled tyres	I512860			I512861			I512862
Solid tyres	I512863			I512864			I512865

* XLT Swing ** Not for XLT Max

Safety

Description	Art.no.	
	left	right
Standard brake	I512662	
Brake kit 22" + 24" DB, XLT Std	I536468	I536469
Brake kit 22" + 24" DB, XLT Max	I539477	I539478
Active brake **	I512665	
Performance brake *	I512668	
One arm brake **	I512669	I512670
Anti tip device	I512715	I512716
Anti tip device extended **	I518488	

* Not for XLT Swing and XLT Max ** Not for XLT Max



Standard brake



Active brake



Performance brake

INVACARE® EUROPE:

Belgium & Luxemburg: Invacare nv • Autobaan 22 • 8210 Loppem • www.invacare.be

Danmark: Invacare A/S • Sdr. Ringvej 37 • 2605 Brøndby • www.invacare.dk

Deutschland & Eastern Europe: Invacare Aquatec GmbH • Alemannenstraße 10 • 88316 Isny • www.invacare-aquatec.de

España: Invacare SA • c/ Areny • s/n Polígono Industrial de Celrà • 17460 Celrà (Girona) • www.invacare.es

European Distributor Organisation: Invacare • Kleiststraße 49 • 32457 Porta Westfalica • www.invacare.eu.com

France & Greece: Invacare Poirier SAS • Route de St Roch • 37230 Fondettes • www.invacare.fr

Ireland: Invacare Ireland • Unit 5 Seatown Business Campus • Seatown Rd • Swords • Dublin • www.invacare.ie

Italia: Invacare Mecc San s.r.l. • Via dei Pini 62 • 36016 Thiene (VI) • www.invacare.it

Nederland: Invacare BV • Celsiusstraat 46 • 6716 BZ Ede • www.invacare.nl

Norge & Iceland: Invacare AS • Grensesvingen 9 • Postboks 6230 Etterstad • 0603 Oslo • www.invacare.no

Österreich: Invacare Austria GmbH • Herzog Odilostrasse 101 • 5310 Mondsee • www.invacare.at

Portugal: Invacare Lda • Rua Estrada Velha, 949 • 4465-784 Leça do Balio • www.invacare.pt

Sverige & Suomi: Invacare AB • Fagerstagatan 9 • P.O. Box 66 • 16391 Spånga • www.invacare.se

Switzerland: Invacare AG • Benkenstrasse 260 • 4108 Witterswil • www.invacare.ch

United Kingdom: Invacare Ltd • Pencoed Technology park • Pencoed • Bridgend CF35 5HZ • www.invacare.co.uk



Yes, you can.®

Invacare® International • Route de Cité-Ouest 2 - 1196 Gland - Switzerland
www.invacare.eu.com