









ROLSTOEL HOOGTE INVACARE® SPIN^X

24" Rear wheel

															
H mm	x°	x"	pos.	pos.	pos.	pos.	pos.	pos.	mm	SD40	SD45	SD50			
410	0	24	-	4	2	-	-	-	100	1,2,3,4	1,2,3,4	1,2,3,4			
410	3	24	-	5	2	-	-	-	100	1,2,3,4	1,2,3,4	1,2,3,4			
410	6	-	-	-	-	-	-	-	-	-	-	-			
430	0	24	-	3	1	-	-	-	100	1,2,3,4	1,2,3,4	1,2,3,4			
430	3	24	-	4	1	-	-	-	100	-	1,2,3,4	1,2,3,4			
430	6	24	-	5	1	-	-	-	100	1,2,3,4	1,2,3,4	1,2,3,4			
430	0	24	-	2	-	4	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
430	3	24	-	4	-	4	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
430	6	24	-	5	-	4	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
455	0	-	-	-	-	-	-	-	-	-	-	-			
455	3	24	-	3	-	3	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
455	6	24	-	4	-	3	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
455	0	24	-	1	-	3	-	-	140	1,2,3,4	1,2,3,4	1,2,3,4			
455	3	24	-	2	-	3	-	-	140	1,2,3,4	1,2,3,4	1,2,3,4			
455	6	24	-	3	-	3	-	-	140	1,2,3,4	1,2,3,4	1,2,3,4			
480	0	-	-	-	-	-	-	-	-	-	-	-			
480	3	24	-	1	-	1	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
480	6	24	-	2	-	1	-	-	120	1,2,3,4	1,2,3,4	1,2,3,4			
480	0	24	5	-	-	-	-	4	120	1,2,3,4	1,2,3,4	1,2,3,4			
480	3	-	-	-	-	-	-	-	-	-	-	-			
480	6	24	-	2	-	-	-	4	120	1,2,3,4	1,2,3,4	1,2,3,4			
480	0	24	4	-	-	-	-	4	140	1,2,3	1,2,3,4	1,2,3,4			
480	3	24	5	-	-	-	-	4	140	1,2,3	1,2,3,4	1,2,3,4			
480	6	-	-	-	-	-	-	-	-	-	-	-			
480	0	24	4	-	-	1	-	-	150	1,2,3,4	1,2,3,4	1,2,3,4			
480	3	24	5	-	-	1	-	-	150	1,2,3,4	1,2,3,4	1,2,3,4			
480	6	-	-	-	-	-	-	-	-	-	-	-			
505	0	24	3	-	-	-	-	1	100	1,2,3,4	1,2,3,4	1,2,3,4			
505	3	24	4	-	-	-	-	1	100	1,2,3,4	1,2,3,4	1,2,3,4			
505	6	24	5	-	-	-	-	1	100	1,2,3,4	1,2,3,4	1,2,3,4			
505	0	24	3	-	-	-	-	2	120	1,2,3,4	1,2,3,4	1,2,3,4			
505	3	24	4	-	-	-	-	2	120	1,2,3,4	1,2,3,4	1,2,3,4			
505	6	24	5	-	-	-	-	2	120	1,2,3,4	1,2,3,4	1,2,3,4			
505	0	24	2	-	-	-	-	2	140	1,2,3	1,2,3,4	1,2,3,4			
505	3	24	3	-	-	-	-	2	140	1,2,3	1,2,3,4	1,2,3,4			
505	6	-	-	-	-	-	-	-	-	-	-	-			
505	0	24	3	-	-	-	-	3	150	1,2,3	1,2,3,4	1,2,3,4			
505	3	24	4	-	-	-	-	3	150	1,2,3	1,2,3,4	1,2,3,4			
505	6	24	5	-	-	-	-	3	150	1,2,3	1,2,3,4	1,2,3,4			
505	0	24	3	-	-	1	-	-	180	-	1,2	1,2,3,4			
505	3	24	4	-	-	1	-	-	180	-	1,2	1,2,3,4			
505	6	-	-	-	-	-	-	-	-	-	-	-			
505	0	24	3	-	-	-	-	4	180	1	1,2,3	1,2,3,4			
505	3	24	4	-	-	-	-	4	180	1	1,2,3	1,2,3,4			
505	6	24	5	-	-	-	-	4	180	1	1,2,3	1,2,3,4			
505	0	24	2	-	-	1	-	-	200	-	1,2	1,2,3,4			
505	3	24	3	-	-	1	-	-	200	-	1,2	1,2,3,4			
505	6	-	-	-	-	-	-	-	-	-	-	-			
530	0	-	-	-	-	-	-	-	-	-	-	-			
530	3	24	2	-	-	-	-	1	140	1,2,3	1,2,3,4	1,2,3,4			
530	6	24	4	-	-	-	-	1	140	1,2,3	1,2,3,4	1,2,3,4			
530	0	24	1	-	-	-	-	1	150	1,2,3	1,2,3,4	1,2,3,4			
530	3	24	2	-	-	-	-	1	150	1,2,3	1,2,3,4	1,2,3,4			
530	6	24	3	-	-	-	-	1	150	1,2,3	1,2,3,4	1,2,3,4			
530	0	24	1	-	-	-	-	2	180	1	1,2,3	1,2,3,4			
530	3	24	2	-	-	-	-	2	180	1	1,2,3	1,2,3,4			
530	6	24	3	-	-	-	-	2	180	1	1,2,3	1,2,3,4			
530	0	24	1	-	-	-	-	3	200	-	1,2	1,2,3,4			
530	3	24	2	-	-	-	-	3	200	-	1,2	1,2,3,4			
530	6	24	3	-	-	-	-	3	200	-	1,2	1,2,3,4			

- = Niet mogelijk