

Spare Parts Catalogue

Version: 2017/04/28a NL



Invacare® **Spin^x**



Yes, you can.®

General Information

For any order of spare parts the following information is required:

- Product type*
- Product name & model*
- Serial Number*
- Part-Number*
- Part-Description*

Exploded view drawings and code numbers are not contractual and subject to modification, in relation to Invacare® continuous product improvement policy to meet end users needs.

The reproduction or use of this work, or of any substantial part of this work, in any form or by any electronic, mechanical or other means, including photocopying, scanning and recording, and in any information storage and retrieval system, including electronic transmission, is forbidden without the prior written permission of Invacare® International Sarl.



Yes, you can.®

TABLE OF CONTENTS

(Chassis)

(Chassis, SpinX)	5
(Rear wheel attachment)	6
(Rear wheel attachment, passive)	7
(Castor attachment)	8

(Seat)

(Seat upholstery)	9
(Seat upholstery, adj.)	10
(Seat upholstery, padded)	11
(Seat cushion)	12

(Backrest)

(Backrest tube, fixed)	14
(Backrest tube, adj.)	15
(Backrest tube with external push handles)	16
(Backrest upholstery, fixed)	17
(Backrest upholstery, adj.)	18
(Backrest MISTRAL2, tension adj. (->2012.04))	19

(Legrests)

(Legrest 80° Alu2003, Fixed)	20
(Legrest 90° Alu2003, Fixed)	21
(Legrest AA)	22
(Amputee legrest, Alu2003)	23
(Amputee legrest)	24

(Footrests)

(Footrest, fixed, Alu2003)	25
(Footrest, AA, "Alu2003")	26
(Footrest 90°, Alu2003)	28
(Footrest one piece 80°)	29
(Footrest one piece 90°)	30
(Calf strap, Alu2003 legrest)	31
(Heel strap)	32

(Siderests)

(Armrest, Flip-up)	33
(Armrest Rio, Height Adj.)	34
(Armrest attachment)	35
(Sideguard)	36
(Mudguard)	37

(Front Wheels)

(Castor wheels)	38
-----------------------	----

(Rear Wheels)



Yes, you can.®

TABLE OF CONTENTS

(Wheel 24"-22-20)	39
(Kit Handrim 24" Stainless 4 Fixations)	41
(Wheel 26-25")	42
(Wheel 24" SPORTY)	43
(Wheel 24"-22" DB)	44
(Kit Handrim 24" Stainless 4 Fixations)	41

(Brakes)


(Brake)	46
(Brake extender)	47
(One arm brake)	48
(Drumbrake, assistant manoeuvred SpinX)	49
(Drumbrake assistant manoeuvred (=> 2010.11))	50

(Safety Options)

(Pelvic belt, SpinX)	51
(Anti-tip device)	52
(Reflector)	53


(Others Options & Accessories)

(Backrest stay)	54
(Step tube)	55
(Crutch holder)	56
(Pump)	57
(Table tray)	58




1	1507973	(STOP BLOCK, HANGER CPL., PAIR , SPINX)	1
2	1507974	SpinX frame dop tbv. zitting set van 4	1
3	1507975	SpinX zitframe zadel pp	1
4	1507976	(COVER, CHASSI CPL., PAIR SPINX)	1
5	1507977	SpinX voorwiellagerhuis bev.plaat pp	1
6	1508141	SpinX kruisfr.vergrendeling zb 38-43 ps	1
7	1508142	SpinX kruisfr.vergrendeling zb 45,5 t/m 50,5 ps	1
8	SP1584734	(Screw kit, cross brace, SpinX)	1

(Rear wheel attachment)





1	1488853	SpinX achterw.bev.0 gr. camber cpl	1
2	1483923	SpinX achterw.bev.3 gr. camber cpl	1
3	5750166	Rea achterwiel as QR incl.asbehuizing cpl	1
4	5750165	Rea achterwiel as QR ps	1


(Rear wheel attachment, passive)



SPS1491236

1	1529770	(Kit, passive rear wheel attachment, pair)	1
2	1491236	SpinX asplaat verlengd tbv.amp.steun ps	1
3	1446117	(BRAKE BODY LE, USER BRAKE, 16 WHEEL)	1
4	1446116	(BRAKE BODY RI, USER BRAKE, 16 WHEEL)	1
5	SP5221280	(Handle grip for parking brake)	1
6	5750165	Rea achterwiel as QR ps	1
7	5750166	Rea achterwiel as QR incl.asbehuizing cpl	1

					
cm	40	45	50		
35.5	1	2			
38	3	4	5		
40.5	6	7	8		
43	9	10	11		
45.5	12	13	14		
48	15	16	17		
50.5	18	19	20		




SP2519006

1	5106214	Activ zitting verst.35,5 x40 cpl AD/SP	1
2	5106216	\$ctiv zitting verst.36x4 5 cpl Ad/Sp	1
3	5106337	Activ zitting verst.38x4 0 cpl Ad/Sp	1
4	5106338	\$ctiv zitting verst.38x4 5 cpl Ad/Sp	1
5	1491251	(Upholstery Adj. Cpl., Seat 38/50)	1
6	5104005	Activ zitting verst.41x4 0 cpl Ad/Sp	1
7	5106217	Activ zitting verst.41x4 5 cpl Ad/Sp	1
8	1491252	(Upholstery Adj. Cpl., Seat 40.5/50)	1
9	5106339	Activ zitting verst.43x4 0 cpl Ad/Sp	1
10	5106340	Activ zitting verst.43x4 5 cpl Ad/Sp	1
11	1491253	Activ zitting verst.43x 50 cpl Spx/Sp	1
12	5104006	Activ zitting verst 45,5 x 40 cpl Ad/Sp	1
13	5106218	Activ zitting verst.46x4 5 cpl Ad/Sp	1
14	1491254	Activ zitting verst.46x 50 cpl	1
15	1491249	Activ zitting verst.48x 40 cpl	1
16	1491250	Activ zitting verst.48x 45 cpl spx/Sp	1
17	1491284	Activ zitting verst.48x 50 cpl	1
18	1440017	\$ctiv zitting verst.50,5 x40 cpl	1
19	1434717	Activ zitting verst.50,5 x45 Ad/Sp	1
20	1488065	Activ zitting verst.50,5 x50 cpl Sp3	1

(Seat upholstery, adj.)



cm	Wheelchair		
	40	45	50
38	1	2	3
40.5	4	5	6
43	7	8	9
45.5	10	11	12
48	13	14	15
50.5	16	17	18




SPS106447

1	5106448	\$ctiv zitting contour d.v.38x40 Ad/Sp	1
2	5106454	Activ zitting contour d.v.38x45 Ad/Sp	1
3	1485558	(Seat 38/50, tension adj.)	1
4	5106449	\$ctiv zitting contour d.v.41x40 Ad/Sp	1
5	5106455	Activ zitting contour d.v.41x45 Ad/Sp	1
6	1485559	SpinX zitting Contour DV 41x50	1
7	5106450	\$ctiv zitting contour d.v.43x40 Ad/Sp	1
8	5106456	Activ zitting Contour d.v.43x45 SpinX/Ad/Sp	1
9	1485560	Activ zitting 43x50 nasp +d.i.spinX/Spirea	1
10	5106451	Activ zitting contour d.v.46x40 Ad/Sp	1
11	5106457	Activ zitting contour d.v.46x45 Ad/Sp	1
12	1485561	Activ zitting 46x50 nasp +d.i.spinX/Spirea	1
13	1485551	(Seat 48/40, tension adj.)	1
14	1485552	SpinX zitting contour DV 48x45	1
15	1485557	SpinX zitting contour DV 48x50	1
16	1435834	Activ zitting verst.50,5 x40	1
17	1435835	Activ zitting verst.50,5 x45	1
18	1485562	Activ zitting contour d.v.50x50 Ad/Sp	1

(Seat upholstery, padded)

			
cm	40	45	50
38	1	2	3
40.5	4	5	6
43	7	8	9
45.5	10	11	12
48	13	14	15
50.5	16	17	18




SPINX®

1	1486496	(SEAT 38/40 PADDED,CPL TR33)	1
2	1486497	(SEAT 38/45 PADDED,CPL TR33)	1
3	1488068	(SEAT 38/50 PADDED,CPL TR33)	1
4	1486498	(SEAT 40.5/40 PADDED,CPL TR33)	1
5	1486499	(SEAT 40.5/45 PADDED,CPL TR33)	1
6	1488069	(SEAT 40.5/50 PADDED,CPL TR33)	1
7	1486500	(SEAT 43/40 PADDED,CPL TR33)	1
8	1486501	(SEAT 43/45 PADDED,CPL TR33)	1
9	1488070	SpinX zitting padded 45x50 TR33	1
10	1486502	(SEAT 45.5/40 PADDED,CPL TR33)	1
11	1486503	(SEAT 45.5/45 PADDED,CPL TR33)	1
12	1488071	(SEAT 45.5/50 PADDED,CPL TR33)	1
13	1488066	(SEAT 48/40 PADDED,CPL TR33)	1
14	1488067	SpinX zitting 48x45 cpl gepolsterd tr33	1
15	1478606	(SEAT 48/50 PADDED,CPL TR33)	1
16	1486504	(SEAT 50.5/40 PADDED,CPL TR33)	1
17	1487441	Rea zitting gepolsterd 50,5x45 TR33 SpX/Sp2 cpl	1
18	1478607	(SEAT 50.5/50 PADDED,CPL TR33)	1

(Seat cushion)

	SW	SD	TR18	TR32	TR35
Spirea ² Focus	380	400	6	38	25
		450	7	39	26
	405	400	9	40	27
		450	10	41	28
	430	400	12	42	29
		450	13	43	30
	455	400	15	44	31
		450	16	45	32
	480	400	18	46	33
		450	19	47	34
505	400	21	48	35	
	450	22	49	36	
530	450	24	50	37	
555	450	24	50	37	

	SW	SD	400	450	500
Spin ^x	330	1	2	-	-
XLT	355	3	4	5	-
Spirea ⁴ NG	380	6	7	8	-
	405	9	10	11	-
	430	12	13	14	-
	455	15	16	17	-
	480	18	19	20	-
TR18	505	21	22	23	-



	Färg	Colour	Farbe	Couleur
TR18	jemima	jemima	Jemima	Jemima
TR32	plysch grå PRIMA	plush grey PRIMA	Plüsch grau PRIMA	velours gris PRIMA
TR35	plysch svart	plush black	Plüsch schwarz	velours noir


SPS510179

1	5226506-TR18	\$ea zitkussen 34x40x5cm zwart	1
2	5226507-TR18	\$ea zitkussen 34x45x5cm zwart	1
3	5104255-TR18	Rea zitkussen 36x40x5cm zwart	1
4	5104257-TR18	Rea zitkussen 36x45x5cm zwart	1
5	1512809	\$CTIV ZITKUSSEN 40,5x50 ZWART PS	1
6	5104317-TR18	Rea zitkussen 38x40x5cm zwart	1
7	5104319-TR18	Rea zitkussen 38x45x5cm zwart	1
8	1495830	Rea zitkussen 38x50x5cm zwart	1
9	5103179-TR18	Rea zitkussen 41x40x5cm zwart	1
10	5104259-TR18	Rea zitkussen 41x45x5cm zwart	1
11	1495831	Rea zitkussen 40x50x5cm zwart	1
12	5104262-TR18	Rea zitkussen 43x40x5cm zwart	1
13	5104264-TR18	Rea zitkussen 43x45x5cm zwart	1
14	1495832	Rea zitkussen 43x50x5cm zwart	1
15	5103180-TR18	Rea zitkussen 46x40x5cm zwart	1
16	5104266-TR18	Rea zitkussen 46x45x5cm zwart	1
17	1495833	Rea zitkussen 46x50x5cm zwart	1
18	5226509-TR18	Rea zitkussen 48x40x5cm zwart	1
19	5226510-TR18	Rea zitkussen 48x45x5cm zwart	1
20	1495834	Rea zitkussen 48x50x5cm zwart	1
21	5104269-TR18	Rea zitkussen 50x40x5cm zwart	1
22	5104271-TR18	Rea zitkussen 50x45x5cm zwart	1
23	1495835	Rea zitkussen 50x50x5cm zwart	1
24	1435429	Rea zitkussen 53x45x5cm zwart	1
25	5104317-TR35	(Seat Cushion 380x400, Plush Black, Standard)	1
26	5104319-TR35	(Seat Cushion 380x450, Plush Black, Standard)	1
27	5103179-TR35	(Seat Cushion 410x400, Plush Black, Standard)	1
28	5104259-TR35	(Seat Cushion 410x450, Plush Black, Standard)	1
29	5104262-TR35	(Seat Cushion 430x400, Plush Black, Standard)	1
30	5104264-TR35	(Seat Cushion 430x450, Plush Black, Standard)	1
31	5103180-TR35	(Seat Cushion 460x400, Plush Black, Standard)	1

(Seat cushion)

32	5104266-TR35	(Seat Cushion 460x450, Plush Black, Standard)	1
33	5226509-TR35	(Seat Cushion 480x400, Plush Black, Standard)	1
34	5226510-TR35	(Seat Cushion 480x450, Plush Black, Standard)	1
35	5104269-TR35	(Seat Cushion 500x400, Plush Black, Standard)	1
36	5104271-TR35	(Seat Cushion 500x450, Plush Black, Standard)	1
37	1435445-TR35	(Seat Cushion 530x450, Plush Black, Standard)	1
38	5104317-TR32	\$ea zitkussen 38x40x5cm grijs velours	1
39	5104319-TR32	\$ea zitkussen 38x45x5cm grijs velours	1
40	5103179-TR32	\$ea zitkussen 41x40x5cm grijs velours	1
41	5104259-TR32	\$ea zitkussen 41x45x5cm grijs velours	1
42	5104262-TR32	\$ea zitkussen 43x40x5cm grijs velours	1
43	5104264-TR32	\$CTIV ZITKUSSEN 43x45+ HOES GRIJS PRIMA AD/SP	1
44	5103180-TR32	\$ea zitkussen 46x40x5cm grijs velours	1
45	5104266-TR32	\$ea zitkussen 46x45x5cm grijs velours	1
46	5226509-TR32	\$ea zitkussen 48x40x5cm grijs velours	1
47	5226510-TR32	\$ea zitkussen 48x45x5cm grijs velours	1
49	5104271-TR32	\$ea zitkussen 50x45x5cm grijs velours	1
50	1435445	\$ea zitkussen 53x45x5cm grijs velours	1


(Backrest tube, fixed)



SPS1491207


1	1507970	RUGBUIZEN VAST SPINX PP	1
2	1491207	SpinX rugb.vast lang ps cpl	1
3	SP1429301	Rea duwh.opb.kort ex.bev ps	1
4	1543736	(FIXATION, BACKREST UPHOLSTERY)	1

Invacare® SpinX | (Backrest)
(Backrest tube, adj.)



1	1507981	(BACKREST TUBE ADJ.LOW PAIR)	1
2	1507984	SpinX rugbuis ihv RH41-49 pp	1
3	1507982	(BACKREST TUBE ADJ.LOW LEFT)	1
4	1507983	(BACKREST TUBE ADJ.LOW RIGHT)	1
5	1507985	(BACKREST TUBE ADJ. HIGH LEFT)	1
6	1507986	(BACKREST TUBE ADJ. HIGH RIGHT)	1
7	1488860	SpinX rugbuis kort incl. duwh.ps	1
8	1489044	SpinX rugbuis met duwhandvat vast rh 41-49 ps	1
9	1485415	SpinX rugb.bev.opkl.tbv. rh 36-40 pp	1
10	1485416	SpinX rugb.bev.opkl.tbv. rh 41-49 pp	1
11	1491211	SpinX rugbuis opn.rh36 -40 links ps	1
12	1491210	SpinX rugbuis opn.rh36 -40 rechts ps	1
13	1489022	SpinX rugbuis opn.rh41- 49 links ps	1
14	1489021	SpinX rugbuis opn.rh41- 49 rechts ps	1
15	5226192	(WIRE BACKREST LOCK 36 SPINX SW 35,5)	1
16	5226193	Rea rug ontgr.touwteje zb38 ps XLT/SpinX	1
17	5226194	Rea rug ontgr.touwteje zb41 ps xlt/SpinX	1
18	5226195	Rea rug ontgr.touwteje zb43 ps XLT/SpinX	1
19	5226196	Rea rug ontgr.touwteje zb46 ps XLT/SpinX	1
20	5226411	Rea rug ontgr.touwteje zb48 ps XLT/SpinX	1
21	5226197	Rea rug ontgr.touwteje zb50 ps xlt/SpinX	1
22	1543736	(FIXATION, BACKREST UPHOLSTERY)	1
23	SP1429301	Rea duwh.opb.kort ex.bev ps	1
24	SP1584148	(Label pair, backrest height)	1


(Backrest tube with external push handles)



1	1491208	SpinX rugb.tbv.opb.duwh. kort	1
2	SP1574655	(Backrest tube for additional push handle)	1
3	1445825	SpinX duwh.opb.kort ps	1
4	1434743	Focus duwh.opb.ps	1
5	1434841	Focus duwh.opb.sterknop ps	1
6	SP1429301	Rea duwh.opb.kort ex.bev ps	1


Invacare® SpinX | (Backrest)
(Backrest upholstery, fixed)

355	1
380	2
405	3
430	4
455	5
480	6
505	7




1	1495839	(Upholstery, Backrest Fixed 35.5 Cpl, Jemima)	1
2	1495840	Focus rugbekl.zb38 zwart jemima tr18	1
3	1495841	(Upholstery, Backrest Fixed 40.5 Cpl, Jemima)	1
4	1495842	Focus rugbekl.zb43 zwart jemima tr18	1
5	1495843	Focus rugbekl.zb45,5 zwart jemima tr18	1
6	1495844	Focus rugbekl.zb48 zwart jemima TR18	1
7	1495845	Focus rugbekl.zb50,5 zwart jemima TR18	1

(Backrest upholstery, adj.)



1	1542540-TR18	Focus rug nasp.+hoes zb38 zwart jemima	1
2	1542541-TR18	Focus rug nasp.+hoes zb40,5 zwart jemima	1
3	1542542-TR18	Focus rug nasp.+hoes zb43 zwart jemima	1
4	1542543-TR18	Focus rug nasp.+hoes zb45,5 zwart jemima	1
5	1542544-TR18	Focus rug nasp.+hoes zb48 zwart jemima	1
6	1542545-TR18	Focus rug nasp.+hoes zb50,5 zwart jemima	1
7	1540578-TR18	Rughoes zb38 tbv.spanbanden jemima Rea/Poir	1
8	1540579-TR18	Rughoes zb41 tbv.spanbanden jemima Rea/Poir	1
9	1540580-TR18	Rughoes zb43 tbv.spanbanden jemima Rea/Poir	1
10	1540581-TR18	Rughoes zb46 tbv.spanbanden jemima Rea/Poir	1
11	1540582-TR18	Rughoes zb48 tbv.spanbanden jemima Rea/Poir	1
12	1540583-TR18	Rughoes zb50 tbv.spanbanden jemima Rea/Poir	1
13	1435415	Focus rugbekl.nasp.ex. hoes zb38-43 ps	1
14	1435416	Focus rugbekl.nasp.ex. hoes zb48-53 ps	1


(Backrest MISTRAL2, tension adj. (->2012.04))



1435761

1	1435761	Focus rug Mistral2 zb38 grijs dartex cpl	1
2	1478566	Rea rug nasp.+Mistral kussen cpl zb40,5	1
3	1435762	Focus rug Mistral2 zb43 grijs dartex cpl	1
4	1478567	Rea rug nasp.+Mistral kussen cpl zb45,5	1
5	1435763	Focus rug Mistral2 zb48 grijs dartex cpl	1
6	1478568	Rea rug nasp.+Mistral kussen cpl zb50,5	1
7	1435765	Focus rug Mistral2 zb38 grijs velours cpl	1
8	1478569	(BACKREST "MISTRAL2", 40,5, TR32)	1
9	1435766	Focus rug Mistral2 zb43 grijs velours cpl	1
10	1478570	(BACKREST "MISTRAL2", 45,5 TR32)	1
11	1435767	Focus rug Mistral2 zb48 grijs velours cpl	1
12	1478571	(BACKREST "MISTRAL2", 50,5, TR32)	1
13	1435764	Focus 150 rug Mistral2 zb53 grijs dartex cpl	1
14	1435768	Focus 150 rug Mistral2 zb53 grijs velours cpl	1
15	1431668	Rea rugkussen 355x510 Mistral dartex	1
16	1424651	Bellis rugkussen zb41 Mistral II grijs dartex	1
17	1424652	Bellis rugkussen zb46 Mistral II grijs dartex	1
18	1424653	Bellis rugkussen zb51 Mistral II grijs dartex	1
19	1431669	(Cushion, Backrest Mistral2 355X510, Plush Grey)	1
20	1435404	Focus rug Mistral zb40,5 grijs velours ps	1
21	1435405	(Cushion, Backrest Mistral2 455X510, Plush Grey)	1
22	1435406	(Cushion, Backrest Mistral2 505X510, Plush Grey)	1
23	1435867	Focus rughoes Mistral2 zb38	1
24	1435868	Focus rughoes Mistral2 zb43	1
25	1435869	Focus rughoes Mistral2 zb48	1
26	1435870	Focus rughoes Mistral2 zb53-55,5	1
27	1435415	Focus rugbekl.nasp.ex. hoes zb38-43 ps	1
28	1435416	Focus rugbekl.nasp.ex. hoes zb48-53 ps	1


(Legrest 80° Alu2003, Fixed)



SFS1493057

1	1493055-0501	Beenst.bovendl.80° rood ps A3/Ass/Bel/Comf/Focus	1
1	1493055-0507	Beenst.bovendl.80° blauw ps A3/Ass/Bel/Comf/Focus	1
1	1493055-0509	Beenst.bovendl.80° nw ps A3/Ass/Bel/Comf/Focus	1
1	1493055-0514	Beenst.bovendl.80° antzw ps cpl A3/As/Bel/Com/Foc	1
1	1493055-0515	Beenst.bovendl.80° grijs A3/Comf/Foc ps	1
1	1493055-0517	Beenst.bovendl.80	1
1	1493055-0518	(HANGER ALU2003, FIXED 80 DEG, ORANGE)	1
2	1424571	BEENST.BOVDL.VAST SCHROEF M10x23 A3/FOC PS	1


1





SPS140206


1	1493056-0501	Beenst.bovendl.90° nw ps rood A3/Ass/Bel/Comf/Foc	1
2	1424571	BEENST.BOVELDL.VAST SCHROEF M10x23 A3/FOC PS	1
1	1493056-0507	Beenst.bovendl.90° nw ps blauw A3/Ass/Bel/Com/Foc	1
1	1493056-0509	Beenst.bovendl.90° nw ps zwart A3/Ass/Bel/Com/Foc	1
1	1493056-0514	Beenst.bovendl.90° nw ps antra.A3/Ass/Bel/Com/Foc	1
1	1493056-0515	Beenst.bovendl.90° nw ps grijs A3/Ass/Bel/Com/Foc	1
1	1493056-0517	Beenst.bovendl.90° nw ps groen A3/Ass/Bel/Com/Foc	1
1	1493056-0518	Beenst.bovendl.90° nw ps oranje A3/Ass/Bel/Com/Fo	1
1	1424571	BEENST.BOVELDL.VAST SCHROEF M10x23 A3/FOC PS	1

(Legrest AA)



1	1524990	Rea kuitst.10' tbv.alu beenst.links ps	1
2	1524991	Rea kuitst.10' tbv.alu beenst.rechts ps	1
3	1524992	(Hanger (+ Calfpad 11.5), Le, Alu2008 Legrest)	1
4	1524993	(Hanger (+ Calfpad 11.5), Ri, Alu2008 Legrest)	1
5	1524994	Rea kuitst.25' tbv alu beenst.links ps	1
6	1524995	(Hanger (+ Calfpad 12), Ri, Alu2008 Legrest)	1
7	1525027	(Hanger (+ Calfpad 15), Le, Alu2008 Legrest)	1
8	1525028	(Hanger (+ Calfpad 15), Ri, Alu2008 Legrest)	1
9	1525029	(Hanger (+ Calfpad 18), Le, Alu2008 Legrest)	1
10	1525030	(Hanger (+ Calfpad 18), Ri, Alu2008 Legrest)	1
11	1508945	Beenst.bovendl.comf.hv li Rea/Poir	1
12	1508948	Beenst.bovendl.comf hv re Rea/Poir	1
13	1518616	Rea kuitst.zb36 ps	1
14	1518617	Rea kuitst.11,5' tbv alu beenst.ps	1
15	1518531	Rea kuitst.16'tbv.alu beenst.ps	1
16	1518532	Rea kuitst.19'tbv.alu beenst.ps	1
17	1518533	Rea kuitst.22'tbv.alu beenst.ps	1
18	5106406	(Pad, Calf Support 10)	1
19	5104185	\$ctiv kuitst.pad small ps	1
20	5104186	\$ctiv kuitst.pad zb46 ps	1
21	5104187	\$ctiv kuitst.pad zb50 ps	1
22	1429215	Rea kuitst.handwiel ps	1
23	1508943	Beenst.bovendl.comf.hv. stelknop ps Rea/Poir	1
24	1428336	Activ beenst.alu.afdekkap foampad ps A3/As/Bel	1
25	1424571	BEENST.BOVENDL.VAST SCHROEF M10x23 A3/FOC PS	1
26	1445811	Beenst.bovendl.hv.kniepa d li cpl A3/Aass/Bel/Foc	1
27	1445812	Beenst.bovendl.hv.kniepa d re cpl A3/Aass/Bel/Foc	1
28	1429326	Rea beenst.bovendl.hv.sterknop schr.M10x23 ps	1

	SW (cm)		
Focus, Focus 150	38-45,5	1	2
	48-55,5	3	4
Spirea ²	38-43	1	2
Spin ^x , XLT swing	45,5-50,5	3	4




SPS1429181

1	1429192-0509	Amp.st.11 li cpl Focus	1
2	1429191-0509	Amp.st.11 re cpl Focus	1
3	1429211-0509	Amp.st.14-17 li cpl Focus	1
4	1429212-0509	Amp.st.14-17 re cpl Focus	1
5	5104667	\$ctiv amp.st.pad ps. 160x55x220	1
6	5104668	Activ amp.st.pad ps. Ad/Ass/Bel/Sp	1

Invacare® SpinX | (Legrests)
(Amputee legrest)


1 2
👤 👤



062941526

1	1442991	Amp.st.li cpl A3/Spx/Bel	1
2	1442990	Amp.st re cpl A3/Spx/Bel	1

(Footrest, fixed, Alu2003)




1	1428357	Rea voetpl.zb38-40+inst. 24 alu li Bel/Com/Foc/Sp	1
1	1428358	Rea voetpl.zb38-40+inst. 24 alu re Bel/Com/Foc/Sp	1
2	1428361	Rea voetpl.zb43-46+inst. 24 alu.li Bel/Com/Foc/Sp	1
2	1428362	Rea voetpl.zb43-46+inst. 24 alu.re Bel/Com/Foc/Sp	1
3	1428365	Rea voetpl.zb44-50+inst.24 alu.li Bel/Com/Foc/Sp	1
3	1428366	Rea voetpl.zb44-50+inst.24 alu.re Bel/Com/Foc/Sp	1
4	1435756	Voetpl.27+inst.b.41-53 kort li Rea	1
4	1435757	Voetpl.27+inst.b.41-53 kort re Rea	1
5	1428359	Rea voetpl.zb38-40+inst. 36 alu.li Bel/Com/Foc/Sp	1
5	1428360	Rea voetpl.zb38-40+inst. 36 alu.re Bel/Com/Foc/Sp	1
6	1428363	Rea voetpl.zb43-46+inst. 36 alu.li Bel/Com/Foc/Sp	1
6	1428364	Rea voetpl.zb43-46+inst. 36 alu.re Bel/Com/Foc/Sp	1
7	1428367	Rea voetpl.zb44-50+inst. 36 alu.li Bel/Com/Foc/Sp	1
7	SP1428368	(FOOTREST 24 "ALU2003",LOWER PART,FIXED,LONG,RIGHT)	1
8	1435659	Voetpl.27+inst.b.41-53 ang li A3/As/Bel/Com/Foc	1
8	1435660	Voetpl.27+inst.b.41-53 ang re A3/As/Bel/Com/Foc	1
9	1434561-0010	Voetpl.insteekb.kort A3/Focus ps	1
9	1434563-0010	Rea insteekbuis lang voor aa/fx ps	1
10	1428340	Voetpl.vast zb38-40 li A3/Foc	1
10	1428341	Voetpl.vast zb38-40 re A3/Foc	1
11	1424574	Rea voetpl.zb43-46 alu. li.Bel/Comf/Foc/Sp	1
11	1424575	Rea voetpl.zb43-46.re. Bel/Comf/Foc/Sp ex pijp	1
12	1428346	Rea voetpl.zb48-50,5 li	1
12	1428347	Rea voetpl.zb48-50,5 re	1
13	1428374	Rea voetpl.zb49-60,5 li ps	1
13	1428375	Rea voetpl.zb49-60,5 re ps	1
14	1430433	\$OETPL.STD.BEV.KIT CPL A3/FOC	1

(Footrest, AA, "Alu2003")

	mm	24 cm	36cm	
Rea Azalea				
340 0	4	11	18	
340 +5	4	11	18	
340 +10	5	12	19	
390 0	4	11	18	
390 +5	5	12	19	
390 +10	6*	13	18	
440 0	5	12	19	
440 +5	6*	13	18	
440 +10	7*	14	19	
490 0	6*	13	18	
490 +5	7*	14	19	
490 +10	7*	14	19	
Rea Bellis				
390 - 430	5	12	19	
440 - 480	7*	14*	19	
490 - 530	7*	14*	19	
Rea Clematis				
390	3	10	17	
440	5	12	19	
490	7*	14*	19	

	mm	24 cm	36cm	
Rea Focus				
380	2	9	16	
430	4	11	18	
480	6*	13*	18	
SpinX, XLT Swing				
330 - 355	1	8	15	
380 - 405	2	9	16	
430 - 455	4	11	18	
480	6*	13*	18	
505	7*	14*	19	
Rea Spirea², Rea Spirea⁴NG				
355	1	8	15	
380	2	9	16	
405	3	10	17	
430	4	11	18	
455	5	12	19	
480	6*	13*	18	
505	7*	14*	19	





09/2013/19

1	1485430	Rea voetpl.zb33-35,5 HV+ inst.24 alu li Sp/XLT/Sw	1
1	1485524	Rea voetpl.zb33-35,5 HV+ inst.24 alu re Sp/XLT/Sw	1
2	1429159	Rea voetpl.zb36-45 hv+in st.24 alu.li As/Co/Fo/Sp	1
2	1429160	Rea voetpl.zb36-45 hv+in st.24 alu.re As/Co/Fo/Sp	1
3	1429161	\$ellis voetpl.zb39 hv+inst.buis 24 alu.li	1
3	1429162	\$ellis voetpl.zb39 hv+inst.buis 24 alu.re	1
4	1429163	Rea voetpl.zb39-43hv+ins t.24 alu.li as/Bel/Co/Sp	1
4	1429164	Rea voetpl.zb39-43hv+ins t.24 alu.re.as/Bel/Co/Sp	1
5	1429165	Rea voetpl.zb44-50 hv+in st.24 alu.li Bel/Com/Sp	1
5	1429166	Rea voetpl.zb44-50 hv+in st.24 alu.re Bel/Com/Sp	1
6	1429167	Voetpl.hv.19+inst.b.32-4 4 li A3/Ass/Bel/Comf/Foc	1
6	1429168	Voetpl.hv.19+inst.b.32-4 4 re A3/Ass/Bel/Comf/Foc	1
7	1429169	Rea voetpl.zb50,5 hv+in st.24 alu.li As/Co/Spi	1
7	1429170	Rea voetpl.zb50,5 hv+in st.24 alu.re As/Co/Spi	1
8	1485428	Rea voetpl.zb33-35,5 HV+ inst.36 alu li Sp/XLT/Sw	1
8	1485429	Rea voetpl.zb33-35,5 HV+ inst.36 alu re Sp/XLT/Sw	1
9	1429171	Rea voetpl.zb36-45 hv+in st.36 alu.li As/Co/Fo/Sp	1
9	1429172	Rea voetpl.zb36-45 hv+in st.36 alu.re As/Co/Fo/Sp	1
10	1429173	Bellis voetpl.zb41 hv+in st.buis 36cm.alu.li	1
10	1429174	Bellis voetpl.zb41 hv+in st.buis 36cm.alu.re	1
11	1429175	Rea voetpl.zb39-43hv+ins t.36 alu.li as/Bel/Co/Sp	1
11	1429176	Rea voetpl.zb39-43hv+ins t.36 alu.re as/Bel/Co/Sp	1
12	1429177	Rea voetpl.zb44-51 hv+in st.36 alu.li Bel/Com/Sp	1
12	1429178	Rea voetpl.zb44-51 hv+in st.36 alu re Bel/Com/Sp	1
13	1429179	Voetpl.hv.19+inst.b.41-5 3 li A3/Ass/Bel/Comf/Foc	1
13	1429180	Voetpl.hv.19+inst.b.41-5 3 re A3/Ass/Bel/Comf/Foc	1
14	1429181	Activ voetpl.hv.20+inst. buis 41-53 li.A3/Ass/Bel	1
14	1429182	Activ voetpl.hv.20+inst. buis 41-53 re.A3/Ass/Bel	1
15	5106372	Activ voetpl.zb38 std.ex inst.buis ps	1
16	5106374	Activ voetpl.std.zb43 ex inst.b.ps A3/Ad/Ass/Bel	1
17	5106375	Activ voetpl.std.zb46 ex inst.b.ps.A3/Ad/Ass/Bel	1

18	5106594	Activ voetpl.d.i.zb48 ex inst.buis ps	1
19	5106376	Activ voetpl.zb46 std.ex inst.buis ps	1
20	1434561-0010	Voetpl.insteekb.kort A3/Focus ps	1
21	1434563-0010	Rea insteekbuis lang voor aa/fx ps	1
22	1430839	Voetpl.verst.zb38,5-45,5 bev.kit A3/Focus	1
23	1430840	Voetpl.verst.zb38,5-45,5 bev.kit A3/Focus	1

(Footrest 90°, Alu2003)


Rea Bellis			
S (39 - 43)	2	6	11
M (44 - 48)	3	7	12
L (49 - 53)	4	8	13
Rea Azalea, Rea Comfort			
34-39	2	6	11
44	3	7	12
49	4	8	13
Rea Clematis			
39	1	5	10
44	2	6	11
49	3	7	12
Rea Focus			
38	1	5	10
43	2	6	11
48	3	7	12
SpinX, XLT Swing, Rea Spirea²			
38 - 40,5	1	5	10
43 - 45,5	2	6	11
48 - 50 5	3	7	12



SPS1446130

1	1446130	Rea voetpl.zb38-40 inst. kort tbv.beenst.90gr li	1
1	1446131	Rea voetpl.zb38-40 inst. kort tbv.beenst.90gr re	1
2	1446132	Rea voetpl.zb43-45 inst. kort tbv.beenst.90gr li	1
2	1446133	Rea voetpl.zb43-45 inst. kort tbv.beenst.90gr re	1
3	1446134	Rea voetpl.zb48-50 inst. kort tbv.beenst.90gr li	1
3	1446135	Rea voetpl.zb48-50 inst. kort tbv.beenst.90gr re	1
4	1446159	(FOOTREST LOWER PART 18, AA, FOOTR. 90, SHORT, LH)	1
4	1446160	(FOOTREST LOWER PART 18, AA, FOOTR. 90, SHORT, RH)	1
5	1446136	Rea voetpl.zb38-40 inst. lang tbv.beenst.90gr li	1
5	1446137	Rea voetpl.zb38-40 inst. lang tbv.beenst.90gr re	1
6	1446138	Rea voetpl.zb43-45 inst. lang tbv.beenst.90gr li	1
6	1446139	Rea voetpl.zb43-45 inst. lang tbv.beenst.90gr re	1
7	1446140	Rea voetpl.zb48-50 inst. lang tbv.beenst.90gr li	1
7	1446141	Rea voetpl.zb48-50 inst. lang tbv.beenst.90gr re	1
8	1446161	(FOOTREST LOWER PART 18, AA, FOOTR. 90, LONG, LH)	1
8	1446162	(FOOTREST LOWER PART 18, AA, FOOTR. 90, LONG, RH)	1
9	1440035	Rea insteekbuis tbv.d.i. voetplaat lang 346 ps	1
9	1446094	(TUBE, LEGREST, SHORT (221))	1
10	5106371	Activ voetpl.zb36 std.ex inst.buis ps	1
11	5106373	Activ voetpl.zb41 std.ex inst.buis ps	1
12	5106375	Activ voetpl.std.zb46 ex inst.b.ps.A3/Ad/Ass/Bel	1
13	5106376	Activ voetpl.zb46 std.ex inst.buis ps	1
14	1497913	(SCREW KIT FOOTREST 90 ALU 2004)	1

SH (cm)	SW (cm)					
XLT swing SpinX Spirea^{2/3/4}	38	40,5	43	45,5	48	50,5
38 - 48	1	2	3	4	5	6
50.5 - 53	7	8	9	10	11	12
Focus	38	40,5	43	45,5	48	50,5
36.5 - 46.5	1	1	3	3	5	5
49 - 54	7	7	9	9	11	11




SPS 1496945

1	1496945	(One piece footboard, short tube, SW380 mm, depth and angle adjustable, hanger 80°)	1
2	1496946	Voetpl.1-dlg zb40,5+inst buis kort tbv.80° Fo/A4	1
3	1496947	Voetpl.1-dlg zb43+inst.buis kort tbv.80° Fo/A4	1
4	1496948	Voetpl.1-dlg ZB45,5+inst buis kort tbv 80Fo/A4	1
5	1496949	Voetpl.1-dlg zb48+inst.buis kort tbv.80° Fo/A4	1
6	1496953	Voetpl.1-dlg zb50,5+inst buis kort tbv.80° Fo/A4	1
7	1496954	(One piece footboard, long, SW380, depth and angle adjustable, hanger 80°)	1
8	1496955	Voetpl.1-dlg zb40,5+inst buis lang tbv.80° Fo/A4	1
9	1496982	Rea voetpl.1-dlg 80gr. zb43 + 45,5 afh.v.zh cpl	1
10	1497004	Rea voetpl.1-dlg D.I.80° zb45,5 zh50,5-53 cpl ps	1
11	1497005	Rea voetpl.1-dlg D.I. zb48 cpl	1
12	1497006	Rea voetpl.1-dlg d.i.80° zb50,5 cpl	1

(Footrest one piece 90°)


	SH (mm)						SW (mm)					
XLT Swing, Spin ^X , Spirea ^{NG}	380	405	430	455	480	505	380	405	430	455	480	505
410 - 480	1	2	3	4	5	6						
505 - 530	11	12	13	14	15	16						
Focus	380	405	430	455	480	505	380	405	430	455	480	505
365 - 465	1	1	3	3	5	5						
490 - 540	11	11	13	13	15	15						
Focus 150	530	580					530	580				
440 - 490	7	7										
490 - 540	17	17										



SP51491248


1	1491247	(FOOTREST ONE PIECE 197, ANGLE ADJ., SHORT)	1
2	1491248	\$ea voetpl.1-dlg inst. buis kort zb40,5	1
3	1493267	\$ea voetpl.1-dlg d.i.90	1
4	1493268	Rea voetpl.1-dlg d.i.90 gr.zb45,5 kort inst	1
5	1493269	(FOOTREST ONE PIECE 297, ANGLE ADJ., SHORT)	1
6	1493270	(FOOTREST ONE PIECE 322, ANGLE ADJ., SHORT)	1
7	1495821	(FOOTREST ONE PIECE 347, ANGLE ADJ., SHORT)	1
8	1495822	(FOOTREST ONE PIECE 372, ANGLE ADJ., SHORT)	1
9	1495823	(FOOTREST ONE PIECE 397, ANGLE ADJ., SHORT)	1
10	1495824	(FOOTREST ONE PIECE 422, ANGLE ADJ., SHORT)	1
11	1493349	Rea voetpl.1-dlg 90gr. zb38 H/Di lang cpl	1
12	1493350	(FOOTREST ONE PIECE 222, ANGLE ADJ., LONG)	1
13	1493351	(FOOTREST ONE PIECE 247, ANGLE ADJ., LONG)	1
14	1493358	\$ea voetpl.1-dlg 90gr.zb45,5 cpl	1
15	1493359	(FOOTREST ONE PIECE 297, ANGLE ADJ., LONG)	1
16	1493377	Rea voetpl.1-dlg 90gr. zb50,5 H/Di lang cpl	1
17	1495825	(FOOTREST ONE PIECE 347, ANGLE ADJ., LONG)	1
18	1495826	(FOOTREST ONE PIECE 372, ANGLE ADJ., LONG)	1
19	1495827	(FOOTREST ONE PIECE 397, ANGLE ADJ., LONG)	1
20	1495828	(FOOTREST ONE PIECE 422, ANGLE ADJ., LONG)	1

SPS1429214



(Heel strap)

1




SP55106433

5321579


Rea hielband set cpl zb40-50+bev pp

1




1	1435616	Focus arml.rio kort+t- stuk li	1
2	1435615	Focus arml.rio kort+t- stuk re	1
3	1435618	Focus arml.rio lang+t- stuk tbv.flipup armst.li	1
4	1435617	Focus arml.rio lang+t- stuk tbv.flipup armst.re	1
5	1435614	Focus armst.h.v.zijschot +bevestiging cpl li	1
6	1435613	Focus armst.h.v.zijschot +bevestiging cpl re	1
7	5106533	Rea inbus m12x1 (fijn) 150/703/Ad/Bel/Sp	1
8	1431731	Focus armst.h.v.opkl. bev.punt tbv.zijschot ps	1
9	1424569	Focus armst.h.v.schroef kxpt phil k50x20 ps	1
10	1431724	Focus armst.oph.bev. achter sch.pen ps	1
11	1435492	\$pirea armst.opkl.bev.ps	1
12	1435493	Focus armst.opkl.bev. achter cpl.ps	1
13	1485272	SpinX armsteun opkl.bev. achter cpl li	1
14	1485246	SpinX armsteun opkl.bev. achter cpl re	1
15	1512698	XLT armst.opkl.bev. achter pp cpl	1
16	SP1574539	Rea armlegger hoes kort zwart plush ps	1
17	SP1574609	Rea armlegger hoes gepolst. lang zw. plush	1

(Armrest Rio, Height Adj.)



1	1435729	Focus armst.rio ihv.kort pp	1
2	1435731	Focus armst.rio ihv.lang pp	1
3	1435461	Focus armst.rio h.v.kort li	1
4	1435487	Focus armst.rio lang ihv arml.li	1
5	1435460	Focus armst.rio h.v.kort re	1
6	1435486	Focus armst.rio lang ihv arml.re	1
7	1434676	Rea arml.rio kort ps Foc/Sp3/Sp3/a3	1
8	1508064	Rea arml.rio lang ps Foc/Sp3/Sp3	1
9	1434725	Focus kled.besch.tbv.Rio arml.	1
10	5724189	\$ea kled.besch.schroef pt-fx h k40x14 zn ps	1
11	1435388	Rea kled.besch.geleiderail ps	1
12	1440027	Rea schroef k50x25 tbv. rio arml.ps	1
14	SP1574539	Rea armlegger hoes kort zwart plush ps	1
15	SP1574609	Rea armlegger hoes gepolst. lang zw. plush	1
13	1560734	Rea armst.vergrendel/ stokhouder veer ps	1





SPS1434680

1508028

SpinX armst.bev blok ps

1



(Sideguard)



SPS146127

1	1514484	(SIDE GUARD, PAIR)	1
2	1491238	SpinX kled.besch.li cpl	1
3	1491237	SpinX kled.besch.re cpl	1
4	1434725	Focus kled.besch.tbv.Rio arml.	1
5	1435388	Rea kled.besch.geleiderail ps	1
6	1484543-0004	Rea kled.besch.buis ps SpinX/XLT2008	1
7	1488504	Rea stop tbv kled.besch bev.ps	1
8	5723007	\$ea schroef m6x10 tbv.t- stuk ps	1
9	5724189	\$ea kled.besch.schroef pt-fx h k40x14 zn ps	1
10	1560734	Rea armst.vergrendel/ stokhouder veer ps	1


	SH (cm)	L	R
	41 - 48	1	2
	50,5 - 53	3	4



SPS 1484562

1	1484563	SpinX spatbord hoog li carbon look zh50,5-53 ps	1
2	1484562	SpinX spatbord hoog re carbon look zh50,5-53 ps	1
3	1484737	SpinX spatbord laag li carbon look	1
4	1484564	SpinX spatbord laag re carbon look	1

(Castor wheels)




	H/W	8	9	13		H/W	15	
	110/60					75/40		
	150/60	180x40 (7") soft	180x40 (7") pneu.	140x40 (5,5") soft			100x27 (4") mh	
	H/W	6	7	10	11	12	14	15
	110/40							
	150/40	200x32 (8") pneu.	200x27 (8") mh	150x32 (6") pneu.	150x30 (6") soft	150x30 (6") mh	125x27 (5") mh	100x27 (4") mh

ISPS1491244

1	1493468	(FORK 75 N F CPL, SPINX, PAIR)	1
2	1493439	(FORK 110 F N CPL, SPINX, PAIR)	1
3	1493449	(FORK 110 W F CPL, SPINX, PAIR)	1
4	1493393	SpinX voorw.vork smal 150 grijs pp	1
5	1493438	Rea voorw.vork 150x60 pp	1
6	5322005	Wiel 8'PN.200x32 ps Rea gebr.5320680	1
7	5322004	Wiel 8'semi-pn.200x27 ps gebr.5320679	1
8	5750781	Wiel 7'flexel 180x40 ps Rea	1
9	5750782	Wiel 7'pn.180x40 ps Rea gebr.5750788	1
10	5322006	Wiel 6'pn.150x32 ps	1
11	5750806	WIEL 6"FLEXEL 150x30 PS DIO	1
12	1555797	Wiel 6'semi-pn.150x30+as ps Rea	1
13	SP1601613	(CASTOR 140x40 (5,5") SOFT)	1
14	1491228	Wiel 5'softroll 125x30 zwart ex.as Rea	1
15	5750890	Wiel 4'soft 100x27 ps	1

(Wheel 24"-22-20)



1	1524030	Wiel 24'prof+21'alu ex. as ps Rea	1
2	1524186	(Wheel 24"x1"3/8 HD Notched Tyre, Handrim, Plastic)	1
3	1530764	Wiel 24' pn+hoepel staal Rea ps	1
4	1531508	Wiel 24'prof+foam coated hoepel Rea 2010	1
5	1531563	(Wheel 24"x1"3/8 HD Notched Tyre, Handrim, Blakgrip)	1
6	1524029	Wiel 24'vis+alu ex.as ps Rea >2010	1
7	1524163	(Wheel 24"x1"3/8 HD Thin Tyre, Handrim, Plastic)	1
8	1530763	(Wheel 24"x1"3/8 HD Thin Tyre, Handrim, Stainless)	1
9	1531462	Wiel 24' visprof.hoepel foam (2010) Rea	1
10	1531538	(Wheel 24"x1"3/8 HD Thin Tyre Profile, Handrim, Blackgrip)	1
11	1524050	Wiel 24'semi-pn.lijnprofiel alu hoepel cpl	1
12	1524188	(Wheel 24"x1"3/8 HD Solid , Handrim, Plastic)	1
13	1530762	(Wheel 24"x1"3/8 HD Solid , Handrim,Stainless)	1
14	1531531	(Wheel 24"x1"3/8 HD Solid , Handrim, Foam)	1
15	1531564	(Wheel 24"x1"3/8 HD Solid , Handrim, Blackgrip)	1
16	1524150	Wiel 24'semi-pn blokprofiel alu hoepel cpl rea	1
17	1524189	(Wheel 24"x1"3/8 HD Punctureproof, Handrim, Plastic)	1
18	1531129	(Wheel 24"x1"3/8 HD Punctureproof, Handrim, Stainless)	1
19	1531532	(Wheel 24"x1"3/8 HD Punctureproof, Handrim, Foam)	1
20	1531565	(Wheel 24"x1"3/8 HD Punctureproof, Handrim, Blackgrip)	1
21	1524049	(Wheel 24"x1" HD Thin Tyre, Handrim, Aluminum)	1
22	1524187	(Wheel 24"x1" HD Thin Tyre, Handrim, Plastic)	1
23	1531063	(Wheel 24"x1" HD Thin Tyre, Handrim, Stainless)	1
24	1531533	Wiel 24' low prof (1') met foam hoepel Rea	1
25	1531566	(Wheel 24"x1" HD Thin Tyre, Handrim, Blackgrip)	1
26	1532143	Wiel 24'+marathon buba alu hoepel ps rea	1
27	1532183	(Wheel 24"x1" HD Black Marathon, Handrim, Plastic)	1
28	1532144	(Wheel 24"x1" HD Black Marathon, Handrim, Inox)	1
29	1532184	(Wheel 24"x1" HD Black Marathon, Handrim, Foam)	1
30	1532185	Wiel 24'+marathon buba black grip hoepel ps Rea	1
31	1523949	Wiel 22' prof + alu hoepel rea	1

(Wheel 24"-22-20)


32	1524024	(Wheel 22"x1"3/8 HD Notched Tyre, Handrim, Plastic)	1
33	1523948	Wiel 22'vis+19 alu hoepel ex.as 2010	1
34	1524023	(Wheel 22"x1" HD Thin Profile Tyre, Handrim, Plastic)	1
35	1524021	Wiel 22'semi pneu + alu hoepel Rea	1
36	1524027	(Wheel 22"x1"3/8 HD Solid , Handrim, Plastic)	1
37	1524022	(Wheel 22"x1"3/8 HD Punctureproof, Handrim, Aluminum)	1
38	1524028	(Wheel 22"x1"3/8 HD Punctureproof, Handrim, Plastic)	1
39	1524020	Wiel 22x1'pneu.incl.alu hoepel ps	1
40	1524026	(Wheel 22"x1" HD Thin Tyre, Handrim, Plastic)	1
41	1523944	(Wheel 20"x1"3/8 HD Notched Tyre, Handrim, Aluminum)	1
42	1523945	(Wheel 20"x1"3/8 HD Notched Tyre, Handrim, Plastic)	1
43	1523885	Wiel 20' vis+18 alu.ex as ps Rea	1
44	1523902	(Wheel 20"x1"3/8 HD Thin Profile Tyre, Handrim, Plastic)	1
45	5750019	Buba 24'prof.ps Rea	1
46	5750014	Buba 24'vis.ps Rea	1
47	5750015	Buba 24'semi-pn.ps Rea	1
48	5207046	\$iba 24" z.prof ps Dio	1
49	5750160	\$UBA 24"x1 VIS.NR.2369GPS.REA	1
50	1512742	Buba 24'x1' (25-540) zwart marathon Rea	1
51	1531455	Buba 22' prof ps Rea	1
52	1531456	Buba 22' non marking HP rea	1
53	5750172	Buba 22'semi-pn.ps Rea	1
55	1531454	(Pneumatic Tyre Thin Profile 22"x1" (25-501))	1
56	5750017	\$uba 20" prof.ps Rea	1
57	5750012	(Pneumatic Tyre Thin 20"x1" 3/8")	1
58	1531577	Hoepel 21'alu Rea 2010 rechte bev.	1
59	1531578	Hoepel 21'tbv 24' wiel plastic coated 2010	1
60	1531305	(Kit Handrim 24" Stainless 4 Fixations)	1
61	1531459	Hoepel 21'foam coated ps Rea	1
62	1531460	Hoepel 21' black grip ps Rea	1
63	1531579	Hoepel 19'tbv.22'wiel alu rechte bev.ps	1
64	1531580	Hoepel 19' tbv 22' wiel plastic gecoat 2010	1
65	1530271	(Kit Handrim 20"Aluminum 4 Fixations)	1
66	1530272	(Kit Handrim 20" Plastic 4 Fixations)	1
67	5102807	SPAAKBESCH.24'TRANSP.PS	1
68	5102800	Spaakbesch.22 transp.ps	1
69	5102795	SPAAKBESCH.20'.TRANSP.PS	1
70	5750036	Biba 24'ps Rea	1
70	5750161	\$iba 24" x1 nr.2939G.psRea	1
71	5750035	Biba 22'ps Rea	1
72	5750034	\$iba 20" ps Dio	1

(Kit Handrim 24" Stainless 4 Fixations)


3854	1528924	()	4
3854	1529757	(Handrim Stainless, 24" Wheel, Std2010)	1
3854	1531457	(Spacer 22mm, Handrim)	4
3854	1531458	(Spacer 5mm, Handrim)	4

(Wheel 26-25")

25"	1
26"	2




1	1508583	Wiel 25'vis+alu.ex.as Sp/XLT	1
2	1508584	Wiel 26'vis+alu.ex.as Sp/XLT	1
3	1517344	(TYRE 25X1" LOW PROFILE (25-559))	1
4	1517345	Buba 26'vis.ps Rea	1
5	1517346	(TUBE 25X1" CAR VALVE (25-559))	1
6	1517359	Biba 26'ps Rea	1



SPX-48324

1	1483924	Wiel 24'sport+21'alu ps cpl Rea	1
2	1532186	Wiel 24'sport+21'alu bl.marathon ps cpl R	1
3	1561767	Wiel 24' sport+21'hoepel staal ps cpl Rea	1

(Wheel 24"-22" DB)




951534416

1	1524437	Wiel 24' pneu alu hoepel tbv tr.rem Rea	1
2	1524514	Wiel 24'profiel PL coate d tbv.trommelrem ps Rea	1
3	1531080	(Wheel 24"x1"3/8 Hub Brake Notched tyre, Handrim,Stainless)	1
4	1524436	Wiel 24' met alu hoepel vis tbv trom rem Rea ps	1
5	1524513	(Wheel 24"x1"3/8 Hub Brake Thin Tyre,Handrim,Plastic)	1
6	1531064	(Wheel 24"x1"3/8 Hub Brake Thin Tyre,Handrim,Stainless)	1
7	1524447	Wiel 24' buba semipn alu hoepel tbv tr.rem Rea	1
8	1524515	(Wheel 24"x1"3/8 Hub Brake Solid,Handrim, Plastic)	1
9	1531127	(Wheel 24"x1"3/8 Hub Brake Solid,Handrim, Stainless)	1
10	1524512	(Wheel 24"x1"3/8 Hub Brake Punctureproof,Handrim,Aluminum)	1
11	1524516	(Wheel 24"x1"3/8 Hub Brake Punctureproof,Handrim,Plastic)	1
12	1531128	(Wheel 24"x1"3/8 Hub Brake Punctureproof,Handrim,Stainless)	1
13	1536613	(Wheel 24"x1" Hub Brake Thin Tyre Punctureproof,Handrim, Stainless)	1
14	1524290	(Wheel 22"x1"3/8 Hub Brake Notched Tyre,Handrim,Aluminum)	1
15	1524361	(Wheel 22"x1"3/8 Hub Brake Notched Tyre,Handrim,Plastic)	1
16	1524289	(Wheel 22"x1"3/8 Hub Brake Thin Tyre,Handrim,Aluminum)	1
17	1524337	(Wheel 22"x1"3/8 Hub Brake Thin Tyre,Handrim,Plastic)	1
18	1524291	(Wheel 22"x1"3/8 Hub Brake Solid,Handrim,Aluminum)	1
19	1524388	(Wheel 22"x1"3/8 Hub Brake Solid,Handrim,Plastic)	1
20	1524292	(Wheel 22"x1"3/8 Hub Brake Punctureproof,Handrim, Aluminum)	1
21	1524435	(Wheel 22"x1"3/8 Hub Brake Punctureproof,Handrim,Plastic)	1
22	1536612	(Wheel 22"x1" Hub Brake Thin Tyre Punctureproof,Handrim,Stainless)	1
23	5750019	Buba 24'prof.ps Rea	1
24	5750014	Buba 24'vis.ps Rea	1
25	5750015	Buba 24'semi-pn.ps Rea	1
26	5207046	\$iba 24" z.prof ps Dio	1
27	5750160	\$UBA 24"x1 VIS.NR.2369GPS.REA	1
28	1531455	Buba 22' prof ps Rea	1
29	1531456	Buba 22' non marking HP rea	1
30	5750172	Buba 22'semi-pn.ps Rea	1
32	1531454	(Pneumatic Tyre Thin Profile 22"x1" (25-501))	1


33	1531577	Hoepel 21'alu Rea 2010 rechte bev.	1
34	1531578	Hoepel 21'tbv 24' wiel plastic coated 2010	1
35	1531305	(Kit Handrim 24" Stainless 4 Fixations)	1
36	1531579	Hoepel 19'tbv.22'wiel alu rechte bev.ps	1
37	1531580	Hoepel 19' tbv 22' wiel plastic gecoat 2010	1
38	1536546	(Kit Handrim 22"Stainless 4 Fixations)	1
39	5102807	SPAARBESCH.24'TRANSP.PS	1
40	5102800	Spaakbesch.22 transp.ps	1
41	5750036	Biba 24'ps Rea	1
42	5750035	Biba 22'ps Rea	1

(Brake)



SPX1492299

1	1542407	(Brake Cpl., Le, Incl attachment)	1
2	1542408	(Brake Cpl., Ri, Incl attachment)	1
3	1415485	Rea rem malta li ex.bev.Focus/SpinX	1
4	1415486	Rea rem malta re ex.bev.Focus/SpinX	1
5	1492299	SpinX rem bev.set op frame ps excl.rem	1







5203340

Remverl.univ.ps Ad/Foc/ Rea/Sp/XLT/Zoom

1

(One arm brake)


1 2
L R



525190200

1	1492107	SpinX rem 1-handsbed.li	1
2	1492000	SpinX rem 1-handsbed.re	1
3	1514787	(ATTACHMENT PAIR ONE ARM BRAKE, SPIREA/SPIN-X)	1
4	1518563	Rea rem 1-handsbed ex kabel ex bev links	1
5	1518564	Focus rem 1-hands bed re ex kabel ex erm links	1
6	1514742	Focus remkabel cpl tbv. 1-handsrem ps	1
7	5203340	Remverl.univ.ps Ad/Foc/ Rea/Sp/XLT/Zoom	1



(Drumbrake, assistant manoeuvred SpinX)



SFS153466


1	1536464	SpinX trommelrem cpl ex. wiel pp <2010-11	1
2	1536465	(BRAKE KIT 24" DB LEFT, CPL)	1
3	1536466	(BRAKE KIT 24" DB RIGHT, CPL)	1
4	1507940	(Kit, attachm. DB SpinX, Aktiv)	1
5	1536721	(KIT, WIRE, DB, L=830)	1
6	5100925	\$CHTERW.AS EDCO AFDEK-DOP PS	1
7	5750192	Rea QR assen tbv trommelrem ps	1
8	1536252	rea trommelrem naaf rechts	1
9	1536253	rea trommelrem naaf links	1
10	1536764	(DECAL, CABLE ASSEMBLY, BRAKE)	1

(Drumbrake assistant manoeuvred (=> 2010.11))



1	1485421	\$pinX trommelrem cpl ex. wielen pp	1
2	1488856	(DRUMBRAKE KIT, CPL, LE)	1
3	1488855	(DRUMBRAKE KIT, CPL, RI)	1
4	1507940	(Kit, attachm. DB SpinX, Aktiv)	1
5	1497616	(Kit, Wire, Db, L=780)	1
6	5100925	\$CHTERW.AS EDCO AFDEK-DOP PS	1
7	5221363	\$ea kant/rugverst.kabel klemdeel adapter ps	1
8	5750199	(Drumbrake Unit, Lh)	1
9	5750705	(Drumbrake Unit, Rh)	1
10	5750192	Rea QR assen tbv trommelrem ps	1
11	1517109	(KIT, SPACERS, BRAKE, 10 MM, PAIR)	1

1




1 1543662

SpinX heupgordel cpl ps


1

(Anti-tip device)




1	1483884	SpinX anti-tip li laag gemonteerd zh41-45,5 ps	1
1	1483885	SpinX anti-tip re laag gemonteerd zh41-45,5 ps	1
2	1483998	SpinX anti-tip re hoog gemonteerd zh48-53 ps	1
2	1483999	SpinX anti-tip li hoog gemonteerd zh48-53 ps	1

	 20-26"	 16"
Azalea	1	2
Dahlia	1	2
Clematis	1	-
Focus	4	6
Spin ^x	1	-
Spirea ² , Spirea ^{4NG}	4	5
XLT, XLT Max	3	-



1	SP5206655	(Reflector kit (Wheel Rear Spoked))	1
2	SP5254249	(Reflector kit (Wheel rear 16") Azalea Clematis)	1
3	5206640	Rea reflectorkit cpl. Comf/xlt	1
4	5104965	\$eflector set cpl >04/00	1
5	1435760	\$pirea2 reflector kit	1
6	1546559	(Reflector Kit Dahlia (wheel 16"))	1


(Backrest stay)



1	5106188-C9005	\$ctiv rugfix.zb36 cpl	1
2	5106399-C9005	Activ rugfix.zb38 cpl	1
3	5221696-C9005	Rea rugfix.zb40	1
4	5106331-C9005	Activ rugfix.zb43 cpl	1
5	5106190-C9005	Activ rugfix.zb46 cpl	1
6	1435484-9005	Activ rugfix.zb48 cpl	1
7	5106191-C9005	Activ rugfix.zb50 cpl	1
8	1435490-9005	Focus rugfix.zb53+bev. mat.cpl	1
9	5106299	(Sticker Warning (Se/Dk/Fi))	1
10	5106300	Activ rugfix.waarschuingssticker NL/DU/FR	1
11	5106301	(Sticker Warning (Gb/No/Sp))	1

(Step tube)


1



1	1493346	SpinX trapdop ps cpl	1
2	5716299	(PLASTIC CAP K 22 BLACK)	1
3	5716296	\$ea afdekdop plast.zwart 22x19,5 150/703/Ass/Comf	1

(Crutch holder)

1




SPS5106390

5106390

Rea stokhouder adapt/ SpinX ps

1




SPS575072€

5750729

Rea pump model 97 ps

1

(Table tray)



1	SP1570747	(TABLE TRAY SMALL INCL. MOUNTING KIT)	1
2	SP1570748	Rea werkblad medium incl montage cpl	1
3	SP1570749	Rea werkblad large incl. montage cpl	1
4	SP1570775	(TABLE TRAY X-LARGE INCL. MOUNTING KIT)	1
5	SP1567517	(TABLE TRAY SMALL)	1
6	SP1567518	(TABLE TRAY MEDIUM)	1
7	SP1567519	(TABLE TRAY LARGE)	1
8	SP1568862	(TABLE TRAY X-LARGE)	1
9	SP1570745	Dahlia werkblad montage set cpl	1
10	SP1572460	Rea werkblad vergrendel set	1
11	SP1570746	(TABLE CUSHION, SMALL)	1
12	SP1571269	(TABLE CUSHION, MEDIUM)	1
13	SP1571270	(TABLE CUSHION, LARGE)	1
14	SP1571271	(TABLE CUSHION, X-LARGE)	1



Yes, you can.®



Belgium & Luxemburg:

Invacare nv, Autobaan 22, B-8210 Loppem
Tel: (32) (0)50 83 10 10, Fax: (32) (0)50 83 10 11
belgium@invacare.com / www.invacare.be

Denmark:

Invacare A/S, Sdr. Ringvej 37, DK-2605 Brøndby
Tel: (45) (0)36 90 00 00, Fax: (45) (0)36 90 00 01
denmark@invacare.com / www.invacare.dk

Deutschland:

Invacare Aquatec GmbH, Alemannenstraße 10, D-88316 Isny
Tel: (49) (0)75 62 7 00 0, Fax: (49) (0)75 62 7 00 66
info@invacare-aquatec.com / www.invacare-aquatec.de

European Distributor Organisation:

Invacare, Kleiststraße 49, D-32457 Porta Westfalica
Tel: (49) (0)57 31 754 540, Fax: (49) (0)57 31 754 541
edo@invacare.com / www.invacare.eu.com

España:

Invacare SA, c/Areny s/n, Polígon Industrial de Celrà, E-17460 Celrà (Girona)
Tel: (34) (0)972 49 32 00, Fax: (34) (0)972 49 32 20
contactsp@invacare.com / www.invacare.es

France:

Invacare Poirier SAS, Route de St Roch, F-37230 Fondettes
Tel: (33) (0)2 47 62 64 66, Fax: (33) (0)2 47 42 12 24
contactfr@invacare.com / www.invacare.fr

Ireland:

Invacare Ireland Ltd, Unit 5 Seatown Business Campus, Seatown Road, Swords, County Dublin - Ireland
Tel: (353) 1 810 7084, Fax: (353) 1 810 7085
ireland@invacare.com / www.invacare.ie

Italia:

Invacare Mecc San s.r.l., Via dei Pini 62, I-36016 Thiene (VI)
Tel: (39) 0445 38 00 59, Fax: (39) 0445 38 00 34
italia@invacare.com / www.invacare.it

Nederland:

Invacare BV, Celsiusstraat 46, NL-6716 BZ Ede
Tel: (31) (0)318 695 757, Fax: (31) (0)318 695 758
nederland@invacare.com, csede@invacare.com / www.invacare.nl

Norge:

Invacare AS, Grensesvingen 9, Postboks 6230, Etterstad, N-0603 Oslo
Tel: (47) (0)22 57 95 00, Fax: (47) (0)22 57 95 01
norway@invacare.com, island@invacare.com / www.invacare.no

Österreich:

Invacare Austria GmbH, Herzog Odilostrasse 101, A-5310 Mondsee
Tel: (43) 6232 5535 0, Fax: (43) 6232 5535 4
info-austria@invacare.com / www.invacare.at

Portugal:

Invacare Lda, Rua Estrada Velha, 949, P-4465-784 Leça do Balio
Tel: (351) (0)225 1059 46/47, Fax: (351) (0)225 1057 39
portugal@invacare.com / www.invacare.pt

Sverige & Suomi:

Invacare AB, Fagerstagatan 9, S-163 91 Spånga
Tel: (46) (0)8 761 70 90, Fax: (46) (0)8 761 81 08
sweden@invacare.com, finland@invacare.com / www.invacare.se

Switzerland:

Invacare AG, Benkenstrasse 260, CH-4108 Witterswil
Tel: (41) (0)61 487 70 80, Fax: (41) (0)61 487 70 81
switzerland@invacare.com / www.invacare.ch

United Kingdom:

Invacare Limited, Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ
Switchboard Tel: (44) (0)1656 776200, Fax: (44) (0)1656 776201
Customer services Tel: (44) (0)1656 776222, Fax: (44) (0)1656 776220
uk@invacare.com / www.invacare.co.uk