

Spare Parts Catalogue

Version: 2010/08/11a EN



Invacare® DRAGON



Yes, you can.®

General Information

For any order of spare parts the following information is required:

- Product type*
- Product name & model*
- Serial Number*
- Part-Number*
- Part-Description*

Exploded view drawings and code numbers are not contractual and subject to modification, in relation to Invacare® continuous product improvement policy to meet end users needs.

The reproduction or use of this work, or of any substantial part of this work, in any form or by any electronic, mechanical or other means, including photocopying, scanning and recording, and in any information storage and retrieval system, including electronic transmission, is forbidden without the prior written permission of Invacare® International Sarl.



Yes, you can.®

TABLE OF CONTENTS

Information

User Manuals Dragon	5
User Manuals REM24SD	6
User Manuals Shark	7
User Manuals Miscellaneous	8

Chassis

Anti-Tipper	9
Chassis	10
Shrouds	11

Seat

Seat Unit	12
Seatplate Dragon Start	13
Childrenseat	14
Seat Hight Adjusting	15
Seat Angle Adjusting	16
Seat Cushion Adults	17
Seat Cushion For Junior Seat	18

Backrest

Backrest Adjustment Pneumatic ==>11/2004	19
Backrest Adjustment 01/2005 ==>	20
Backrest Adjustment Setable	21
Settable Backrestrecliner Junior Seat	22
Upholstery	23

Legrests

Legrest Support Manuall	24
Legrest Support Electr. / Vari A 07.2005 ==>	25
Hanger manually height adjustable	26
Hanger Standard 80°	27
Vari A Hanger RH/LH 07.2005 ==>	28
ADM Hanger RH/LH Electr. 07.2005 ==>	29
Calf Plate Holder LH/RH 07.2005 ==>	30
Legrests Dragon Junior	31

Footrests

Footplate flip up 01.08.10 ==>	32
Footplate flip up, angle-adjustable 01.08.10 ==>	33
Legrests, manual	34
Footrest Details	35

Siderests

Side Panel for Standard-Seat ==> 07/31/08	36
Side Panel for Standard-Seat 08/01/08 ==>	37
Armrest Junior Seat	38



Yes, you can.®

TABLE OF CONTENTS

Front Wheels

Front Wheels	39
--------------------	----

Rear Wheels

Rear Wheel	40
------------------	----

Propulsion

Motor & Declutching 6/8 km/h	41
Motor & Declutching 6 km/h 16.07.08 ==>	42
Motor & Declutching 8 km/h 16.07.08 ==>	43

Controls Drive

Electronic	44
------------------	----

Control Adjustment

REM24SD Remote	45
Remote holder for accompanying person	46
Remote REM550 & Holder	47
Keypads Shark2	48
Dual Control	49
Remote Brackets	50
Control Unit/ Joystickbox	51
Remoteholder	52

Cables & Charger

Batterys & Charger	53
Cables	54

Safety Options

Front Lights	55
Rear Light	56
Belts	57

Seating Options

Headrest	58
Headrest with Cheeksupport	59

Other Options & Accessories

Options	60
Tray	61
Kerb Climber	62

CS ČESKY	FR FRANÇAIS	PL POLSKI
DA DANSK	HR HRVATSKI	PT PORTUGUÊS
DE DEUTSCH	HU MAGYAR	SK SLOVENČINA
EN ENGLISH	IT ITALIANO	SL SLOVENŠČINA
ES ESPAÑOL	LT LIETUVIŲ KALBA	SV SVENSKA
ET EESTI KEEL	NL NEDERLANDS	RO ROMÂNĂ
FI SUOMEN KIELI	NO NORSK	JA 日本語

Pos.	Part number	Description	Valid from → until	Qty
NL	1442476	User Manual Dragon Dutch		1
EN	1442477	User Manual Dragon English		1
FI	1442478	User Manual Dragon Finnish		1
FR	1442479	User Manual Dragon French		1
DE	1442480	User Manual Dragon German		1
IT	1442481	User Manual Dragon Italian		1
NO	1442482	User Manual Dragon Norwegian		1
PT	1442483	User Manual Dragon Portuguese		1
ES	1442484	User Manual Dragon Spanish		1
SV	1442485	User Manual Dragon Swedish		1
SL	1443130	User Manual Dragon Slovenian		1
SK	1443131	User Manual Dragon Slovak		1
PL	1443132	User Manual Dragon Polish		1
HU	1443133	User Manual Dragon Hungarian		1
CZ	1443134	User Manual Dragon Czech		1
LT	1493185	User Manual Dragon Lithuanian		1

User Manuals REM24SD

CS ČESKY	FR FRANÇAIS	PL POLSKI
DA DANSK	HR HRVATSKI	PT PORTUGUÊS
DE DEUTSCH	HU MAGYAR	SK SLOVENČINA
EN ENGLISH	IT ITALIANO	SL SLOVENŠČINA
ES ESPAÑOL	LT LIETUVIŲ KALBA	SV SVENSKA
ET EESTI KEEL	NL NEDERLANDS	RO ROMÂNĂ
FI SUOMEN KIELI	NO NORSK	JA 日本語

Pos.	Part number	Description	Valid from → until	Qty
DE	1436142	User Manual REM24SD German		1
FR	1436143	User Manual REM24SD French		1
EN	1436144	User Manual REM24SD English		1
SV	1436145	User Manual REM24SD Swedish		1
DA	1436146	User Manual REM24SD Danish		1
NO	1436147	User Manual REM24SD Norwegian		1
NL	1436148	User Manual REM24SD Dutch		1
ES	1436149	User Manual REM24SD Spanish		1
PT	1436150	User Manual REM24SD Portuguese		1
IT	1436385	User Manual REM24SD Italian		1
FI	1436386	User Manual REM24SD Finnish		1

CS ČESKY	FR FRANÇAIS	PL POLSKI
DA DANSK	HR HRVATSKI	PT PORTUGUÊS
DE DEUTSCH	HU MAGYAR	SK SLOVENČINA
EN ENGLISH	IT ITALIANO	SL SLOVENŠČINA
ES ESPAÑOL	LT LIETUVIŲ KALBA	SV SVENSKA
ET EESTI KEEL	NL NEDERLANDS	RO ROMÂNĂ
FI SUOMEN KIELI	NO NORSK	JA 日本語

Pos.	Part number	Description	Valid from → until	Qty
NL	1436157	User Manual Shark Dutch		1
SV	1436154	User Manual Shark Swedish		1
ES	1436158	User Manual Shark Spanish		1
DE	1436151	User Manual Shark German		1
IT	1436383	User Manual Shark Italian		1
PT	1436159	User Manual Shark Portuguese		1
DA	1436155	User Manual Shark Danish		1
FI	1436384	User Manual Shark Finnish		1
EN	1436153	User Manual Shark English		1
NO	1436156	User Manual Shark Norwegian		1
FR	1436152	User Manual Shark French		1
1	1436151	User Manual Shark German		1
2	1436153	User Manual Shark English		1
4	1436157	User Manual Shark Dutch		1
5	1436158	User Manual Shark Spanish		1
6	1436159	User Manual Shark Portuguese		1
9	1436156	User Manual Shark Norwegian		1
10	1436155	User Manual Shark Danish		1
11	1436383	User Manual Shark Italian		1
8	1436154	User Manual Shark Swedish		1


User Manuals Miscellaneous

CS	ČESKY	FR	FRANÇAIS	PL	POLSKI
DA	DANSK	HR	HRVATSKI	PT	PORTUGUÊS
DE	DEUTSCH	HU	MAGYAR	SK	SLOVENČINA
EN	ENGLISH	IT	ITALIANO	SL	SLOVENŠČINA
ES	ESPAÑOL	LT	LIETUVIŲ KALBA	SV	SVENSKA
ET	EESTI KEEL	NL	NEDERLANDS	RO	ROMÂNĂ
FI	SUOMEN KIELI	NO	NORSK	JA	日本語

Pos.	Part number	Description	Valid from → until	Qty
DE	1419316	User Manual Handheld Programmer German		1
EN	1419317	User Manual Handheld Programmer English		1
DE	Q99081	User Manual Table- & Middle Control REM41-S German		1
FR	Q99082	User Manual Chin Control, French		1
DE	Q99083	User Manual Chin Control, German		1
ED	Q99086	User Manual Chin Control, English		1

Anti-Tipper


1440628B.dwg




Pos.	Part number	Description	Valid from → until	Qty
1	1440743-9005	Anti-Tipper Support Tube Black		1
2	M14906	Anti-Tipper Wheel VG80/24-GI10-SS		1
3	SP1440686	Anti-Tipper Complete		1
4	1442593	Mounting Kit Anti-Tip		1

Chassis

1442035.dwg




Pos.	Part number	Description	Valid from → until	Qty
1	SP1442035	Chassis, blue, cpl.		1
1	SP1448642	Chassis, black, cpl.		1
1	SP1536283	Chassis, silver, cpl.		1
1	SP1536284	Chassis, transparent blue, cpl.		1
1	SP1536285	Chassis, transparent red, cpl.		1
1	SP1536286	Chassis, transparent purpur, cpl.		1
2	1419254	Angular Ball Bearing 7202B 2Rs		1
3	1418547-9005	Bearing Plate LH (black)		1
3	1418544-9005	Bearing Plate RH (black)		1
4	1492498	Battery Frame ACS		1
5	SP1440745	Luggage Carrier cpl.		1
7	SP1492647	Taxi Fixation cmpl.		1
-	1442551	Fixing Kit Luggage Carrier cpl.		1
-	1442601	Battery Clamp Fixing Kit		1

Shrouds


Pos.	Part number	Description	Valid from → until	Qty
1	SP1442599	Battery Cover Rear Complete		1
2	1440684-9006	Battery Cover Front		1
3	1440769	Hand Wheel		1
4	SP1442592	Facing for Chassis LH cpl.		1
5	SP1442596	Facing for Chassis RH cpl.		1

Seat Unit



Pos.	Part number	Description	Valid from → until	Qty
2	1442655-9005	Seat Frame For Electrical Tilt		1
3	1419999-9005	Seat Plate		1
4	P31298	Reflector, Yellow		1
5	U21087	Velcro Fastener Strip 38 X 400mm Adhesive		1
6	SP1423205	Seat Mounting Kit (One-Piece Version)		1


1450311.dwg




Pos.	Part number	Description	Valid from → until	Qty
1	SPP31298	Reflector, yellow		1
2	U21087	Velcro Fastener Strip 38 X 400mm Adhesive		1
3	-	HEX NUT M8		1
4	-	Hex Head Screw M8 X 20		1
5	1442753-9005	Bracket f.siederest lh Dragon Junior		1
6	1442754-9005	Bracket f.siederest rh Dragon Junior		1
7	1442761-9005	Seat plate cpl. std Dragon Junior		1

Childrenseat

1444900X.dwg




Pos.	Part number	Description	Valid from → until	Qty
1	SP1446520	Actuator ACS2 / Shark		1
2	1434609	KIT INDEXEUR DE R-J+FIX.A2000 JUNIOR		1
3	5009551	Plastic Upholstery Washer		1
4	F368016.99	Locking Bolt - Quick Release		1
5	SPP31298	Reflector, yellow		1
6	1443083	Part not for sale		1
7	-	HEX NUT M8		1
8	-	Nut, M8 DIN 985		1
9	-	Srew M8X35 DIN 933		1
10	-	Srew M8X35 DIN 933		1
11	-	Countersunk Screw W. Allen Key M8*16 10.9		1
12	5824282	M6X16 ISO7380 SCREW		1
13	SP1443076	Seat Cushion 360mm For Junior Seat		1
14	1443069-0002	Thread Sheet Dragon Junior		1
15	-	Lock Washer 10mm		1
16	-	Tube End Plug SI 22		1
17	P22414	CABLE STRAP L=142 MM; BLACK		1
18	P22417	CABLE FASTENER 7.9 MM (CORD CLIP); BLACK		1
19	D16108	Countersunk Self-Tapping Cross Head Screw 4,8 X 22		1
20	G00107	Scheibe 8, 4X17X1, 6 DIN 125		1
21	J50999-0002	Bolt		1
22	-	Bush Rd. 12 X 1.8 X 23		1
23	SP1448718	Upholstery cpl For Junior Seat		1
24	1442761-9005	Seat plate cpl. std Dragon Junior		1
25	SP1448703	Side part LH cpl. Dragon Junior		1
26	SP1448702	Side part RH cpl. Dragon Junior		1
27	SP1448712	Backresttube		1
28	SP1445023	Support bracket		1



Pos.	Part number	Description	Valid from → until	Qty
1	1440733-9005	Plate Black Right		1
1	1440763-9005	Plate Black Left		1
2	1440721-9005	Fixing Plate		1
3	1442594	Mounting Kit Lifter And Rear Seat Fixation		1

Seat Angle Adjusting



Pos.	Part number	Description	Valid from → until	Qty
-	SP1440716	Manual Seat Angle Adjusting Complete		1
-	SP1442796	Electric Seat Angle Adjuster ; cpl.		1
1	1440744-9005	Cross Tube	→ 20050530	1
1	1449085-9005	Tube	20050601 →	1
2	1440718-9005	Top Fixing Plate		1
3	SP1446521	Actuator La 2825S-75-24-002		1
4	1442800-9005	Tilt Suspension		1
5	1442595	Mounting Kit Tilt And Front Seat Fixation		1
6	1446515	Mounting Kit Actuator Seat Angle		1

1420460DR.dwg

STANDARD

KONTUR

LAGUNA

FLO-TECH LITE VISCO

1	41cm	black
2	46cm	
3	51cm	

4	43cm	black
---	------	-------

5	41,5cm
---	--------



6	40,5 x 43cm
---	-------------

Pos.	Part number	Description	Valid from → until	Qty
1	1420460	Cushion Standard SD 41cm		1
2	1420461	Cushion Standard SD 46cm		1
3	1420462	Cushion Standard SD 51cm		1
4	B260005	Seat Cushion Contour		1
5	1435857	CUSHION, SEAT, "MISTRAL2" 415/490, DARTEX GREY		1
6	SP1501110	Flotech Lite Visco Seat Cushion (40.5 X 43cm)		1


Seat Cushion For Junior Seat

1443076.dwg


STANDARD



Pos.	Part number	Description	Valid from → until	Qty
1	SP1443076	Seat Cushion 360mm For Junior Seat		1


Backrest Adjustment Pneumatic ==>11/2004

Pos.	Part number	Description	Valid from → until	Qty
a	-	Axle Bolt M8 X 45		1
c	G00107	Scheibe 8, 4X17X1, 6 DIN 125		1
b	-	Nut, M8 DIN 985		1
-	SP1436509	Backrest Adjustment Pneumatic For Standard Seat Complete		1
d	Z131009A	Bush Rd. 18 X 3 X 22		1
f	-	Washer 8,4 X 25 X 2		1
e	-	Bush Rd. 12 X 1.8 X 23		1
g	L01041	Cover Cap For Wheel Invacare P452E		1
1	1442845-9005	Lower Receiving Bracket		1
2	1440395-9005	Backtube		1
3	N368003	Release Lever Stabilus		1
4	N368007	Back Adjuster Pneumatic cpl.		1
5	F367005	Fixing Kit For Seat Mounting		1

Backrest Adjustment 01/2005 ==>


Pos.	Part number	Description	Valid from → until	Qty
-	SP1444993	Back Adjuster Laguna Pneumatic, Without Back Cushion, cpl.		1
-	SP1443268	Electrical Back Adjuster Laguna, Without Back Cushion, cpl.		1
-	SP1443267	Electrical Back Adjuster Contur, Without Back Cushion, cpl.		1
-	SP1444991	Back Adjustment Standard, Back Height 54CM, Pneumatically, cpl with backbelt		1
-	SP1443269	Backrest Adjustment Pneumatic For Standard Seat Complete		1
-	SP1443270	Back Adjustment Standard, Back Height 54cm, Electzrically, Cpl with backbelt		1
-	SP1444990	Backrest Adjustment Pneumatic For Standard Seat Complete		1
-	SP1444992	Back Adjuster Kontur Pneumatic, Without Back Cushion, cpl.		1
1	1442866-9005	Back Tube Standard, Back Height 54 cm, Infinitely Adjustable		1
1	1442833-9005	Back Tube Contour, Infinitely Adjustable		1
1	1442864-9005	Back Tube Standard, Back Height 48 cm, Infinitely Adjustable		1
2	1442874-9005	Back Tube Laguna, Infinitely Adjustable		1
3	1442845-9005	Lower Receiving Bracket		1
4	N368007	Back Adjuster Pneumatic cpl.		1
5	N368003	Release Lever Stabilus		1
6	SP1446520	Actuator ACS2 / Shark		1
7	1446702	Backtube Fixation Kit Laguna		1
7	1446699	Backtube Fixation Kit		1
8	1446700	Mounting Kit Actuator		1

Backrest Adjustment Setable




Pos.	Part number	Description	Valid from → until	Qty
1	1434878-9005	Backtube Contour For Storm ³ Euro		1
1	1419996-9005	Back Tube 480 mm		1
1	1420465-9005	Back Tube 540 mm		1
2	SP1423208	Std. Backrest Assembly Kit cpl.		1
3	SP1420444	Star Grip M8 cpl.		1
4	SP1444989	Back Adjuster Kontur, Without Back Cushion, Progressively Adjustable cpl.		1
4	SP1444988	Back Adjuster Standard, Back Height 54 cm, With Backrest Upholstery,		1
4	SP1444987	Back Adjuster Standard, Back Height 48 cm, With Backrest Upholstery,		1

Settable Backrestrecliner Junior Seat




Pos.	Part number	Description	Valid from → until	Qty
-	-	Backrest Dragon Junior cpl.		1
1	5009551	Plastic Upholstery Washer		1
2	-	Nut, M8 DIN 985		1
3	-	Srew M8X35 DIN 933		1
4	-	Countersunk Screw W. Allen Key M8*16 10.9		1
5	1443069-0002	Thread Sheet Dragon Junior		1
6	1442925-9005	Back Adjustment Dragon Junior		1
7	D16108	Countersunk Self-Tapping Cross Head Screw 4,8 X 22		1
8	G00107	Scheibe 8, 4X17X1, 6 DIN 125		1
9	-	Bush Rd. 12 X 1.8 X 23		1
10	SP1448718	Upholstery cpl For Junior Seat		1
11	1442980-9005	Back Tube cpl.Std Kid. Dragon Junior		1
12	SP1420444	Star Grip M8 cpl.		1

1444993.dwg




Pos.	Part number	Description	Valid from → until	Qty
1	1422062	Sling Back with removable cushion 480 mm		1
1	1422063	Sling Back with removable cushion 540 mm		1
2	SPF365037	Back Cushion For Storm³ Euro cpl.		1
3,4	F367006	Fixing Kit For Backrest Comfort And Contur		1
5	1427119	CUSHION,BACKREST LAGUNA 402X550,DARTEX GREY		1
6	1420072	BACKREST SHELL "LAGUNA" 410,CPL		1
3	Z16001A	Clamp		1
4	D16108	Countersunk Self-Tapping Cross Head Screw 4,8 X 22		1

1450310.dwg




Pos.	Part number	Description	Valid from → until	Qty
1	SP1450310	Fixation for Legrest LH Dragon Junior Complet		1
1	SP1450309	Fixation for Legrest RH Dragon Junior Complet		1
2	SP1441526	Leg Support Carrier LH cpl.		1
2	SP1441527	Leg Support Carrier RH cpl.		1


Legrest Support Electr. / Vari A 07.2005 ==>

Pos.	Part number	Description	Valid from → until	Qty
1	SP1449737	Footrest/Legrest receptacle, manuell, LH, 8cm		1
1	SP1449739	Footrest/Legrest receptacle, powered, LH, 8cm		1
1	SP1449735	Footrest/Legrest receptacle, powered, LH, 5cm		1
1	SP1449733	Footrest/Legrest receptacle, manuell, LH, 5cm		1
2	SP1449738	Footrest/Legrest receptacle, powered, RH, 8cm		1
2	SP1449736	Footrest/Legrest receptacle, manuell, RH, 8cm		1
2	SP1449734	Footrest/Legrest receptacle, powered, RH, 5cm		1
2	SP1449732	Footrest/Legrest receptacle, manuell, RH, 5cm		1
3	SP1450228	Legrest Cable		1


Hanger manually height adjustable



Pos.	Part number	Description	Valid from → until	Qty
1	1513006	HANGER ALU2008, AA, RIGHT, CPL		1
1	1513007	HANGER ALU2008, AA, LEFT, CPL		1
2	1516314	CALF PAD 15 (165x170x80), CPL		1
2	1518532	Calf Support 15		1
3	SP1442853	Lower Legrest With Angle Adjustable Footplate, With Heel Strap. Width		1
3	SP1442854	Lower Legrest With Angle Adjustable Footplate, With Heel Strap. Width		1
4	1430402	HEELSTRAP COMPLETE		1




Pos.	Part number	Description	Valid from → until	Qty
1	1493055-0509	HANGER ALU2003, FIXED 80 DEG, CPL		1

Vari A Hanger RH/LH 07.2005 ==>


Pos.	Part number	Description	Valid from → until	Qty
-	SP1527820	Rework-Kit VariA-Legrest Version 1 ==> Version 2	→ 20090401	1
1	SP1449348	Hanger Vari A RH cpl.		1
1	SP1449349	Hanger Vari A LH cpl.		1
2	SPL32906	Excentric Clamping Lever		1
3	SP1434955	Cover Legrest cpl.		1
6	6301333	Legrest Calf Panel		1
7	SP1450447	Calf Pad Holder Left Turnable, cpl.		1
7	SP1450446	Calf Pad Holder Right Turnable, cpl.		1

ADM Hanger RH/LH Electr. 07.2005 ==>




Pos.	Part number	Description	Valid from → until	Qty
1	6301333	Legrest Calf Panel		1
2	SP1450446	Calf Pad Holder Right Turnable, cpl.		1
2	SP1450447	Calf Pad Holder Left Turnable, cpl.		1
3	SPZ173178	Spring		1
4	SPZ174134	Floating Cylinder With Pin		1
5	SP1450410	Knee Upholstery LH		1
6	SP1450448	Knee Upholstery RH		1
7	SP1449293	Ad-Cover cpl.		1
8	SPZ151006	Actuator With Circuit Board, Preassembled		1
9	SPN12196	Spring		1
10	SP1449408	Hanger ADP RH cpl.		1
10	SP1449407	Hanger ADP LH cpl.		1

Calf Plate Holder LH/RH 07.2005 ==>




Pos.	Part number	Description	Valid from → until	Qty
1	SP1450447	Calf Pad Holder Left Turnable, cpl.		1
1	SP1450446	Calf Pad Holder Right Turnable, cpl.		1
2	SP1487409	Calf Plate Holder for Recaroseat RH		1
2	SP1487410	Calf Plate Holder for Recaroseat LH		1
3	SP1449021	Assembly Kit, Calf Plate Holder		1
4	SP1450226	CALF PLATE HINGE, RH		1
4	SP1450225	CALF PLATE HINGE, LH		1



Pos.	Part number	Description	Valid from → until	Qty
1	1449018-0601	Legrest cpl. RH Junior		1
2	1449019-0601	Legrest cpl. LH Junior		1

Footplate flip up 01.08.10 ==>


1529644X.iso



Pos.	Part number	Description	Valid from → until	Qty
1	SP1541575	Footplate flip up, LH, cpl.		1
2	SP1541574	Footplate flip up, RH, cpl.		1
3	SP1529645	FOOTPLATE LH		1
4	SP1529644	FOOTPLATE RH		1
5	SP1543045	Plug-in tube		1
6	SP1423202	HEEL STRAP CPL.		1


Footplate flip up, angle-adjustable 01.08.10 ==>

1529645X.iso




Pos.	Part number	Description	Valid from → until	Qty
1	SP1541580	Footplate flip up, angle adjustable, LH, cpl.		1
2	SP1541579	Footplate flip up, angle adjustable, RH, cpl.		1
3	SP1529645	FOOTPLATE LH		1
4	SP1529644	FOOTPLATE RH		1
5	SP1423204	FOOT PLATE TUBE, ANGLE ADJUSTMENT LH CPL.		1
6	SP1423203	FOOT PLATE TUBE, ANGLE ADJUSTMENT RH CPL.		1
7	SP1423202	HEEL STRAP CPL.		1

Legrests, manual




Pos.	Part number	Description	Valid from → until	Qty
1	SP1442852	Lower Legrest Standard With Heel Strap, Width 43-45,5cm, Left		1
1	SP1442851	Lower Legrest Standard With Heel Strap, Width 43-45,5cm, Right		1
2	SP1442853	Lower Legrest With Angle Adjustable Footplate, With Heel Strap. Width		1
2	1428368	FOOTREST 24ALU2003,LOWER PART,FIXED,LONG,RIGHT		1
2	SP1442854	Lower Legrest With Angle Adjustable Footplate, With Heel Strap. Width		1
2	1428367	FOOTREST 24ALU2003,LOWER PART,FIXED,LONG,LEFT		1
3	1508945	HANGER, CPL LEFT		1
3	1513007	HANGER ALU2008, AA, LEFT, CPL		1
3	1513006	HANGER ALU2008, AA, RIGHT, CPL		1
3	1508948	HANGER, CPL RIGHT		1
4	1518532	Calf Support 15		1
4	1516314	CALF PAD 15 (165x170x80), CPL		1
5	1493055-0509	HANGER ALU2003,FIXED 80 DEG,CPL		1
5	1424577-0509	HANGER STANDARD, USABLE RIGHT AND LEFT ==>08/2006 - ==>08/2006		1
6	SP1441527	Leg Support Carrier RH cpl.		1
6	SP1441526	Leg Support Carrier LH cpl.		1
7	1430402	HEELSTRAP COMPLETE		1
8	1440874	Part not for sale		1
9	1430117	LEGRESTS STORM 3 EURO CPL. (PAIR)		1

Footrest Details




Pos.	Part number	Description	Valid from → until	Qty
1	SP1423202	HEEL STRAP CPL.		1
2	Z173217	AT 2000 footplate		1
2	Z173218	AT 2000 footplate		1
3	SP1423204	FOOT PLATE TUBE, ANGLE ADJUSTMENT LH CPL.		1
3	SP1423203	FOOT PLATE TUBE, ANGLE ADJUSTMENT RH CPL.		1
4	5106376	FOOTPLATE 18 AA, CPL.		1
5	5321579	HEEL LOOP, PAIR		1

Side Panel for Standard-Seat ==> 07/31/08


Pos.	Part number	Description	Valid from → until	Qty
b	1450095	Assembly Kit, Side Part Plate, Soft		1
a	1450296	Fixing Kit Side Part Plate		1
1	1442618	Armrest Pad With Assembly Kit - LH cpl.		1
1	1442619	Armrest Pad With Assembly Kit - RH cpl.		1
2	1442021-9005	Arm Pad Bracket RH		1
2	1442020-9005	Arm Pad Bracket LH		1
3	1442012-9005	Side Part Adjustment RH		1
3	1442015-9005	Side Part Adjustment LH		1
4	L10024	Wing Thumb Screw M6 X 16		1
5	1442024-9005	Armrest Support Tube LH		1
5	1442025-9005	Armrest Support Tube RH		1
6	1433292	Skirt-Guard RH		1
6	1433291	Skirt-Guard LH		1
7	1450102	Armrest Pad, Soft, Right With Assembly Kit		1
7	1450100	Armrest Pad, Soft, Left With Assembly Kit		1
8	SP1430417	Bag For Side Panel LH cpl.		1
8	SP1430416	Bag For Side Panel RH cpl.		1
9	1448749-9005	Base Plate cpl. Dragon		1
10	-	Washer DIN 5,3 9021		1
11	-	Hexagon Bolt DIN 912 M5X16		1
12	L10024	Wing Thumb Screw M6 X 16		1
13	SP1448754	Side Part Soft, Right cpl.		1
13	SP1448753	Side Part Soft, Left, cpl.		1
14	SP1442017	Side Part RH cpl.		1
14	SP1442016	Side Part LH cpl.		1

Side Panel for Standard-Seat 08/01/08 ==>




Pos.	Part number	Description	Valid from → until	Qty
1	SP1526440	Sidepart RH		1
2	SP1526441	Sidepart LH		1
3	SP1528846	Armpad, long, LH (376mm)		1
3	SP1528845	Armpad, long, RH (376mm)		1
4	SP1443004	Armpad LH (326mm)		1
4	SP1449043	Armpad RH (326mm)		1
5	SP1541505	Skirt guard, RH		1
5	SP1541506	Skirt guard, LH		1
6	SP1541507	T-Adapter		1
7	SP1541508	Insertion tube, RH		1
7	SP1541509	Insertion tube, LH		1
8	SP1541510	Siderestbracket, RH		1
8	SP1541511	Siderestbracket, LH		1

Armrest Junior Seat




Pos.	Part number	Description	Valid from → until	Qty
1	SP1448702	Side part RH cpl. Dragon Junior		1
2	5009885	Armrest Pad Comet		1
3	D36112	Self-Tapping Cross Head Screw 4,8 X 16		1
4	SP1448703	Side part LH cpl. Dragon Junior		1
5	1425601	Armrest Skirtguard 1		1
6	L10024	Wing Thumb Screw M6 X 16		1
7	D90184	Screw M6X12 DIN 7380		1

Front Wheels



Pos.	Part number	Description	Valid from → until	Qty
-	M12017	Castor Inner Tube 200x50, kpl		1
-	1446709	TYRE 220 X 50, GREY		1
1	1442604	Front Fork cpl.		1
1	SP1523435	Front Fork cpl. (Frame Melone Yellow)		1
2	1441999	Wheel 220X50 Pneumatic Grey		1
2	1442607	Wheel 220X50 Puncture Proof (Solid Inlay) Grey		1

Rear Wheel




Pos.	Part number	Description	Valid from → until	Qty
5	SP1440764	Pneumatic Rear Wheel Complete		1
5	SP1442858	Puncture-Proof Rear Wheel Complete (Solid Inlay)		1
5	SP1488561	Rearwheel Pneumatic, cpl. Dragon Start		1
5	SP1488562	Rearwheel Puncture Proof, cpl. Dragon Start		1
1	M12004	Inner Tube (for Tyre - 12 1/2" x 2 1/4")		1
1	23009262	Puncture Proof Infill (Complete with tyre)		1
2	M11005	Tyre 12 1/2 X 2 1/4		1
4	1440294	Wheel Fixing Kit Wheel 12-1/2		1
3	1442602	Wheel Rim Assembly Kit Black, cpl.		1
3	SP1523436	Rim Assembly-Kit		1


This diagram illustrates the assembly of a hydraulic cylinder with a lever arm. The components are numbered as follows:

- 1**: Hydraulic cylinder body.
- 1a**: Piston rod end cap.
- 2**: Flange mounting bracket.
- 3**: Piston rod end fitting.
- 4**: Main lever arm.
- 5**: Lever arm bracket.
- 6**: Lever arm bracket (alternative view).
- 7**: Lever arm bracket (alternative view).
- 8**: Piston rod end seal.
- 9**: Lever arm bracket (alternative view).
- 10**: Lever arm bracket (alternative view).

The diagram shows the assembly sequence with dashed lines indicating the alignment of parts. Red dashed boxes highlight specific sub-assemblies: the flange (2) and seal (8) on the cylinder end; the lever arm bracket (5) and its associated fittings (3, 4, 6, 7, 9, 10).

Motor & Declutching 6 km/h 16.07.08 ==>

Pos.	Part number	Description	Valid from → until	Qty
i	-	For 6km/h Version Effective from 16.07.2008 new Software.		1
-	SP1530172	Rework-Kit Schmid ==> Hua-Feng (without motors)		1
1	1440643-9005	Disengager		1
2	1430139	KNUCKLE JOINT		1
3	1509571	Gas spring		1
4	1442019	Starlock Cap 8 mm		1
5	1430138	Angle Joint M6 R		1
6	5006742	PULL-LATCH SMALL		1
7	SP1509652	Retro-Fitting Kit Declutching Mechanism (Lever, 2 Gassprings)		1
8	1436633	Motor 6km/h LH <i>From 16.07.08 new software</i>		1
8	1436632	Motor 6km/h RH <i>From 16.07.08 new software</i>		1
9	1509053-0002	Wheel Hub		1


Motor & Declutching 8 km/h 16.07.08 ==>

Pos.	Part number	Description	Valid from → until	Qty
-	SP1530172	Rework-Kit Schmid ==> Hua-Feng (without motors)		1
1	1440643-9005	Disengager		1
2	1430139	KNUCKLE JOINT		1
3	1509571	Gas spring		1
4	1442019	Starlock Cap 8 mm		1
5	1430138	Angle Joint M6 R		1
6	5006742	PULL-LATCH SMALL		1
7	SP1509652	Retro-Fitting Kit Declutching Mechanism (Lever, 2 Gassprings)		1
8	1440652	Drive Unit 8km/h Right		1
8	1440653	Drive Unit 8km/h Left		1
9	1509053-0002	Wheel Hub		1




20050801 →
20050801 →

REM24SD Remote




Pos.	Part number	Description	Valid from → until	Qty
1	1436541	ACS REM24SD Joystick		1
2	SP1497416	casing bottom part for REM24SD		1
3	SP1497417	casing upper part for REM24SD with keypad		1
4	SP1497418	Internal Cable Loom for REM 24SD		1
5	SP6001771	Joystick upper part		1
6	SP1497419	Joystick for REM 24SD		1

Remote holder for accompanying person





Pos.	Part number	Description	Valid from → until	Qty
1	SP1422937	Assembly Kit F. Std. Attendant Seat Control LH .		1
1	SP1422936	Assembly Kit RH F. Std. Attendant Seat Control		1
2	SP1422818	Assembly Kit For Attendantcontrol (Standard/Kontour/Flex)		1
3	-	Handknob		1
4	Z12008A.P40	Gripping Clip With Hexagon Hole, Ral 9005 Black		1



Pos.	Part number	Description	Valid from → until	Qty
1	SP1523444	Remote REM550		1
2	1440440-9005	JOYSTICK BOX SUPPORT SHARK, WITHOUT SCREWS		1
3	SP1498222	Installation-Kit For Control Unit Support Swing-Away Mechanism REM 24 SD		1
4	SP1535082	Joystickknob, cpl.		1

Keypads Shark2


1537479X.dwg



Pos.	Part number	Description	Valid from → until	Qty
1	SP1537479	Keypad DK-REMD01		1
2	SP1537480	Keypad DK-REMD11		1
3	SP1537481	Keypad DK-REMD21		1


Dual Control

6001425.dwg




Pos.	Part number	Description	Valid from → until	Qty
-	6001425	ACS Dual Control		1
1	P41105	Bus Cable Socket		1
2	GME61927	Fixing Clamp For ACS Bus Connector		1
3	P10304	ACS Joystick With Connecting Cable		1
4	Z183001A	Rubber Gaiter For ACS Remote		1

Remote Brackets




Pos.	Part number	Description	Valid from → until	Qty
-	SP1422818	Assembly Kit For Attendantcontrol (Standard/Kontour/Flex)		1
1	SPS183040	Remote Bracket Standard		1
2	1442777	Remote Support Shark, cpl. With Screws		1
3	SP1421796	Swing-Out Height-Adjustable Remote Holder LH		1
3	SP1421795	Swing-Out Height-Adjustable Remote Holder RH		1
4	SP1436542	Joystick Box Support Std. For REM 24SD, With Screws		1
5	SP1436543	Joystick Box Support Swing Away REM 24SD & Shark, With Screws		1
6	SP1496249	Remote Bracket, Height Adjustable And Swing Away LH		1
6	SP1496248	Remote Bracket, Height Adjustable And Swing Away RH		1
7	SP1536180	Montagekit for swing-away remote		1
7	SP1498222	Installation-Kit For Control Unit Support Swing-Away Mechanism REM 24 SD		1




Pos.	Part number	Description	Valid from → until	Qty
	1444717	Shark-Attendant Control, Switchable		1
	1484911	Shark2 Buscable L=2000mm	20051101 →	1
-	SP1419695	Magnet Key REM24 With Pouch cpl.		1
-	P10062	ACS Programming Kit (Wizard Software)		1
	1449299	Shark Retrofit Kit Operation Accompanying Person, Switchable		1
-	6001450	ACS-Hand Programmer, Instruccions In French, German, Spanish		1
1	AT.1448743	replacement part for Remote Shark REMB11 (drive and adjustments) incl. BUS cable 1.5m	20060501 →	1
1	AT.1441612	replacement part for Remote Shark REMA05 (drive only) incl. BUS cable 1.5m	20060501 →	1
1	1441612	Remote Shark With Bus Cable 1,5m	→ 20051201	1
2	1484904	Remote DK-REMD21, Shark2 (Drive, Seat Adjustment, Light)	20051101 →	1
2	1449177	Remote Shark DK-REMD11 (with actuators)	20051101 →	1
2	1492567	Remote Shark DK-REMD01 (drive only)	20051101 →	1
3	AT.1440541	Comando Luzes Shark (Not longer available. Please order exchange-kit)		1
3	SP1522494	Exchange-Kit For Shark1 Lightmodule (At.1440541) Standard-Remote-Holder		1
4	SP1440312	Protective Bumber For Shark Remote cpl.		1
5	1442691	Guard Bow For Shark With Lights Remote cpl.		1
6	1441551	Programming Adapter Shark (Unconditionally Necessary To Program Shark)		1
7	1436541	ACS REM24SD Joystick		1
8	6001425	ACS Dual Control		1
9	SP1422828	Protective Bumber F. REM 24 Remote cpl.		1

Remoteholder



Pos.	Part number	Description	Valid from → until	Qty
	SP1422818	Assembly Kit For Attendantcontrol (Standard/Kontour/Flex)		1
1	SPS183040	Remote Bracket Standard		1
2	1442777	Remote Support Shark, cpl. With Screws		1
3	SP1421796	Swing-Out Height-Adjustable Remote Holder LH		1
3	SP1421795	Swing-Out Height-Adjustable Remote Holder RH		1
4	SP1436542	Joystick Box Support Std. For REM 24SD, With Screws		1
5	SP1436543	Joystick Box Support Swing Away REM 24SD & Shark, With Screws		1

Batterys & Charger



Pos.	Part number	Description	Valid from → until	Qty
-	SP1523354	Battery Fixation-Kit (not shown)		1
-	1440758	Battery Cable Shark		1
-	1442552	Battery Cable Shark With Light		1
-	1444668	Battery Cable ACS		1
-	SP1532880	Batterycable ACS2		1
6,7	SP1423696	Charger 24V/8A With Continental Mains Cable, cpl.		1
6,8	SP1423697	Charger 24V/8A With Uk Mains Cable, cpl.		1
1	1432944	Battery Mk Sldg 60Ah (Gel)		1
1	SP1522107	AGM Battery 50Ah		1
2	1421655	Battery Strap		1
3	P52996	Battery Terminal Cover		1
4	1440759	Strip Fuse Holder With 40A Fuse Shark		1
4	1444669	Strip Fuse Holder With Fuse 50A		1
4	1432561	Strip Fuse Holder with 80A Fuse <i>Only for Shark II with 90Ah Powermodule</i>		1
5	1442600	Strip Fuse 40A (Shark)		1
5	P12191	Fuse 50A		1
6	SPP51981	Gel Battery Charger Unit		1
7	P21079	Charger Power Cable (Shockproof Plug)		1
8	P21080	UK 3 Pin Mains Lead		1
8	1466783	Power Cable, 3-pole, black, 6A/250VAC (Switzerland)		1
9	5813987	Battery Charger 8A Gi Popular (only Sweden)		1
10	1492498	Battery Frame ACS		1

Cables

GSM63003DR.dwg

ACS-Buskabel

Länge / length (mm)	Art.-Nr.	Pos.
300	-	1
1100	-	2
1600	-	3

Shark-Buskabel

Länge / length (mm)	Art.-Nr.	Pos.
1000	-	4
2000	-	5
Y-Bus-Cable	-	6

Kabel für Licht- und Aktuatormodule (CLAM)

Cable for Light- and Actuatomoduls (CLAM)


Licht / Light	Sitzwinkelverstellung Seatangleadjustment	Rückenverstellung Backrestadjustment	Sitzhöhenverstellung Seathightadjustment	Beinstützenverstellung Legrestadjustment	Art.-Nr.	Pos.
-	-	x	x	x	-	7
x	x	x	-	x	-	8
x	-	x	x	x	-	9

Batterie kabel / Battery cable

Licht / Light	mit Sicherung with Fuse	Art.-Nr.	Pos.
24V	-	-	10
24V	50A	-	11
24V (Shark)	-	-	12
24V (Shark)	40A	-	13


Pos.	Part number	Description	Valid from → until	Qty
1	GSM63003	ACS Bus Cable L=300 mm		1
2	GSM63011	ACS Bus-Cable 1,1M		1
3	GSM63016	ACS Bus-Cable 1,6M		1
4	1448922	Shark Buscable 1000mm		1
5	1484911	Shark2 Buscable L=2000mm		1
6	1450314	Shark Y-Bus Cable		1
7	1445554	Harness Shark For Electrical Adjustments		1
8	1444670	Harness ACS For Electrical Adjustments And Lightening		1
9	1488951	Harness Shark For Ligth And Electrical Adjustments		1
10	1444668	Battery Cable ACS		1
11	1444669	Strip Fuse Holder With Fuse 50A		1
12	1440758	Battery Cable Shark		1
13	1440759	Strip Fuse Holder With 40A Fuse Shark		1

Front Lights




Pos.	Part number	Description	Valid from → until	Qty
-	1442734	Headlamp RH 24V cpl.		1
-	1442735	Headlamp LH 24V cpl.		1
1	1433737	Glass For Front Light		1
2	P31108	Light Bulb 24V/10W		1
3	1419271	Light Bulb 24V/21W Yellow, turn indicator		1
4	1428082	Swivel Mechanism RH		1
4	1428081	Swivel Mechanism LH		1
5	1442685	Mounting Kit Front Light		1

Rear Light




Pos.	Part number	Description	Valid from → until	Qty
b	-	Insulating Sheath		1
e	-	Hex Nut M5 Self-Locking		1
d	-	Washer 5.3 DIN 125		1
a	-	Screw DIN 7991 M8X20 Black Galvanised		1
-	SP1441571	Light Support With Rearlight LH cpl.		1
c	-	Cable Fixation Small		1
-	SP1441565	Light Support With Rearlight RH cpl.		1
1	1442736	Rear Light 24V, cpl.		1
2	1419271	Light Bulb 24V/21W Yellow, turn indicator		1
3	P31191	Lamp 24V/5W		1
4	5807542	REFLECTOR, RED		1
5	1441573-9005	Light Support Rear LH		1
5	1441569-9005	Light Support Rear RH		1
6	1441568	Fixing For Reflector		1
7	1419274	Reflector Red, Round		1
8	1442603	Tail Light Fixing Kit		1

Belts

Pos.	Part number	Description	Valid from → until	Qty
1	1422869	Positioning Strap With Velcro cpl. for Junior-Seat		1
2	1422845	Positioning Strap With Snap cpl. For Standard-Seat		1
2	5370088	Belt Dragon Start		1
3	SP1422848	Holding Strap For Standard And Contour Seat H-Form cpl.		1

Headrest


1422865DR.dwg



Pos.	Part number	Description	Valid from → until	Qty
1	SP1422925	Headrest Standard/Contour Black cpl.		1
2	SP1422865	Headrest Cushion Contour Black, Wo Bracket		1
3	SP5313308	Headrest Cushion Junior, Grey Vinyl, cpl.		1
4	SP5313262	Headrest Cushion Junior, Grey Vinyl, Wo Bracket		1
5	1440900	Headrest Comfort (black)		1
6	1440903	Headrest Pad Comfort (black)		1
7	SP1422547	Headrest Bracket For Contour-/Flex Seat cpl.		1
8	SP1422914	Headrest Bracket For Standard Seat cpl.		1

Headrest with Cheeksupport


1446560X.dwg



Pos.	Part number	Description	Valid from → until	Qty
1	1446560	Cheek Support, Sidewise Adj.		1
1	SP1541658	HEAD/CHEEK SUPP. TR26 CPL. SIDEWISE ADJUSTABLE FOR EASY ADAPT		1
2	1444210	HEAD REST ATTACHM. CPL.		1
3	1444304	CUSHION, CHEEK SUPPORT		1
4	1445691	ATTACHMENT HEADREST CPL		1
5	1444183-0509	POLE NECKREST, CPL		1
5	SP1529572	Headrestholder for EasyAdapt		1
6	5716497	PLASTIC CAP		1
7	1448661	Stop cpl., Neck- & Headrest		1
8	1445838	SKENA HUVUDSTÖD KPL.		1
9	1485538	Part not for sale		1

Options


1422955.dwg



Pos.	Part number	Description	Valid from → until	Qty
-	1449280	USB Wizard Dongle Kit, EDL-USB, cpl.		1
1	SP1422955	Rear Mirror		1
2	SP1423435	Rear Mirror Support Rail cpl.		1
3	SP1422881	Rear Mirror Complete With Support		1
4	SP1430416	Bag For Side Panel RH cpl.		1
4	SP1430417	Bag For Side Panel LH cpl.		1
5	SP1422844	Cane Holder Std. cpl., RH/LH		1
6	5370157	Hour counter		1
6	SP1492110	Cabele for hour counter		1


Tray

Z177036.dwg



Pos.	Part number	Description	Valid from → until	Qty
-	SPN377001	Tray RH Foldaway With Holder With Cut-Out For Remote LH, Junio-Seat		1
-	SPN377002	Tray LH Foldaway With Holder With Cut-Out For Remote RH, Junior-Seat		1
-	SP1422954	Swing-Out Tray Without Cut-Out, cpl.		1
-	SP1422952	Swing-Out Tray With Cut-Out F. Control LH cpl.		1
-	SP1422953	Swing-Out Tray With Cut-Out F. Control RH cpl.		1
1	Z177016	Tray-RH Control, Transparent		1
1	1415083	Tabletop Cut RH		1
1	Z177017	Tray-LH Control, Transparent		1
1	Z17024A	Tray 660 X 540 X 8, Transparent		1
1	1415084	Tabletop Cut LH		1
2	-	Tube End Plug SI 22		1
3	-	Sks Nut Cap M12		1
4	-	Handknob		1
5	-	Bush Rd. 19 X 3.5 X 25		1
6	-	Tube End Plug E 15 X 15		1
7	SPZ177011	Tray holder		1

Kerb Climber



Pos.	Part number	Description	Valid from → until	Qty
	SP1434453	Curbclimber With Front Fixation cpl.		1
1	SP1418452	Kerb Climber cpl.		1
2	1434454-9005	Kerbstone Raiser Receptacle With Fixing Brackets		1
3	N10499	Gas Spring		1
4	1433733	Kerb Climber With Shoe		1
5	N02201	Bush BB 1517 DU		1
6	L01039	Tube End Plug F 38 X 20		1



Yes, you can.®



Belgium & Luxemburg:

Invacare nv, Autobaan 22, B-8210 Loppem
Tel: (32) (0)50 83 10 10, Fax: (32) (0)50 83 10 11
belgium@invacare.com

Danmark:

Invacare A/S, Sdr. Ringvej 37, DK-2605 Brøndby
Tel: (45) (0)36 90 00 00, Fax: (45) (0)36 90 00 01
denmark@invacare.com

Deutschland:

Invacare Aquatec GmbH, Alemannenstraße 10, D-88316 Isny
Tel: (49) (0)75 62 7 00 0, Fax: (49) (0)75 62 7 00 66
info@invacare-aquatec.com

European Distributor Organisation:

Invacare, Kleiststraße 49, D-32457 Porta Westfalica
Tel: (49) (0)57 31 754 540, Fax: (49) (0)57 31 754 541
edo@invacare.com

España:

Invacare SA, c/Areny s/n, Polígon Industrial de Celrà, E-17460 Celrà (Girona)
Tel: (34) (0)972 49 32 00, Fax: (34) (0)972 49 32 20
contactsp@invacare.com

France:

Invacare Poirier SAS, Route de St Roch, F-37230 Fondettes
Tel: (33) (0)2 47 62 64 66, Fax: (33) (0)2 47 42 12 24
contactfr@invacare.com

Ireland:

Invacare Ireland Ltd, Unit 5 Seatown Business Campus, Seatown Road, Swords, County Dublin - Ireland
Tel: (353) 1 810 7084, fax: (353) 1 810 7085
ireland@invacare.com

Italia:

Invacare Mecc San s.r.l., Via dei Pini 62, I-36016 Thiene (VI)
Tel: (39) 0445 38 00 59, Fax: (39) 0445 38 00 34
italia@invacare.com

Nederland:

Invacare BV, Celsiusstraat 46, NL-6716 BZ Ede
Tel: (31) (0)318 695 757, Fax: (31) (0)318 695 758
nederland@invacare.com, csede@invacare.com

Norge:

Invacare AS, Grensesvingen 9, Postboks 6230, Etterstad, N-0603 Oslo
Tel: (47) (0)22 57 95 00, Fax: (47) (0)22 57 95 01
norway@invacare.com, island@invacare.com

Österreich:

Mobitec Mobilitätshilfen GmbH, Herzog Odilostrasse 101, A-5310 Mondsee
Tel: (43) 6232 5535 0, Fax: (43) 6232 5535 4
office@mobitec-austria.com, austria@invacare.com

Portugal:

Invacare Lda, Rua Estrada Velha, 949, P-4465-784 Leça do Balio
Tel: (351) (0)225 1059 46147, Fax: (351) (0)225 1057 39
portugal@invacare.com

Sverige & Suomi:

Invacare AB, Fagerstagatan 9, S-163 91 Spånga
Tel: (46) (0)8 761 70 90, Fax: (46) (0)8 761 81 08
sweden@invacare.com, finland@invacare.com

Switzerland:

Mobitec Rehab AG, Benkenstrasse 260, CH-4108 Witterswil
Tel: (41) (0)61 487 70 80, Fax: (41) (0)61 487 70 81
office@mobitec-rehab.ch, switzerland@invacare.com

United Kingdom:

Invacare Limited, Pencoed Technology Park, Pencoed, Bridgend CF35 5HZ
Switchboard Tel: (44) (0)1656 776200, Fax: (44) (0)1656 776201
Customer services Tel: (44) (0)1656 776222, Fax: (44) (0)1656 776220
UK@invacare.com