


Customer:  
 Customer Code:  
 Phone:

Contact:  
 Commission:

ALIZE stool with arms H2240/1



Pos.	Qty	Part Number	Designation
1	1	1410052	Seat
2	6	1410491	Form Washer
3	6	C 1001	Screw cap
4	6	-	M8 x 50 Screw
5	4	1410505	ANTI FOLDING CLIP
6	4	1423382	Plastic Int Tube
7	4	1410068	Adjustable feet
8	4	1410495	Frame Tip

**Shower Stool with arms H2240/1**



Alize Range

**Edition  
 04-02**