






Quick Instruction for Use for Invacare® Cushions

Matrix® and Flo-Tech® Range



www.invacare.eu.com



		kg	kg	mm	mm	mm	°C
Matrix®	Libra	2.0 kg	136; 227 kg	355 - 560 mm	305 - 560 mm	85 mm	80 °C
	Contour Visco NG	2.5 kg	80 - 140 kg	355 - 555 mm	405 - 530 mm	100 mm	90 °C
Flo-Tech®	Lite	0.9 kg	127 kg	255 - 510 mm	255 - 535 mm	80 mm	80 °C
	Lite Visco	1.6 kg	32 - 127 kg	255 - 510 mm	255 - 535 mm	80 mm	80 °C
	Contour	1.5 kg	152; 191 kg	355 - 610 mm	330 - 560 mm	115 mm	80 °C
	Contour Visco	1.8 kg	51 - 127 kg	355 - 510 mm	330 - 535 mm	115 mm	80 °C
	Plus	2.4 kg	152; 191 kg	355 - 610 mm	330 - 560 mm	115 mm	80 °C
	Image	1.9 kg	140 kg	355 - 510 mm	330 - 535 mm	80 mm	80 °C
	Solution	3.5 kg	140 kg	355 - 610 mm	330 - 560 mm	80 mm	80 °C
	Solution Xtra	3.5 kg	140 kg	355 - 610 mm	330 - 560 mm	80 mm	80 °C



Invacare UK Operations Limited
Pencoed Technology Park
Pencoed Bridgend CF35 5AQ
UK

Invacare Portugal Unipessoal, Lda
Rua Estrada Velha 949
4465-784 Leça do Balio
Portugal



Yes, you can.