



C E R T I F I C A T E

No. Z1 03 05 11048 005

Holder of Certificate: Ulrich Alber GmbH + Co. KG

Vor dem Weißen Stein 21
72461 Albstadt
Germany

Certification Mark:



Product: Wheelchair for handicapped

The product was tested on a voluntary basis and complies with the essential requirements.
The certification mark shown above can be affixed on the product. See also notes overleaf.

Test report no.: 70046021



Date, 2003-06-02

Page 1 of 2

A handwritten signature in blue ink, appearing to read 'J. Meyer', is written over the bottom portion of the certification mark.

CERTIFICATE
No. Z1 03 05 11048 005



Model(s): **Electrically driven wheelchair**
Type: Adventure A10

Parameters:

class of use:	class B, rear-wheel drive
total height:	935 mm
total length:	1172 mm
total width:	680 mm
max. load:	140 kg
electro-motor:	brushless, electronically commutate d.c. motors
driving controller:	Adventure A10
batteries:	Frank Electronic GmbH
	type: MP22-12 C (4 x 12 V, 22 Ah)
battery charger:	Soneil, model: 2412SRF-B, 24 V DC, 6 A

Tested according to: EN 12184:1999

Production Facility(ies): 11048