



Product Service

# CERTIFICATE

No. Z1 05 09 15775 039

**Holder of Certificate:** AQUATEC GmbHAlemannenstr. 10  
88316 Isny  
GERMANY**Production Facility(ies):** 15775**Certification Mark:****Product:** Sanitary Aids for Disabled Persons**Model(s):** Toilet- and shower chair  
Type: Aquatec OCEAN VIP

<b>Parameters:</b>	max. width:	576 mm
	max. depth:	1058 mm
	seat angle adjustment:	0° - +35°
	seat width:	500 mm
	seat height:	500, 550, 600 mm
	seat depth:	452 mm
	distance between the armrests:	457 mm
	wheel size (front/ rear):	ø 125 x 30 mm (all brakable)
	weight:	21,5 kg
	frame material:	stainless steel / plastic
	max. load:	150 kg

**Tested according to:** EN 12182:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. See also notes overleaf.

**Test report no.:** 70102645**Date,** 2005-09-22

Page 1 of 1

