


DOLOMITE MAXI +



- | | | |
|----|----------|---|
| 1 | 13565-25 | SL II-wheel complete in front fork (pcs) |
| 2 | 12590-25 | SL II-wheel complete in rear fork (pcs) |
| 3 | 12609 | SL II-wheel kit, complete (pcs) |
| 4 | 70050 | Cover (pcs) |
| 5 | 12620 | Brake flap, complete DABSII (pcs) grey |
| 6 | 13096 | Rear fork, DABSII (pcs) |
| 7 | 12708-M | Cross bar spacer |
| 8 | 12627-10 | Folding hinge kit (pair) |
| 9 | 12707-M | Folding lock - kit (pcs) |
| 10 | 12665-M | Knob and screw for heightadjustment (pair) |
| 11 | 80060 | Wire (pcs) |
| 12 | 80094 | Wire cover L=820 mm (pcs) |
| 13 | 13238-L | Handle kit DABS, left. Incl. grip and wire (pcs) |
| 13 | 13238-R | Handle kit DABS, right. Incl. grip and wire (pcs) |
| 14 | 12677 | Soft grip, L=220 mm (pair) |
| 15 | 12509 | Anatomic grips, incl. Cover (pair) |
| 16 | 13464-L | Brake handle kit, left incl. wire (pcs) |
| 16 | 13464-R | Brake handle kit, right incl. wire (pcs) |
| 17 | 12721-M | Flip-up seat, compl w. attachmentparts (pcs) |
| 18 | 70043 | Cable holder (pcs) |
| | 82249 | Operating instructions (pcs) |