



# CERTIFICATE

No. Z1 03 05 11048 005

**Holder of Certificate:** Ulrich Alber GmbH + Co. KG

Vor dem Weißen Stein 21  
72461 Albstadt  
Germany

**Certification Mark:**



**Product:** Wheelchair for handicapped

The product was tested on a voluntary basis and complies with the essential requirements.  
The certification mark shown above can be affixed on the product. See also notes overleaf.

**Test report no.:** 70046021



**Date,** 2003-06-02

Page 1 of 2

A handwritten signature in blue ink, appearing to read 'J. Meyer', is written over the bottom portion of the page.

**CERTIFICATE**

No. Z1 03 05 11048 005



**Model(s):** **Electrically driven wheelchair**  
**Type: Adventure A10**

**Parameters:**

class of use:	class B, rear-wheel drive
total height:	935 mm
total length:	1172 mm
total width:	680 mm
max. load:	140 kg
electro-motor:	brushless, electronically commutate d.c. motors
driving controller:	Adventure A10
batteries:	Frank Electronic GmbH type: MP22-12 C (4 x 12 V, 22 Ah)
battery charger:	Soneil, model: 2412SRF-B, 24 V DC, 6 A

**Tested according to:** EN 12184:1999

**Production Facility(ies):** 11048